

# ISOLATION GOWN

LEVEL 1 LEVEL 2



## Technical Data Sheet

### Description

Level 1 and level 2 isolation gowns available in yellow and blue.  
Fluid resistant, full length, knit cuff, latex free and disposable.

### General Information:

<b>Size</b>	
GAH-GISO1YL (Level 1)	Large
GAH-GISO1YX (Level 1)	X-Large
GAH-GISO2BL (Level 2)	Large
GAH-GISO2BX (Level 2)	X-Large
<b>Material</b>	Spunbonded Polypropylene
Level 1	Laminated Spunbonded
Level 2	Polypropylene
<b>Color</b>	Yellow
Level 1	
Level 2	Blue
<b>Additional Details</b>	No pockets
<b>Packaging</b>	10/Bag, 5 Bags/ Case QTY 50



Level 1



Level 2

### TESTING: Level 1 Gown:

<b>Liquid Barrier</b>	ANSI/AAMI PB70
<b>Water Resistance:</b> Impact Penetration Test	AATCC 42
<b>Water Resistance:</b> Hydrostatic Pressure Test	AATCC 127

### TESTING: Level 2 Gown:

<b>Liquid Barrier</b>	ANSI/AAMI PB70
<b>Water Resistance:</b> Impact Penetration Test	AATCC 42
<b>Water Resistance:</b> Hydrostatic Pressure Test	AATCC 127

The above is typical physical property test data applicable to the sample reference name and should not be construed as a specification.

**Disclaimer:** The information supplied in this document is for guidance only and should not be construed as a warranty. All implied warranties are expressly disclaimed, including without limitation any warranty of merchantability and fitness for use. All users of the material are responsible for assuring that it is suitable for their needs, environment and end use. All data is subject to change as High-Tech Conversions deems appropriate.

