

# smiths medical

## portex®



### BLUselect® Suctionaid® Tracheostomy Tubes

Additional Lumen For Effective  
Suctioning Above The Cuff



# BLUselect® Suctionaid® Tracheostomy Tubes

Each Suctionaid® tube is equipped with an additional lumen to permit effective removal of secretions above the cuff, minimising the risk of aspiration. The suction lumen can be used to apply continuous or intermittent suction (thumb valve included) above the cuff. Change kits include a J-tip guidewire to help ensure the correct placement of the replacement tube.



## Features include:

- Integral 15 mm connector allows for ventilation with or without inner cannula in place
- Thermosensitive material with sufficient initial rigidity for insertion which, when heated to body temperature then conforms to the individual's upper respiratory tract ensuring minimal trauma
- Special obturator clip design minimises obturator movement during tube insertion
- Easy to read laser printing on flanges
- Color coding indicating size on tubes and packages
- Not manufactured with DEHP plasticizers
- SOFT-SEAL® cuff — low pressure, high volume cuff for minimal trauma while providing effective seal:
  - When deflated, the cuff fits smoothly to the tube minimising discomfort during insertion and removal
  - 30 mm cuff resting diameter on most common adult sizes (8,9 and 10 mm ID)
- Markings on pilot balloon provide easy-to-read with relevant tube and cuff dimensional information

## Ordering Information

Cuffed BLUselect® Suctionaid® Tracheostomy Tubes				
Reference Code	Tube Size	Tube ID	Tube OD	Tube Length
101/875/060	6.0 mm	6.0 mm	9.2 mm	64.5 mm
101/875/070	7.0 mm	7.0 mm	10.5 mm	70.0 mm
101/875/080	8.0 mm	8.0 mm	11.9 mm	75.5 mm
101/875/090	9.0 mm	9.0 mm	13.3 mm	81.0 mm
101/875/100	10.0 mm	10.0 mm	14.0 mm	87.5 mm

Cuffed tubes include an obturator, two inner cannulae, suction control valve, Velcro® trach tube holder, disconnect wedge, cleaning brush, and patient label.

Cuffed BLUselect® Suctionaid® Tracheostomy Tube Change Kits				
101/885/070	7.0 mm	7.0 mm	10.5 mm	70.0 mm
101/885/080	8.0 mm	8.0 mm	11.9 mm	75.5 mm
101/885/090	9.0 mm	9.0 mm	13.3 mm	81.0 mm

Cuffed tube change kits include an obturator, two inner cannulae, J-tip guide wire, suction control valve, Velcro® trach tube holder, disconnect wedge, cleaning brush, and patient label.

Inner Cannulas For BLUselect® Suctionaid® Tracheostomy Tubes					
Description	Inner Canula ID	Units Per Case	Reference Code	Units Per Cse	Reference Code
Inner Cannula for 6.0 mm Tube	5.0 mm	20	101/856/060	50	101/858/060
Inner Cannula for 7.0 mm Tube	6.0 mm	20	101/856/070	50	101/858/070
Inner Cannula for 8.0 mm Tube	7.0 mm	20	101/856/080	50	101/858/080
Inner Cannula for 9.0 mm Tube	8.0 mm	20	101/856/090	50	101/858/090
Inner Cannula for 10.0 mm Tube	9.0 mm	20	101/856/100	50	101/858/100

Please refer to the Instructions for Use supplied with the product for detailed instructions, warnings, and cautions.

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN ALL GEOGRAPHIES

**Smiths Medical ASD, Inc.**  
6000 Nathan Lane North  
Minneapolis, MN 55442, USA  
Tel: 1-614-210-7300  
Toll-Free USA: 1-800-258-5361  
www.smiths-medical.com

**Smiths Medical International Ltd.**  
1500 Eureka Park, Lower Pemberton  
Ashford, Kent, TN25 4BF  
Tel: +44 (0)845 850 0445

MMSPCA-0578

**Rx  
ONLY**

Find your local contact information at: [www.smiths-medical.com/customer-support](http://www.smiths-medical.com/customer-support)

**smiths medical**

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Portex, BLUselect and the Smiths Medical and design marks are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2019 Smiths Medical. All rights reserved. TR194554GB-082019