

# CADD<sup>®</sup> AMBULATORY TUBING SETS WITH FLOW STOP FREE FLOW PROTECTION

CADD <sup>™</sup> Medication Cassette Reservoirs for CADD <sup>®</sup> -Solis, CADD-Legacy <sup>®</sup> and CADD-Prizm <sup>®</sup> Pumps - Must be used with a CADD <sup>®</sup> extension set							
	Reorder No.	Units	Description	Tubing Length	Priming Vol.	Plasticizer	Previous No.
	21-7301	12/box	50 mL CADD <sup>™</sup> medication cassette reservoir with flow stop, clamp, female Luer and non-vented cap	8 in (20 cm)	0.2 mL	DEHP	21-7001
	21-7302	12/box	100 mL CADD <sup>™</sup> medication cassette reservoir with flow stop, clamp, female Luer and non-vented cap	8 in (20 cm)	0.2 mL	DEHP	21-7002
	21-7300	12/box	100 mL yellow CADD <sup>™</sup> medication cassette reservoir with flow stop, clamp, female Luer and non-vented cap	8 in (20 cm)	0.2 mL	DEHP	21-7100
	21-7308	12/box	250 mL clear CADD <sup>™</sup> medication cassette reservoir with flow stop, clamp, female Luer and non-vented cap	8 in (20 cm)	0.2 mL	TOTM (non-DEHP)	
	21-7309	12/box	250 mL yellow CADD <sup>™</sup> medication cassette reservoir with flow stop, clamp, female Luer and non-vented cap	8 in (20 cm)	0.2 mL	TOTM (non-DEHP)	
	21-7310	12/box	250 mL blue CADD <sup>™</sup> medication cassette reservoir with flow stop, clamp, female Luer and non-vented cap	8 in (20 cm)	0.2 mL	TOTM (non-DEHP)	
	21-7600	12/box	100 mL yellow CADD <sup>™</sup> medication cassette reservoir with NRFit <sup>™</sup> connector, flow stop, female Luer, clamp, and non-vented stopper	8 in (20 cm)	0.2 mL	DEHP	21-7100 21-7300
	21-7609	12/box	250 mL yellow CADD <sup>™</sup> medication cassette reservoir with NRFit <sup>™</sup> connector, flow stop, female Luer, clamp, and non-vented stopper	8 in (20 cm)	0.2 mL	TOTM (non-DEHP)	21-7309
	21-7664-24	50/box	Filling Adapter for Pharmacy use only				
	24-1001-64	10/box	Infusion Adapter to convert female Luer to NRFit <sup>™</sup> female				
	24-1009-64	50/box	Cap, female, yellow				
	24-1004-64	50/box	Cap, male, non-vented, yellow				

CADD <sup>®</sup> Administration Sets for CADD <sup>®</sup> -Solis, CADD-Legacy <sup>®</sup> and CADD-Prizm <sup>®</sup> Pumps							
	Reorder No.	Units	Description	Tubing Length	Priming Vol.	Plasticizer	Previous No.
	21-7321	12/box	CADD <sup>®</sup> administration set with female Luer, flow stop, clamp, and one-way checkvalve with male Luer (for use with a syringe)	102 in (259 cm)	2.7 mL	TOTM (non-DEHP)	21-7021
	21-7322	12/box	CADD <sup>®</sup> administration set with bag spike, flow stop, clamp, and one-way checkvalve with male Luer	78 in (198 cm)	2.3 mL	TOTM (non-DEHP)	21-7022 21-7034
	21-7323	12/box	CADD <sup>®</sup> blue-striped administration set with bag spike, flow stop, clamp, and one-way checkvalve with male Luer	78 in (198 cm)	2.3 mL	TOTM (non-DEHP)	21-7023
	21-7324	12/box	CADD <sup>®</sup> yellow-striped administration set with bag spike, flow stop, clamp, and one-way checkvalve with male Luer	123 in (312 cm)	3.2 mL	TOTM (non-DEHP)	21-7024
	21-7359	12/box	CADD <sup>®</sup> administration set with male Luer, flow stop, clamp, and one-way checkvalve with male Luer	69 in (175 cm)	2.0 mL	TOTM (non-DEHP)	21-7059
	21-7390	12/box	CADD <sup>®</sup> administration set with female Luer, flow stop, 7.6 cm y-extension, one-way checkvalve with female Luer, one-way checkvalve with male Luer, and clamps (for use with a syringe)	102 in with 3 in pigtail (259 cm with 7.6 cm pigtail)	4.0 mL	TOTM (non-DEHP)	21-7090
	21-7391	12/box	CADD <sup>®</sup> administration set with bag spike, flow stop, 7.6 cm y-extension, one-way check valve with female Luer, one-way check valve with male Luer, and clamps	108 in with 3 in pigtail (274 cm with 7.6 cm pigtail)	4.2 mL	TOTM (non-DEHP)	21-7091
	21-7624	12/box	CADD <sup>®</sup> yellow-striped administration set with NRFit <sup>™</sup> connector, bag spike, flow stop, clamp, and one-way checkvalve with male Luer	123 in (312 cm)	3.2 mL	TOTM (non-DEHP)	21-7024 21-7324

CADD <sup>®</sup> Administration Sets with Air-eliminating Filters for CADD <sup>®</sup> -Solis, CADD-Legacy <sup>®</sup> and CADD-Prizm <sup>®</sup> Pumps							
	Reorder No.	Units	Description	Tubing Length	Priming Vol.	Plasticizer	Previous No.
	21-7343	15/box	CADD <sup>®</sup> administration set with bag spike, flow stop, 1.2-micron air-eliminating filter, clamp, needleless injection site, and one-way check valve with male Luer	114 in (290 cm)	15.5 mL	TOTM (non-DEHP)	21-7333
	21-7346	15/box	CADD <sup>®</sup> administration set with bag spike, flow stop, 0.2-micron air-eliminating filter, spiral tubing, clamp, and one-way checkvalve with male Luer	94 in (239 cm)	10.0 mL	TOTM (non-DEHP)	21-7336 21-7036
	21-7363	15/box	CADD <sup>®</sup> administration set with bag spike, flow stop, 1.2-micron air-eliminating filter, spiral tubing, clamp, and one-way checkvalve with male Luer	93 in (236 cm)	8.0 mL	TOTM (non-DEHP)	21-7383 21-7083 21-7071
	21-7349	12/box	CADD <sup>®</sup> yellow-striped administration set with bag spike, flow stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male Luer	130 in (330 cm)	6.2 mL	TOTM (non-DEHP)	21-7339 21-7039
	21-7394	12/box	CADD <sup>®</sup> administration set with bag spike, flow stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male Luer	108 in (274 cm)	3.3 mL	TOTM (non-DEHP)	21-7094
	21-7395	12/box	CADD <sup>®</sup> administration set with female Luer, flow stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male Luer	102 in (259 cm)	3.1 mL	TOTM (non-DEHP)	21-7095
	21-7649	12/box	CADD <sup>®</sup> yellow-striped administration set with NRFit <sup>™</sup> connector, bag spike, flow stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male Luer	130 in (330 cm)	6.2 mL	TOTM (non-DEHP)	21-7639 21-7039 21-7339

CADD <sup>®</sup> High Volume Administration Sets for use with the CADD <sup>®</sup> -Solis and CADD-Prizm <sup>®</sup> Pumps							
	Reorder No.	Units	Description	Tubing Length	Priming Vol.	Plasticizer	Previous No.
	21-7345	15/box	CADD <sup>®</sup> high volume administration set with bag spike, flow stop, 0.2-micron air-eliminating filter, clamp, injection site and one-way checkvalve with male Luer	120 in (305 cm)	20 mL	TOTM (non-DEHP)	21-7355 21-7055 21-7055V
	21-7357	15/box	CADD <sup>®</sup> high volume administration set with bag spike, flow stop, clamp, and one-way checkvalve with male Luer	92 in (234 cm)	11 mL	TOTM (non-DEHP)	21-7057 21-7057V
	21-7364	15/box	CADD <sup>®</sup> high volume administration set with bag spike, flow stop, 1.2-micron air-eliminating filter, Luer-activated needleless injection site, clamp, and one-way checkvalve with male Luer	118 in (300 cm)	18 mL	TOTM (non-DEHP)	21-7375 21-7075 21-7075V
	21-7361	15/box	CADD <sup>®</sup> high volume administration set with bag spike, flow stop, 1.2-micron air-eliminating filter, clamp, injection site, and one-way checkvalve with male Luer	118 in (300 cm)	18 mL	TOTM (non-DEHP)	21-7381 21-7081 21-7389 21-7089
	21-7386	15/box	CADD <sup>®</sup> high volume administration set with bag spike, flow stop, 1.2-micron air-eliminating filter, clamp, and one-way checkvalve with male Luer	104 in (264 cm)	15 mL	TOTM (non-DEHP)	
	21-7384	15/box	CADD <sup>®</sup> yellow-striped high volume administration set with bag spike, flow stop, clamp, and one-way checkvalve with male Luer	123 in (312 cm)	16 mL	TOTM (non-DEHP)	
	21-7365	15/box	CADD <sup>®</sup> yellow-striped high volume administration set with bag spike, flow stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male Luer	136 in (345 cm)	23 mL	TOTM (non-DEHP)	21-7385 21-7085
	21-7684	15/box	CADD <sup>®</sup> yellow-striped high volume administration set with NRFit <sup>™</sup> connector, bag spike, flow stop, clamp, and one-way checkvalve with male Luer	123 in (312 cm)	16 mL	TOTM (non-DEHP)	
	21-7665	15/box	CADD <sup>®</sup> yellow-striped high volume administration set with NRFit <sup>™</sup> connector, bag spike, flow stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male Luer	136 in (345 cm)	23 mL	TOTM (non-DEHP)	21-7685

CADD <sup>®</sup> Extension Sets for CADD <sup>®</sup> -Solis, CADD-Legacy <sup>®</sup> and CADD-Prizm <sup>®</sup> Pumps - Designed to be used with a CADD <sup>™</sup> medication cassette reservoir							
	Reorder No.	Units	Description	Tubing Length	Priming Vol.	Plasticizer	Previous No.
	21-7045	50/box	CADD <sup>®</sup> extension set with male/male Luer, integral anti-siphon valve, and clamp	30 in (76 cm)	1.1 mL	TOTM (non-DEHP)	
	21-7060	12/box					
	21-7046	50/box	CADD <sup>®</sup> extension set with male/male Luer, integral anti-siphon valve, and clamp	45 in (114 cm)	1.4 mL	TOTM (non-DEHP)	
	21-7061	12/box					
	21-7047	50/box	CADD <sup>®</sup> extension set with male/male Luer, integral anti-siphon valve, and clamp	60 in (152 cm)	1.7 mL	TOTM (non-DEHP)	
	21-7062	12/box					
	21-7092	12/box	CADD <sup>®</sup> checkvalve extension set with male/male Luer, integral anti-siphon valve, y-extension with one-way checkvalve with female Luer, and clamp	68 in with 3 in pigtail (173 cm with 7.6 cm pigtail)	2.9 mL	TOTM (non-DEHP)	
	21-7115	12/box	CADD <sup>®</sup> checkvalve extension set with male/male Luer, 0.2-micron air-eliminating filter, integral anti-siphon valve with clamp, and y-extension with one-way checkvalve with female Luer	69 in with 3 in pigtail (175 cm with 7.6 cm pigtail)	3.2 mL	TOTM (non-DEHP)	
	21-7105	12/box	CADD <sup>®</sup> yellow-striped extension set with male/male Luer, integral anti-siphon valve, and clamp	90 in (229 cm)	2.3 mL	TOTM (non-DEHP)	
	21-7044	25/box	CADD <sup>®</sup> extension set with male/male Luer, 0.2-micron air-eliminating filter, and integral anti-siphon valve with clamp	55 in (140 cm)	4.6 mL	TOTM (non-DEHP)	21-7040
	21-7052	50/box	CADD <sup>®</sup> microbore extension set with male Luer, 0.2-micron air-eliminating filter, clamp, and integral anti-siphon valve with male Luer	60 in (152 cm)	1.1 mL	TOTM (non-DEHP)	
	21-7106	50/box	CADD <sup>®</sup> extension set with male/male Luer, 0.2-micron air-eliminating filter, integral anti-siphon valve, and clamp	60 in (152 cm)	1.8 mL	TOTM (non-DEHP)	
	21-7147	12/box	CADD <sup>®</sup> yellow-striped extension set with male/male Luer, 0.2-micron air-eliminating filter, and integral anti-siphon valve with clamp	100 in (254 cm)	5.1 mL	TOTM (non-DEHP)	21-7107
	21-7108	50/box	CADD <sup>®</sup> yellow-striped extension set with male/male Luer, 0.2-micron air-eliminating filter, and integral anti-siphon valve with clamp	90 in (229 cm)	2.4 mL	TOTM (non-DEHP)	
	21-7109	50/box	CADD <sup>®</sup> blue-striped extension set with male/male Luer, 0.2-micron air-eliminating filter, and integral anti-siphon valve with clamp	90 in (229 cm)	2.4 mL	TOTM (non-DEHP)	
	21-1200	12/box	Bag spike set with male Luer lock and clamp	15 in (38 cm)	0.5 mL	TOTM (non-DEHP)	
	21-7605	12/box	CADD <sup>®</sup> yellow-striped extension set with NRFit <sup>™</sup> connector, male/male Luer, integral anti-siphon valve, and clamp	90 in (229 cm)	2.3 mL	TOTM (non-DEHP)	21-7105
	21-7608	12/box	CADD <sup>®</sup> yellow-striped extension set with NRFit <sup>™</sup> connector, male/male Luer, 0.2-micron air-eliminating filter, and integral anti-siphon valve with clamp	90 in (229 cm)	2.4 mL	TOTM (non-DEHP)	21-7108

For more information visit our website at [www.smiths-medical.com](http://www.smiths-medical.com)

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN ALL GEOGRAPHIES

Smiths Medical ASD, Inc.  
6000 Nathan Lane North,  
Minneapolis, MN 55442, USA  
Tel: 1-614-210-7300,  
Toll-Free USA: 1-800-258-5361  
[www.smiths-medical.com](http://www.smiths-medical.com)

Rx  
ONLY  
smiths medical

Find your local contact information at: [www.smiths-medical.com/customer-support](http://www.smiths-medical.com/customer-support)

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. CADD, CADD-Legacy, CADD-Prizm, and the Smiths Medical design mark are trademarks of Smiths Medical. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. Product referenced is CE-Marked. ©2020 Smiths Medical. All rights reserved. IS0254.EN.AM.Rev.A 0820