

Characteristics

This popular black color masks a variety of stains and is a true medical grade exam glove. Engineered with a medium modulus and high tensile strength, this glove is first choice for a wide range of end users including law enforcement and corrections. Textured for excellent grip. Black color. Tested to provide maximum protection from fentanyl exposure.



**NitriDerm®
Ultra Black**

**Nitrile
Series 187**

**Exam Glove
Non-Sterile**



PRODUCT DETAILS

SIZE	ITEM NO.	PACKAGING	DESCRIPTION
XS	187050	100 Gloves/box, 10 boxes/case	Gloves, Exam, Nitrile, Non Sterile, Powder-Free, Chemo, Textured, Black
S	187100	100 Gloves/box, 10 boxes/case	
M	187200	100 Gloves/box, 10 boxes/case	
L	187300	100 Gloves/box, 10 boxes/case	
XL	187350	100 Gloves/box, 10 boxes/case	
XXL	187400	90 Gloves/box, 10 boxes/case	
XXXL	187450	90 Gloves/box, 10 boxes/case	

Product Attributes

- NFPA & UL Certified
- Low Modulus
- Non-Latex
- Textured Finish

Benefits

- Quality Assurance
- Softer, More Comfortable Fit
- No Risk of Latex Allergens
- Improved Wet/Dry Grip

Product Solutions You Trust



INNOVATIVE
HEALTHCARE
CORPORATION

www.ihsolutions.com

Specification Sheet

NitriDerm® Powder-Free **Ultra Black** Nitrile Exam Gloves



- 100 Qty (By Weight)
- Single Use
- Ambidextrous
- Non-Sterile

NitriDerm® Ultra Black is manufactured in compliance with multiple international standards, including the following:

Designation	Standard
ASTM D6319	Standard Specification for Nitrile Examination Gloves for Medical Application
ASTM D5151	Standard Test Method for Detection of Holes in Medical Gloves
ASTM F1671	Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens

Average Length	Average Palm Thickness	Average Finger Thickness
9.5 in ♦ 240 mm	4.0 mil ♦ 0.10 mm	5.5 mil ♦ 0.14 mm

Tensile Strength & Elongation	Before Aging	After Accelerated Aging
Tensile Strength (Mpa)	32	27
ASTM Requirement Min. (Mpa)	14	14
Elongation (%)	650	560
ASTM Requirement Min. (%)	500	400



Intertek

Innovative Healthcare Corporation
is certified to ISO 13485:2003 QMS
for medical devices.



INNOVATIVE
HEALTHCARE
CORPORATION

Product Solutions You Trust

1-800-272-1533