



Prep-Lock™

Tamper Evident Cap

With Male Luer Lock

Provides a secure tamper evident closure to devices that feature an ISO Compatible female Luer connector, such as CADD Cassettes, or a variety of sets. Prep-Lock™ Caps enhance medication safety and serve as an active deterrent to diversion. Tampering is immediately evident when the outer sleeve is broken or missing. They are the perfect compliment to your overall drug security program.

Improve <797> Compliance

Ensure the integrity of your compounds with USP recommended security.

Maintain Sterility

Evidence of access indicates the potential compromise of content's sterility.

Active Diversion Deterrent

Provides a physical barrier to the opportunity for misuse.



To place an order, please call or visit
integratedmedsys.com | 800.755.3800





Prep-Lock™ | Tamper Evident Cap With Male Luer Lock

Order Number	Description	Qty
57-91R	Male Luer Lock Red	Box of 100

Packaged Non-Sterile, Individual Blister, 100 Units Per Box

To place an order, please call or visit
integratedmedsys.com | 800.755.3800

