



Prep-Lock™

Tamper Evident Additive Port Cap

Compatible with B. Braun Excel® IV Bags

The Prep-Lock™ Tamper Evident Additive Port Cap provides remarkable protection and protocol assurance to the medication port of B. Braun Excel IV Bags. They provide visual evidence that medication has been added, tamper-evident protection and more. The simple one-handed installation and considerable contributions to the integrity of IV compounds have made these caps a favorite of pharmacists everywhere.

Diversion Prevention, Protocol Assurance & More

- Strengthens USP <797> Compliance.
- Provides protocol assurance with visual indication medication has been added.
- Protects against tampering or misuse.
- Simple one-handed installation is a superior solution over sticky foil seals.
- Advanced design guarantees once the cap is applied it cannot be removed without damage to the IV Bag.
- Reduces risk of contamination.
- Helps ensure product integrity.
- Helps reduce waste and shortages, by providing confidence that unadministered bags can be returned to the pharmacy.



MADE IN THE USA

IMIWEB.COM | 1.800.344.2554



Prep-Lock™ | Tamper Evident Additive Port Caps Compatible with B. Braun Excel

Order #	Description	Qty
21-02	B. Braun Excel® ·Light Blue	Box of 300
21-02-100	B. Braun Excel® ·Light Blue	Box of 100

Packaged Non-Sterile, 50 Units Per Bag

Technical Factors	
Sterility:	·Packaged non sterile.
DEHP:	·This product does not contain DEHP.
BPA:	·This product does not contain BPA.
Latex:	·This product is not made with natural rubber latex.
Product Material:	·Made with polypropylene.

To place an order, please call or visit

IMIWEB.COM | 1.800.344.2554

