




INNOVATION IN SECURE DRUG DELIVERY

PREP-LOCK™ | PREP-SEAL™ | PREP-FILL™ | RX-VENT™ | RX-TRACT™




At IMI we take great pride in every step of our design and manufacturing process. Quality is at the forefront in everything we do. IMI is certified to meet the requirements of ISO13485. All IMI products are made in the USA at our FDA registered facility under strict quality control.

To Place an Order, Please Call Us at
(800) 344-2554 or Visit www.imiweb.com

MADE IN THE
USA



LEADING THE WAY
IN THE DESIGN AND MANUFACTURE
OF INNOVATIVE MEDICAL DEVICES



TAMPER EVIDENT CAPS FOR LUER LOCK IV SYRINGES

Prep-Lock™ Tamper Evident Caps safeguard the integrity of your medications, by providing a secure, sterile, tamper-evident closure to IV Syringes. Prep-Lock™ Caps enhance medication safety, reduce waste and associated cost, and serve as an active deterrent to diversion.

The Tamper Evident Cap for Luer Lock IV Syringes is also available individually packaged. Perfect for lower volume compounding requirements.

- Help maintain security and sterility, increasing the confidence level of HCPs.

ORDER #	COLOR	QTY
57-24-100	RED	BOX OF 100
57-25-100	BLUE	BOX OF 100
57-26-100	WHITE	BOX OF 100
57-24	RED	BOX OF 1000
57-25	BLUE	BOX OF 1000
57-26	WHITE	BOX OF 1000

PACKAGED STERILE
10 UNITS PER TRAY, SOLID OUTER SLEEVE

ORDER #	COLOR	QTY
57-80	RED	BOX OF 100

PACKAGED STERILE
INDIVIDUAL BLISTER, SOLID OUTER CAP



TAMPER EVIDENT CAPS FOR LUER LOCK IV SYRINGES

Compounded sterile preparations (CSPs) are at their greatest risk when they leave the custody of your pharmacy. The Prep-Lock Tamper Evident Cap greatly reduces the risk of your CSPs being compromised either accidentally or intentionally. It ensures that product integrity is maintained from the time it leaves the sterile hood until it is administered by an authorized clinician.

- Tampering is immediately evident when the outer sleeve is broken or missing.
- Syringes with intact caps can be returned to the pharmacy.

ORDER #	COLOR	QTY
57-50R	RED	BOX OF 1000
57-50B	BLUE	BOX OF 1000
57-50W	WHITE	BOX OF 1000

PACKAGED STERILE
10 UNITS PER TRAY, CLEAR OUTER SLEEVE



TAMPER EVIDENT CAP WITH MALE LUER LOCK

Provides a secure tamper evident closure to devices that feature an ISO Compatible female Luer connector, such as CADD Cassettes, or a variety of sets.

- An active deterrent to potential diversion and misuse.
- Clearly indicates evidence of access, whether malicious or unintentional.
- Maintain Product Sterility. Evidence of access indicates potential compromise of container contents.

ORDER #	COLOR	QTY
57-91R	RED	BOX OF 100

PACKAGED STERILE
INDIVIDUAL BLISTER, SOLID OUTER CAP



TAMPER EVIDENT CAPS FOR ORAL SYRINGES

These Tamper Evident Caps are designed specifically for Oral Dispensers and utilize IMI's patented "Break Away" sleeve technology.

- Makes tampering immediately evident when the colored outer sleeve is broken or missing.
- Featuring retained ring for improved control of disposable components.

ORDER #	TYPE/COLOR	QTY
57-70	BAXTER - RED TOP	BOX OF 1000
57-70-100	BAXTER - RED TOP	BOX OF 100
57-73	BD - YELLOW TOP	BOX OF 1000
57-73-100	BD - YELLOW TOP	BOX OF 100
57-74	COMAR® - WHITE TOP (AVAILABLE 2ND QUARTER 2021)	BOX OF 1000
57-75	NEOMED/AVANOS - ORANGE TOP (AVAILABLE 2ND QUARTER 2021)	BOX OF 1000

PACKAGED NON-STERILE, 100 UNITS PER BAG

*** IMPORTANT:**

Oral Tamper Evident Caps are not a universal fit. You must match your Oral Tamper Evident Cap to the dispenser manufacturer.



TAMPER EVIDENT CAP FOR ENFIT® SYRINGES

Enteral administration is nutrition or drug administered via the GI tract, such as elixirs and suspensions. Prep-Lock Tamper Evident Caps for ENFit® Syringes are designed to provide secure tamper evident closure to the ENFit connector standard (ISO 80369-3) and are the only ENFit Syringe Cap available with built in Tamper Evident Technology.

Compatible with low dose and standard ENFit Syringes, they provide an additional level of security to syringe contents during transportation prior to administration, optimizing the care of neonatal, pediatric and any patients requiring enteral medications.

ORDER #	COLOR	QTY
57-220	PURPLE	BOX OF 1000
57-220-100	PURPLE	BOX OF 100

PACKAGED NON-STERILE, 100 UNITS PER BAG

ENFit® is a GEDSA trademark.

- Cannot be reinstalled after removal.
- Compatible with both standard and low dose ENFit syringes.
- Clearly indicates potential compromise of syringe contents.
- Increases confidence in your enteral medications, by ensuring content's integrity.



TAMPER EVIDENT ADDITIVE PORT CAPS FOR SELECT BAXTER & B.BRAUN IV BAGS

The Prep-Lock™ Tamper Evident Additive Port Cap provides remarkable protection and protocol assurance to the medication port of select Baxter and B.Braun IV Bags. They provide both visual evidence that medication has been added and tamper evident protection. The simple one-handed installation has made these caps a favorite of compounding pharmacists.

The enclosed hinge and hidden lock design ensures that once the Cap is applied to the medication port, it cannot be removed without damage to the IV bag.


ORDER #	COLOR/TYPE	QTY
21-00	RED Baxter Viaflex, Intravia, Aviva, and All-In-One	BOX OF 300
21-00-100	RED Baxter Viaflex, Intravia, Aviva, and All-In-One	BOX OF 100

PACKAGED NON-STERILE, 100 UNITS PER BAG

ORDER #	COLOR/TYPE	QTY
21-01	GREEN B. Braun PAB®/B. Braun E3®	BOX OF 300
21-01-100	GREEN B. Braun PAB®/B. Braun E3®	BOX OF 100

PACKAGED NON-STERILE, 50 UNITS PER BAG

- Simple one-handed installation is a superior solution over sticky foil seals.
- Protects against potential tampering or misuse.



GUARDED LUER CONNECTORS FOR LUER LOCK & LUER SLIP

Prep-Fill™ Adapters and Connectors provide convenient, economical, high-quality solutions to facilitate transferring or mixing of fluids between drug delivery devices or IV sets.

- Expanded wing size means larger finger grips for secure manipulation, helping to reduce touch contamination
- Patented notched wing design enhances first air over the connection port.
- Built-in drip tray design to capture stray fluid droplets during repeater pump procedures.

ORDER #	TYPE/COLOR	QTY
57-400	FEMALE LUER LOCK / FEMALE LUER LOCK - RED	BOX OF 50
57-400W	FEMALE LUER LOCK / FEMALE LUER LOCK - WHITE	BOX OF 50
57-401	FEMALE LUER LOCK / MALE LUER LOCK - YELLOW	BOX OF 50
57-402	FEMALE LUER SLIP / FEMALE LUER LOCK - BLUE	BOX OF 50
57-403	FEMALE LUER SLIP / FEMALE LUER SLIP - VIOLET	BOX OF 50

PACKAGED STERILE, NON-PYROGENIC
INDIVIDUAL BLISTER




LUER CONNECTORS & ADAPTERS

FOR LUER LOCK & LUER SLIP

Prep-Fill™ Adapters and Connectors provide convenient, economical, high-quality solutions to facilitate transferring or mixing of fluids between drug delivery devices or IV sets.

57-02	MALE - MALE LUER SLIP	BOX OF 100
57-13	FEMALE - FEMALE LUER LOCK	BOX OF 100

PACKAGED STERILE, NON-PYROGENIC
INDIVIDUAL BLISTER



CAPS & PLUGS

FOR LUER LOCK, LUER SLIP & NRFIT®

The Prep-Seal™ line of products offers high-quality, sterile closure solutions for a variety of connector types, whether Luer Lock, Luer Slip, or NRFit® IMI has your closure devices covered.

To further increase patient safety and reduce the risk of misconnections a new international standard referred to as NRFit (ISO 80369-6) has been developed. The new standard specifies the requirements for small bore connectors used in neuraxial applications such as spinal and epidural procedures. IMI Prep-Seal™ NRFit Caps provided sterile closure to both male and female NRFit devices.

ORDER #	TYPE/COLOR	QTY
57-11	FEMALE LUER LOCK CAP	BOX OF 100
57-20	MALE LUER SLIP PLUG	BOX OF 100
57-602	FEMALE NRFit® CAP	BOX OF 100
57-601	MALE NRFit® CAP	BOX OF 100

PACKAGED STERILE, NON-PYROGENIC
INDIVIDUAL BLISTER

PREP-SEAL™ | LUER LOCK CAPS




CAPS

FOR LUER LOCK SYRINGES

The Prep-Seal™ Luer Lock Syringe Cap comes in a convenient sterile tray of 10 units. If your compounding needs call for volume capping of syringes, the Prep-Seal Luer Lock Cap is the perfect solution. A unique “keyed” tray bottom lets you install caps quickly and efficiently with a simple twist of the syringe; further enhancing aseptic technique. Offered in 6 unique color options to choose from.

ORDER #	COLOR	QTY
57100G	GREEN	BOX OF 100
57100R	RED	BOX OF 100
57100B	BLUE	BOX OF 100
57100Y	YELLOW	BOX OF 100
57100V	VIOLET	BOX OF 100
57100W	WHITE	BOX OF 100

PACKAGED STERILE, NON-PYROGENIC
10 UNITS PER TRAY



ASPIRATING NEEDLES

SHARP & BLUNT CANNULA

Rx-Tract™ Needles provide high-quality solutions for introducing fluids to deep wells and vials, drawing liquids for reconstruction or pre-fill syringes prior to administration.

- Includes a convenient needle guard.

ORDER #	TYPE	QTY
32-21	19 GAUGE X 3" (7.62 CM) SHARP	BOX OF 100
32-22	19 GAUGE X 5" (12.7 CM) SHARP	BOX OF 100
32-23	19 GAUGE X 8" (20.3 CM) SHARP SPECIAL ORDER ONLY	BOX OF 100
32-34	16 GAUGE X 3.5" (8.8 CM) BLUNT	BOX OF 100
32-31	16 GAUGE X 5" (12.7 CM) BLUNT	BOX OF 100
32-32	16 GAUGE X 8" (20.3 CM) BLUNT	BOX OF 100

STAINLESS STEEL NEEDLE, PACKAGED STERILE
NON-PYROGENIC, INDIVIDUAL BLISTER



VENTING NEEDLES CHEMO-VENT™ / FILTERED

Filtered Venting Needles are designed to vent evacuated containers and to maintain sterility. Less expensive than other venting needles, the non-coring needle easily penetrates all stoppers with minimum risk of damage and hydrophobic filter helps maintain asepsis.

The Chemo-Vent™ filtered venting needle is designed for use when reconstituting and withdrawing cytotoxic anti-neoplastic agents. The non-coring needle easily penetrates all stoppers with minimum risk of damage, while the hydrophobic filter maintains asepsis and prevents toxic, aerosolized materials from leaving the container during filling.

ORDER #	TYPE	QTY
31-50	CHEMO-VENT™ 20 GAUGE X 1.5" (3.8 CM) NON-CORING VENTING NEEDLE WITH .2 MICRON HYDROPHOBIC FILTER	BOX OF 100
31-20	FILTERED 20 GAUGE X 1.5" (3.8 CM) NON-CORING VENTING NEEDLE WITH .2 MICRON HYDROPHOBIC FILTER	BOX OF 100

STAINLESS STEEL NEEDLE, PACKAGED STERILE
INDIVIDUAL BLISTER



INTERNATIONAL MEDICAL INDUSTRIES, INC.
2981 GATEWAY DRIVE
POMPANO BEACH, FLORIDA USA 33069
1 (954) 917-9570

TO PLACE AN ORDER, PLEASE CALL US AT
(800) 344-2554 OR VISIT WWW.IMIWEB.COM



© 2021 International Medical Industries, Inc. All rights reserved.
Companies and products mentioned herein may be trademarks or registered trademarks of their respective trademark owners.