

CLEO 90

INFUSION SET

THE ADVANTAGES ARE CLEAR



All-in-One Design

Includes an introducer and cannula to make insertions simple and allow insertion speed control.

Enclosed Needle Safety

The enclosed needle retracts automatically into the housing after insertion, reducing the risk of accidental needlestick injuries.

Hypoallergenic Dressing

The small, transparent adhesive dressing is designed to breathe and stay in place with less tape. Its release tab makes removal easy too.

Flexible 360° Connection

The Flex AttachTM 360° site tubing connection provides 360 on-site rotation to allow for greater flexibility.

Ordering Details

Item #	Description	UOM (box)
IM53080	Cleo SubQ Infusion Cannula, 6mm	50 Eaches
IM53081	Cleo SubQ Infusion Cannula, 9mm	50 Eaches
IM53082	Cleo SubQ Infusion Set Tubing, 24"	50 Eaches
IM53083	Cleo SubQ Infusion Set Tubing, 31"	50 Eaches
IM53084	Cleo SubQ Infusion Set Tubing, 42"	50 Eaches

**Fast. Simple.
Ready to Order?**

Scan the QR code or visit integratedmedsys.com to buy now or learn more.



Proudly
Distributed by

IMS
INTEGRATED MEDICAL SYSTEMS



Transforming Healthcare, Together. Beyond Distribution.

Solutions. Technology. Service.

IMed[™]
Products

Leading the healthcare industry towards a more efficient and effective future through our commitment to excellence in developing cutting-edge medical products and solutions.

The product details provided herein, including any images or supplements, reflect information from various sources regarding the product; this information has not been independently verified by Integrated Medical Systems ("IMS") or its affiliates. IMS is not liable for any inaccuracies or omissions in the product information, which is provided solely for your convenience.

The information presented here is subject to change without prior notice. Customers should review product packaging and labeling before use or purchase to ensure accuracy, especially when referencing information from IMS's online or printed materials. IMS disclaims any responsibility for the accuracy or completeness of the information reflected herein, regardless of whether any inaccuracies or omissions are due to errors or other factors.