

# EMPTY EVA BAGS

IMed Products™ Empty EVA Bags are used in IV infusion systems to deliver fluids and medications. Made from durable, biocompatible ethylene-vinyl acetate, they are flexible, lightweight, and leak-resistant for safe treatment. Featuring injection sites and patient access spike ports, they enable controlled medication delivery.



### Features & Benefits

- Not Made with Natural Rubber Latex
- DEHP-Free
- Pyrogen-Free
- Non-Toxic
- Made with Ethyl Vinyl Acetate (EVA)
- PVC Ports & Tubing
- Durable & Transparent
- Vented Spikes\*
- Leg Closing Cap\*
- Detachable Transfer Sets\*
- 2-Port, 2-Leg, and 3-Leg Options


\*2-Leg and 3-Leg options only.

### Ordering Details

Item #	Size	Legs	Detachable Transfer Set	Vented Spikes	Leg Closing Cap	Access Port	UOM (case)
<b>2-Port Bags</b>							
IM68040	150mL	-	-	-	-	Y	50 Eaches
IM68045	250mL	-	-	-	-	Y	50 Eaches
IM68050	500mL	-	-	-	-	Y	50 Eaches
IM68055	1000mL	-	-	-	-	Y	50 Eaches
<b>2-Leg Bags</b>							
IM68102	500mL	2-Leg	Y	Y	Y	Y	50 Eaches
IM68100	1000mL	2-Leg	Y	Y	Y	Y	50 Eaches
IM68105	2000mL	2-Leg	Y	Y	Y	Y	50 Eaches
IM68110	3000mL	2-Leg	Y	Y	Y	Y	50 Eaches
<b>3-Leg Bags</b>							
IM68120	1000mL	3-Leg	Y	Y	Y	Y	50 Eaches
IM68125	2000mL	3-Leg	Y	Y	Y	Y	50 Eaches
IM68130	3000mL	3-Leg	Y	Y	Y	Y	50 Eaches
IM68135	4000mL	3-Leg	Y	Y	Y	Y	50 Eaches

**Fast. Simple.  
Ready to Order?**

Scan the QR code or visit [integratedmedsys.com](http://integratedmedsys.com) to buy now or learn more.



Proudly  
Distributed by

**IMS**  
INTEGRATED MEDICAL SYSTEMS



**Transforming Healthcare, Together. Beyond Distribution.**

**Solutions. Technology. Service.**

**IMed**<sup>™</sup>  
Products

Leading the healthcare industry towards a more efficient and effective future through our commitment to excellence in developing cutting-edge medical products and solutions.

The product details provided herein, including any images or supplements, reflect information from various sources regarding the product; this information has not been independently verified by Integrated Medical Systems ("IMS") or its affiliates. IMS is not liable for any inaccuracies or omissions in the product information, which is provided solely for your convenience.

The information presented here is subject to change without prior notice. Customers should review product packaging and labeling before use or purchase to ensure accuracy, especially when referencing information from IMS's online or printed materials. IMS disclaims any responsibility for the accuracy or completeness of the information reflected herein, regardless of whether any inaccuracies or omissions are due to errors or other factors.