

IMed Products® Empty EVA Bags are durable, biocompatible IV infusion bags used to deliver fluids and medications. Made from flexible, lightweight, leak-resistant EVA, they feature injection sites and patient access spike ports for safe, controlled delivery.

## FEATURES & BENEFITS

- Not Made with Natural Rubber Latex
- DEHP-Free
- Pyrogen-Free
- Non-Toxic
- Made with Ethyl Vinyl Acetate (EVA)
- PVC Ports & Tubing
- Durable & Transparent
- Vented Spikes\*
- Leg Closing Cap\*
- Detachable Transfer Sets\*
- 2-Port, 2-Leg, and 3-Leg Options

\* 2-Leg and 3-Leg options only.



## ORDERING DETAILS

ITEM #	SIZE	LEGS	DETACHABLE TRANSFER SET	VENTED SPIKES	LEG CLOSING CAP	ACCESS PORT	UOM (CASE)
<b>2-PORT BAGS</b>							
<b>IM68040</b>	100mL	-	-	-	-	Y	50 Eaches
<b>IM68045</b>	250mL	-	-	-	-	Y	50 Eaches
<b>IM68050</b>	500mL	-	-	-	-	Y	50 Eaches
<b>IM68055</b>	1000mL	-	-	-	-	Y	50 Eaches
<b>2-LEG BAGS</b>							
<b>IM68102</b>	500mL	2-Leg	Y	Y	Y	Y	50 Eaches
<b>IM68100</b>	1000mL	2-Leg	Y	Y	Y	Y	50 Eaches
<b>IM68105</b>	2000mL	2-Leg	Y	Y	Y	Y	50 Eaches
<b>IM68110</b>	3000mL	2-Leg	Y	Y	Y	Y	50 Eaches
<b>3-LEG BAGS</b>							
<b>IM68120</b>	1000mL	3-Leg	Y	Y	Y	Y	50 Eaches
<b>IM68125</b>	2000mL	3-Leg	Y	Y	Y	Y	50 Eaches
<b>IM68130</b>	3000mL	3-Leg	Y	Y	Y	Y	50 Eaches
<b>IM68135</b>	4000mL	3-Leg	Y	Y	Y	Y	50 Eaches

# Strengthened by Nationwide Distribution and Dedicated Customer Support

Proudly  
Distributed by

**IMS**  
INTEGRATED MEDICAL SYSTEMS



**IMed**<sup>®</sup>  
Products

**IMed Products delivers reliable, high-quality medical support solutions at value-driven pricing—designed to ensure consistent supply, trusted clinical performance, and operational ease for outpatient centers, clinics, hospitals, and pharmacies.**

The product details provided herein, including any images or supplements, reflect information from various sources regarding the product; this information has not been independently verified by Integrated Medical Systems ("IMS") or its affiliates. IMS is not liable for any inaccuracies or omissions in the product information, which is provided solely for your convenience.

The information presented here is subject to change without prior notice. Customers should review product packaging and labeling before use or purchase to ensure accuracy, especially when referencing information from IMS's online or printed materials. IMS disclaims any responsibility for the accuracy or completeness of the information reflected herein, regardless of whether any inaccuracies or omissions are due to errors or other factors.

© 2026 IMed® Products. All rights reserved. 02/26 IMS25090 IMed Sell Sheets\_Brand Refresh EVA Bags\_01