

IV ADMINISTRATION SETS

IMed Products™ IV Administration Sets are designed to connect an IV solution bag to a patient’s access site. These products are certified sterile using ethylene oxide. The sets include a pinch or roller clamp, with options of a flow regulator for precise flow rate control.



Features & Benefits

- Not Made with Natural Rubber Latex
- DEHP-Free
- Pyrogen-Free
- Non-Toxic
- Sterilized by Ethylene Oxide
- Drip Chamber with 15u Filter
- Universal Vented/Non-Vented Spike
- Distal Y-Injection Site (needle-free option available)
- Offered in a Variety of Lengths, Priming Volumes, and Connection Types

Ordering Details

| Item # | Length | Priming Volume | Drops Per mL | Y-Injection Site | Flow Regulator | 0.22u In-Line Filter | Clamp Type | Connection Type | UOM (case) |
|---------------------------------------|--------|----------------|--------------|------------------|----------------|----------------------|--------------|--------------------|------------|
| IV ADMINISTRATION GRAVITY SETS | | | | | | | | | |
| IM42075* | 42" | 7mL | 20 | - | - | - | Roller Clamp | Rotating Male Luer | 50 Eaches |
| IM42070 | 96" | 19mL | 20 | - | - | Y | Roller Clamp | Rotating Male Luer | 50 Eaches |
| IM42071 | 96" | 16mL | 10 | Y | - | - | Roller Clamp | Rotating Male Luer | 50 Eaches |
| IM42072 | 96" | 16mL | 15 | Y | - | - | Roller Clamp | Rotating Male Luer | 50 Eaches |
| IM42073 | 96" | 16mL | 20 | Y | - | - | Roller Clamp | Rotating Male Luer | 50 Eaches |
| WITH FLOW REGULATOR | | | | | | | | | |
| IM42031 | 92" | 15mL | 20 | Y | Y | - | Pinch Clamp | Rotating Male Luer | 50 Eaches |
| IM42074 | 92" | 15mL | 20 | Y | Y | - | Pinch Clamp | Fixed Male Luer | 50 Eaches |
| IM42076** | 92" | 15mL | 20 | Y | Y | - | Pinch Clamp | Rotating Male Luer | 50 Eaches |
| IM42069 | 92" | 18mL | 20 | Y | Y | Y | Pinch Clamp | Rotating Male Luer | 50 Eaches |

*includes a bag hanger.

**includes a needle-free y-site.

**Fast. Simple.
Ready to Order?**

Scan the QR code or visit integratedmedsys.com to buy now or learn more.



Proudly
Distributed by

IMS
INTEGRATED MEDICAL SYSTEMS



Transforming Healthcare, Together. Beyond Distribution.

Solutions. Technology. Service.

IMed[™]
Products

Leading the healthcare industry towards a more efficient and effective future through our commitment to excellence in developing cutting-edge medical products and solutions.

The product details provided herein, including any images or supplements, reflect information from various sources regarding the product; this information has not been independently verified by Integrated Medical Systems ("IMS") or its affiliates. IMS is not liable for any inaccuracies or omissions in the product information, which is provided solely for your convenience.

The information presented here is subject to change without prior notice. Customers should review product packaging and labeling before use or purchase to ensure accuracy, especially when referencing information from IMS's online or printed materials. IMS disclaims any responsibility for the accuracy or completeness of the information reflected herein, regardless of whether any inaccuracies or omissions are due to errors or other factors.