

IMed Products® IV Administration Extension Sets are designed to securely extend the connection between an IV fluid source and a patient's access site, offering flexible, sterile solutions for fluid delivery. These sets feature male/female luer connections, pinch clamps for flow control, and optional flow regulator for precise rate adjustment.

FEATURES & BENEFITS

- Not Made with Natural Rubber Latex
- DEHP-Free
- Pyrogen-Free
- Non-Toxic
- Sterilized by Ethylene Oxide
- Pinch Clamp
- Fixed Female Luer Connection
- Fixed or Rotating Male Luer Connection
- Offered in a Variety of Lengths and Attributes
 - In-Line Filters
 - Flow Regulator
 - Y-Injection Site



ORDERING DETAILS

ITEM #	LENGTH	PRIMING VOLUME	IN-LINE FILTER	FLOW REGULATOR	Y-INJECTION SITE	CONNECTION TYPE	UOM (CASE)
IM471000	7"	0.27mL	-	-	-	Fixed Male Luer	50 Eaches
IM471001	12"	0.40mL	-	-	-	Fixed Male Luer	50 Eaches
IM471002	15"	0.48mL	-	-	-	Fixed Male Luer	50 Eaches
IM471067	12"	4.8mL	0.2 Micron	-	-	Rotating Male Luer	50 Eaches
IM471068	12"	4.8mL	0.2 Micron	-	-	Rotating Male Luer	50 Eaches
IM471018	19"	3.5mL	-	Y	Y	Rotating Male Luer	50 Eaches

Strengthened by Nationwide Distribution and Dedicated Customer Support

Proudly
Distributed by

IMS
INTEGRATED MEDICAL SYSTEMS



IMed[®]
Products

IMed Products delivers reliable, high-quality medical support solutions at value-driven pricing—designed to ensure consistent supply, trusted clinical performance, and operational ease for outpatient centers, clinics, hospitals, and pharmacies.

The product details provided herein, including any images or supplements, reflect information from various sources regarding the product; this information has not been independently verified by Integrated Medical Systems ("IMS") or its affiliates. IMS is not liable for any inaccuracies or omissions in the product information, which is provided solely for your convenience.

The information presented here is subject to change without prior notice. Customers should review product packaging and labeling before use or purchase to ensure accuracy, especially when referencing information from IMS's online or printed materials. IMS disclaims any responsibility for the accuracy or completeness of the information reflected herein, regardless of whether any inaccuracies or omissions are due to errors or other factors.

© 2026 IMed® Products. All rights reserved. 02/26 IMS25044 IMed Sell Sheets_Brand Refresh IV Extension Sets_01