





Product Catalog









Integrated Medical Systems' (IMS) private label line, IMed Products™, offers a broad range of high quality and cost-effective medical disposables spanning across many product categories. Strategically cultivated to help offset the continuation of rising costs, IMed Products promotes a standardized approach to products and pricing. We are constantly adding and updating our product line to best align with our customers' evolving wants and needs.

Follow Us for the Latest News & Updates





Scan to View the Entire Product Line:



Table of Contents

Infusion Offerings

Backpacks & Pouches		Kits	
Backpacks	05	IV Start Kits	26
Pouches	06	Dressing Change Kits	30
Elastomeric Pouches	07	Miscellaneous Kits	36
Dressings & Bandages		Custom Kits	39
IMed Guard	08		
Elastic Nets	09	Prep Pads & Wipes	
Cohesive Bandages	10	Alcohol Prep Pads	40
Elastomeric Pumps		Subcutaneous Supplies	
Short Duration	11	Cleo® 90 Infusion Set	41
Long Duration	13	Subcutaneous Dose Kit	42
Chemotherapy	15	Swinger	
End Caps		Syringes Luer Lock Syringes	43
Sterile End Caps	17	· -	13
EVA Bags		Enteral Offerings	
2 Port	18	Backpacks & Pouches	
Gravity Bags	19	Backpacks	45
Extension Sets		Syringes	
Smallbore	20	Enfit Syringes	46
With Micron Filter	21		
With Flow Regulator	22	Cleanroom Offerings	
IV Admin Sets		Disinfectants	
Gravity	23	Sterile Isopropyl Alcohol (IPA)	48
With Flow Regulator	24		
Secondary	25	Respiratory Offerings	
		HME's and Filters	
		Filters	50
		HME's	51



Infusion



IMed Products[™] durable backpacks are economical and stylish for patients receiving outpatient infusion therapies via ambulatory infusion pumps. These backpacks are universal size, for use with CADD, Curlin, and Sapphire pumps; ideal for patients on the go.





Item #	Description	Size	UOM
IM58021	IMed Ambulatory Pump Backpack	2000mL	1/Each
IM58008U	IMed Ambulatory Pump Backpack	3000mL	1/Each

IMed Products[™] durable pouches are economical and stylish for patients receiving outpatient infusion therapies via ambulatory infusion pumps; ideal for patients on the go.










Item #	Description	Size	UOM
IM58005U	IMed Ambulatory Pump Pouch	100mL	1/Each
IM58006U	IMed Ambulatory Pump Pouch	500mL	1/Each
IM58007U	IMed Ambulatory Pump Pouch	1000mL	1/Each
IM58019	IMed Fanny Pump Pouch	500mL	1/Each


IMed Products[™] durable pouches are economical and stylish for patients receiving outpatient infusion therapies via elastomeric pumps; ideal for patients on the go. Perfect for use with EZ-FLOW Elastomeric Pumps!



Item #	Description	Size	UOM
IM58010	IMed Elastomeric Pump Pouch	Holds up to 200mL	1/Each
IM58011	IMed Elastomeric Pump Pouch	Holds up to 500mL	1/Each


The IMed Guard is a sterile, single-use disposable disc infused with the antimicrobial agent chlorhexidine gluconate (CHG) to reduce contamination of the dressing by inhibiting microbial growth and colonization of the dressing.

- Beveled edges to reduce potential for catheter kinking and securement
- Soft material to reduce skin irritation/abrasion
- 38% higher absorbency of foam than market-leading product
- Pyrogen and latex-free
- IMed Guard has undergone applicable in-vitro, biocompatibility, sterility, and performance testing for safety and effectiveness



Item #	Description	UOM
IM36026	IMed Guard, 1" (4mm) Disc with CHG	10/Box

IMed Products™ Elastic Net Retainer Dressings are available in a variety of sizes to provide securement for IV lines/dressing for adults and pediatric patients. These dressing retainers offer compression and support to hold tubing firmly in place without the need for tape or gauze. They eliminate skin stripping and excoriation caused by tape removal and are very effective when used for burn care, surgery, post-op dressing changes, donor wound graft sites, daily dressing changes, EKG stress testing and sleep studies.



Item #	Description	Size	UOM
IMS03	IMed Elastic Net Retainer Dressing	Adult: small hand or foot Child: arm	1 Roll/Box
IMS04	IMed Elastic Net Retainer Dressing	Adult: medium hand or foot Child: medium arm	1 Roll/Box
IMS05	IMed Elastic Net Retainer Dressing	Adult: small arm or foot Child: XL hand, large leg	1 Roll/Box
IMS06	IMed Elastic Net Retainer Dressing	Adult: medium arm, XL hand or foot, small leg	1 Roll/Box
IMS07	IMed Elastic Net Retainer Dressing	Adult: large arm, small thigh Child: small chest or back	1 Roll/Box
IMS08	IMed Elastic Net Retainer Dressing	Adult: XXL arm, leg or shoulder, small chest	1 Roll/Box
IMS09	IMed Elastic Net Retainer Dressing	Adult: head, abdomen, chest, back, XXL leg	1 Roll/Box

IMed Products[™] Cohesive Bandages are an excellent securement device that adapts easily to body contours.

- Self-adherent
- Easily tear with hands
- Secures IVs and other devices
- Slip-resistant



Item #	Description	Size	UOM
IM36085	Tan Cohesive Bandages	2" x 5 Yards	36/Box
IM36086	Assorted Color Cohesive Bandages	2" x 5 Yards	36/Box
IM36087	Tan Cohesive Bandages	3" x 5 Yards	24/Box
IM36088	Assorted Color Cohesive Bandages	3" x 5 Yards	24/Box



EZ-FLOW™ Elastomeric Pumps are a convenient alternative to electronic infusion pumps. The EZ-FLOW™ offers a safe and effective way to deliver medications while allowing patients freedom in their day-to-day lives without the worry of programming and battery power.

- Latex-free
- DEHP-free
- Single-use only
- Filter to eliminate particles (1.2 microns) and air (0.03 microns)
- Available in 50mL/hour 500mL/hour flow rates
- Color-coded in multiple locations for easily identifiable flow rate






Item #	Fill Volume/Flow Rate	UOM
IM050050S	50ML/50ML Hr	24/CASE
IM100050S	100ML/50ML Hr	24/CASE
IM100100S	100ML/100ML Hr	24/CASE
IM100200S	100ML/200ML Hr	24/CASE
IM200100S	200ML/100ML Hr	24/CASE
IM200200S	200ML/200ML Hr	24/CASE
IM250050S	250ML/50ML Hr	24/CASE
IM250100S	250ML/100ML Hr	24/CASE
IM250125S	250ML/125ML Hr	24/CASE
IM250175S	250ML/175ML Hr	24/CASE
IM250250S	250ML/250ML Hr	24/CASE
IM250500S	250ML/500ML Hr	24/CASE
IM270175S	270ML/175ML Hr	24/CASE
IM400100S	400ML/100ML Hr	24/CASE
IM400200S	400ML/200ML Hr	24/CASE
IM500250S	500ML/250ML Hr	24/Case



EZ-FLOW™ Elastomeric Pumps are a convenient alternative to electronic infusion pumps. The EZ-FLOW™ offers a safe and effective way to deliver medications while allowing patients freedom in their day-to-day lives without the worry of programming and battery power.

- Latex-free
- DEHP-free
- Single-use only
- Filter to eliminate particles (1.2 microns) and air (0.03 microns)
- Available in 2mL/hour 20mL/hour flow rates
- Color-coded in multiple locations for easily identifiable flow rate
- Priming cap on long duration and chemotherapy pumps eliminates the use of the pinch clamp and removal of patient end cap when priming





Item #	Fill Volume/Flow Rate	UOM
IM060005L	60ML/5ML Hr	24/CASE
IM100002L	100ML/2ML Hr	24/CASE
IM100005L	100ML/5ML Hr	24/CASE
IM125005L	125ML/5ML Hr	24/CASE
IM270010L	270ML/10ML Hr	24/CASE
IM270001L	270ML/1ML Hr	24/CASE
IM270002L	270ML/2ML Hr	24/CASE
IM270004L	270ML/4ML Hr	24/CASE
IM270005L	270ML/5ML Hr	24/CASE
IM500020L	500ML/20ML Hr	24/CASE



EZ-FLOW™ Elastomeric Pumps are a convenient alternative to electronic infusion pumps. The EZ-FLOW™ offers a safe and effective way to deliver medications while allowing patients freedom in their day-to-day lives without the worry of programming and battery power.

- Latex-free
- DEHP-free
- Single-use only
- Filter to eliminate particles (1.2 microns) and air (0.03 microns)
- Amber colored outer covering for chemotherapy pumps
- Available in 2mL/hour 10mL/hour flow rates
- Color-coded in multiple locations for easily identifiable flow rate
- Priming cap on long duration and chemotherapy pumps eliminates the use of the pinch clamp and removal of patient end cap when priming





Item #	Fill Volume/Flow Rate	UOM
IM100002C	100ML/2ML Hr	12/CASE
IM125005C	125ML/5ML Hr	12/CASE
IM270002C	270ML/2ML Hr	12/CASE
IM270005C	270ML/5ML Hr	12/CASE
IM270010C	270ML/10ML Hr	12/CASE
IM300006C	300ML/6ML Hr	12/CASE

IMed Products[™] Dual-Function Dead End Caps feature a male luer on one end and a female luer on the other, for safe and secure connection.

- Latex-free
- DEHP-free
- PVC-free
- BPA-free
- Meets ISO 80369 standard
- Sterile-packed



Item #	Description	UOM
IM51003	IMed Sterile End Cap, Red	100/Case

IMed Products™ Ethyl Vinyl Acetate (EVA) 2-port bags are made with an EVA material that does not contain harmful plasticizers such as DEHP, which can leach into TPN solutions. EVA is more durable for reduced breakage, more stable for improved safety and more transparent for improved admixture inspection.

- Latex-free
- DEHP-free
- PVC ports



Item #	Description	Size	UOM
IM68040	IMed EVA Bag, 2-Ports	150mL	50/Case
IM68045	IMed EVA Bag, 2-Ports	250mL	50/Case
IM68050	IMed EVA Bag, 2-Ports	500mL	50/Case
IM68055	IMed EVA Bag, 2-Ports	1000mL	50/Case

IMed Products™ Ethyl Vinyl Acetate (EVA) gravity bags are made with an EVA material that does not contain harmful plasticizers such as DEHP, which can leach into TPN solutions. EVA is more durable for reduced breakage, more stable for improved safety and more transparent for improved admixture inspection.


- Latex-free
- DEHP-free
- Detachable 2 or 3 leg transfer sets
- PVC ports and tubing
- Available in 2 or 3 leg versions
- Fill legs have vented spikes
- Cap for closing fill leg
- Available in multiple sizes: 1000mL 4000mL



Item #	Description	Size	UOM
IM68102	IMed EVA Bag, 2 Legs	500mL	50/Case
IM68100	IMed EVA Bag, 2 Legs	1000mL	50/Case
IM68105	IMed EVA Bag, 2 Legs	2000mL	50/Case
IM68110	IMed EVA Bag, 2 Legs	3000mL	50/Case
IM68120	IMed EVA Bag, 3 Legs	1000mL	50/Case
IM68125	IMed EVA Bag, 3 Legs	2000mL	50/Case
IM68130	IMed EVA Bag, 3 Legs	3000mL	50/Case
IM68135	IMed EVA Bag, 3 Legs	4000mL	50/Case

IMed Products[™] Smallbore Extension Sets are equipped with a male and female luer and pinch clamp.

- Latex-free
- DEHP-free








Item #	Description	Length	UOM
IM41000	Smallbore Extension Set	7"	50/Case
IM41001	Smallbore Extension Set	12"	50/Case
IM41002	Smallbore Extension Set	15"	50/Case

IMed Products[™] Extension Set with Micron Filter feature a male/female luer, pinch clamp, and micron filter.


- Latex-free
- DEHP-free



Item #	Description	Length	UOM
IM41067	Extension Set with .22 Micron Filter	12"	50/Case
IM41068	Extension Set with 1.2 Micron Filter	12"	50/Case


IMed Products[™] Extension Set with Flow Regulator features a male and female luer, Y-injection site 7" from distal end, and a priming volume of 3.1 mL.

- Latex-free
- DEHP-free
- Length: 19"
- Range of 5mL-250mL per hour



Item #	Description	UOM
IM41018	19" Extension Set with Flow Regulator	50/Case

IMed Products[™] gravity sets feature a universal vented/non-vented spike, male/female luer, Y-injection site, roller clamp, and a drip chamber.



Item #	Description	UOM
IM42070	IV Admin Set with .22 Micron In Line Filter	50/Box
IM42071	IV Admin Set, 10 Drops/mL	50/Case
IM42072	IV Admin Set, 15 Drops/mL	50/Case
IM42073	IV Admin Set, 20 Drops/mL	50/Case

IMed Products[™] IV Admin Sets with Flow Regulators feature a universal vented/non-vented spike, male luer, pinch clamp, and drip chamber.

- Latex-free
- DEHP-free









Item #	Description	UOM
IM42031	IV Admin Set with Flow Regulator	50/Box
IM42074	IV Admin Set with 5-300mL Flow Regulator	50/Box
IM42076	Needle-Free IV Admin Set with Flow Regulator	50/Box
IM42069	IV Admin Set with Flow Regulator and .2m Filter	50/Box

IMed Products[™] Secondary IV Administration Set features a hanger, universal spike, and roller clamp.

- Latex-free
- DEHP-free



Item #	Description	UOM
IM42075	IMed Secondary IV Admin Set	50/Box

IV START KITS - COMPONENT MATRIX

	IM1200C	IM1205C	IM1207C	IM1252	IM2200C	IM2201
Таре	Transparent	Transpore	Transparent	Transpore	Transparent	Transpore
Transparent Dressing	#1614	#9534	#1614	#9534	#LMP14	#LMP14
Gauze 2 x 2	1	1	2		1	2
Non-Woven Gauze				2		
Alcohol Prep Pad		1		1		1
PVP Prep Pad, 10%				1		
Chloraprep™		Freppp™ 1.5mL	1mL			1mL
Prevantics® Swab	1				1	
Tourniquet	1	1	1	1	1	1
Gloves	Nitrile	Nitrile	Nitrile	Nitrile	Nitrile	Vinyl
UOM	50/Case	50/Case	50/Case	100/Case	50/Case	50/Case

IV Start Kits

IV START KIT WITH PREVANTICS®

ITEM #: IM1200C

KIT CONTENTS:

- 1 Pair of Gloves
- 1 Tourniquet
- 1 Roll of Tape
- 2 Gauze, 2" x 2"
- 1 Prevantics® Swab
- 1 Transparent Dressing (3M Item #1614)
- 50/Case



IV START KIT WITH CHLORAPREP™ FREPP™

ITEM #: IM1205C

KIT CONTENTS:

- 1 Pair of Gloves
- 1 Tourniquet
- 1 Roll of Tape
- 2 Gauze, 2" x 2"
- 1 Alcohol Prep Pad
- 1 ChloraPrep™ Frepp™
- 1 Transparent Dressing (3M Item #9534)
- 50/Case



IV Start Kits

IV START KIT WITH CHLORAPREP™ 1ML ITEM #: IM1207C

KIT CONTENTS:

- 1 Pair of Gloves
- 1 Tourniquet
- 1 Roll of Tape
- 1 Towel
- 2 Gauze, 2" x 2"
- 1 ChloraPrep™, 1mL
- 1 Transparent Dressing (3M Item #1614)
- 50/Case



IV START KIT WITH PVP PREP PAD

ITEM #: IM1252

KIT CONTENTS:

- 1 Pair of Gloves
- 1 Tourniquet
- 1 Roll of Tape
- 2 Gauze, 2"x 2"
- 1 Alcohol Prep Pad
- 1 PVP Prep Pad
- 1 3M[™] Tegaderm[™] HP Transparent Dressing (3M #9534)
- 100/Case



IV Start Kits

IV START KIT WITH PREVANTICS®

ITEM #: IM2200C

KIT CONTENTS:

- 1 Pair of Gloves
- 1 Tourniquet
- 1 Roll of Tape
- 2 Gauze, 2" x 2"
- 1 Prevantics® Swab
- 1 Transparent Dressing (#LMP14)
- 50/Case



IV START KIT WITH CHLORAPREP™ 1ML

ITEM #: IM2201

KIT CONTENTS:

- 1 Pair of Gloves
- 1 Tourniquet
- 1 Roll of Tape
- 2 Gauze, 2"x 2"
- 1 Alcohol Prep Pad
- 1 ChloraPrep™, 1mL
- 1 Transparent Dressing (#LMP14)
- 50/Case



DRESSING CHANGE KITS - COMPONENT MATRIX

	IM1210	IM1211	IM1212	IM1219	IM1222	IM1250	IM46023	IM46061	IM46062
ChloraPrep™		Swabsticks 3/pack 5.25mL	3mL	3mL	3mL	Swabsticks 3/pack 5.25mL	3mL	3mL	3mL
10% PVP I-Prep Swabsticks (3/Pack)	1								
4mm Disc with CHG								1	1
Alcohol Prep Pad		2	2	2	2	1	2	2	2
Skin Prep Pad					1	1			
Statlock PICC Plus with Skin Prep							1	1	
Alcohol Swabsticks (3/Pack)	1						1	1	1
Таре	Transpore	Transpore	Transpore	Transpore	Transpore	Transpore	Transpore	Transpore	Transparent
Gauze 4x4	1	1	1	1	1	1	1	1	1
Gauze 2x2	1	1	1	1	1	1			
Non-Woven Split Gauze	1	1	1		1	1	1	1	1
Non-Adherent Gauze	1	1	1	1	1	1	1	1	1
Dressing	#1626W	#1626W	#1626W	#1658	#1616	#1616	#1616	#1616	#1616
Gloves	Nitrile	Nitrile	Nitrile	Nitrile	Nitrile	Vinyl	Nitrile (2)	Nitrile	Nitrile
Tape Measure	1	1	1	1	1	1	1	1	
Towel	1	1	1	1	1	1	1	1	1
Mask	1	1	1	1	1	1	2	1	1
иом	30/Case	30/Case	30/Case	20/Case	30/Case	30/Case	30/Case	30/Case	30/Case

DRESSING CHANGE KIT WITH PVP SWABS

ITEM #: IM1210

KIT CONTENTS

- 1 Pair of Gloves
- 1 Mask
- 1 Roll of Tape
- 1 Towel
- 1 Tape Measure
- 1 Gauze, 2" x 2"
- 1 Gauze, 4" x 4"
- 1 Non-Adherent Gauze, 2" x 2"
- 1 Split Non-Woven Gauze, 2" x 2"
- 1 Alcohol Swabsticks, 3/Pack
- 1 PVP Swabsticks, 3/Pack
- 1 3M[™] Tegaderm[™] Dressing (3M Item #1626W)
- 30/Case



DRESSING CHANGE KIT WITH CHLORAPREP™ SWABS

ITEM #: IM1211

KIT CONTENTS

- 1 Pair of Gloves
- 1 Mask
- 1 Roll of Tape
- 1 Towel
- 1 Tape Measure
- 1 Gauze, 2" x 2"
- 1 Gauze, 4" x 4"
- 1 Non-Adherent Gauze, 2" x 2"
- 1 Split Non-Woven Gauze, 2" x 2"
- 2 Alcohol Prep Pads
- 1 Chloraprep™ Swabsticks, 3/Pack
- 1 3M[™] Tegaderm[™] Dressing (3M Item #1626W)
- 30/Case



DRESSING CHANGE KIT WITH CHLORAPREP™ 3ML

ITEM #: IM1212

KIT CONTENTS

- 1 Pair of Gloves
- 1 Mask
- 1 Roll of Tape
- 1 Towel
- 1 Tape Measure
- 1 Gauze, 2" x 2"
- 1 Gauze, 4" x 4"
- 1 Non-Adherent Gauze, 2" x 2"
- 1 Split Non-Woven Gauze, 2" x 2"
- 2 Alcohol Prep Pads
- 1 Chloraprep[™], 3mL
- 1 3M[™] Tegaderm[™] Dressing (3M Item #1626W)
- 30/Case



DRESSING CHANGE KIT WITH CHG DRESSING

ITEM #: IM1219

KIT CONTENTS

- 1 Pair of Gloves
- 1 Mask
- 1 Roll of Tape
- 1 Towel
- 1 Tape Measure
- 1 Gauze, 2" x 2"
- 1 Gauze, 4" x 4"
- 1 Non-Adherent Gauze, 2" x 2"
- 2 Alcohol Prep Pads
- 1 Chloraprep[™], 3mL
- 1 Tegaderm CHG Insert
- 1 3M[™] Tegaderm[™] Dressing with CHG (3M Item #1658)
- 20/Case



DRESSING CHANGE KIT WITH

CHLORAPREP™ 3ML

ITEM #: IM1222

KIT CONTENTS

- 1 Pair of Gloves
- 1 Mask
- 1 Roll of Tape
- 1 Towel
- 1 Tape Measure
- 1 Gauze, 2" x 2"
- 1 Gauze, 4" x 4"
- 1 Non-Adherent Gauze, 2" x 2"
- 1 Split Non-Woven Gauze, 2" x 2"
- 2 Alcohol Prep Pads
- 1 Skin Protectant Pad
- 1 Chloraprep[™], 3mL
- 1 3M[™] Tegaderm[™] Dressing (3M Item #1616)
- 30/Case



DRESSING CHANGE KIT WITH CHLORAPREP™ SWABS

ITEM #: IM1250

KIT CONTENTS

- 1 Pair of Gloves
- 1 Mask
- 1 Roll of Tape
- 1 Towel
- 1 Tape Measure
- 1 Gauze, 2" x 2"
- 1 Gauze, 4" x 4"
- 1 Non-Adherent Gauze, 2" x 2"
- 1 Split Non-Woven Gauze, 2" x 2"
- 2 Alcohol Prep Pads
- 1 Skin Prep Pad
- 1 Chloraprep™ Swabsticks, 3/pack
- 1 3M[™] Tegaderm[™] Dressing with CHG (3M Item #1616)
- 30/Case



DRESSING CHANGE KIT WITH STATLOCK™

ITEM #: IM46023

KIT CONTENTS

- 2 Pairs of Gloves
- 2 Masks
- 1 Roll of Tape
- 1 Towel
- 1 Tape Measure
- 1 Gauze, 4" x 4"
- 1 Non-Adherent Gauze, 2" x 2"
- 1 Split Non-Woven Gauze, 2" x 2"
- 2 Alcohol Prep Pads
- 1 Alcohol Swabstick, Triple
- 1 StatLock® PICC Plus with Skin Prep
- 1 ChloraPrep™ 3mL
- 1 3M[™] Tegaderm[™] Dressing (3M #1616)
- 30/Case



KIT CONTENTS

- 1 Pair of Gloves
- 1 Mask
- 1 Roll of Tape
- 1 Towel
- 1 Tape Measure
- 2 Gauze, 4" x 4"
- 2 Alcohol Prep Pads
- 1 Alcohol Swabstick, Triple
- 1 Statlock™ PICC Plus with Skin Prep
- 1 HaloGUARD, 1" 4mm
- 1 ChloraPrep™ 3mL
- 1 Dressing Change Label
- 1 3M[™] Tegaderm[™] Dressing (3M Item #1616)
- 30/Case

DRESSING CHANGE KIT WITH STATLOCK™ & HALOGUARD

ITEM #: IM46061



DRESSING CHANGE KIT WITH CHLORAPREP™ 3ML & HALOGUARD

KIT CONTENTS

ITEM #: IM46062

- 1 Pair of Gloves
- 1 Mask
- 1 Roll of Tape
- 1 Towel
- 1 Tape Measure
- 1 Gauze, 4" x 4"
- 1 Non-Adherent Gauze, 2" x 2"
- 1 Non-Woven Gauze, 2" x 2"
- 1 Split Non-Woven Gauze, 2" x 2"
- 2 Alcohol Prep Pads
- 1 Alcohol Swabsticks, 3/Pack
- 1 ChloraPrep™ 3mL
- 1 HaloGUARD 1" Disc with CHG
- 1 3M[™] Tegaderm[™] Dressing (3M Item #1616)
- 30/Case



MISCELLANEOUS KITS - COMPONENT MATRIX

	IMHYP1000	IMPA1213D	IM1300A
OPTFlow Sub-Q 2-Needle Safety Set	1		
IV Flow Regulator	1		
Gloves	Nitrile	Vinyl	Nitrile (2)
Towel		1	
ChloraPrep™	FREPP™	3mL	
Alcohol Prep Pad	1		
Roll of Tape	Transpore		
Prefilled Saline Syringe		3	
Non-Woven Split Sponge		1	
Gauze 2x2		1	
Chemo Waste Bag			1
Antiseptic Towelette			1
Absorbent Towels			2
Poly-Coated ISO Gown			1
UOM	5/Case	30/Case	20/Case

Miscellaneous Kits

HYPODERMOCLYSIS KIT

ITEM #: IMHYP1000

KIT CONTENTS

- 1 Pair of Gloves
- 1 Roll of Tape
- 2 Alcohol Prep Pads
- 1 OPTFlow® Sub-Q 2-Needle Safety Set. 26G x 9mm
- 1 IV Flow Regulator Set
- 2 ChloraPrep™ Frepp™, 1.5mL
- 5/Case



PORT ACCESS KIT WITH CHLORAPREP™

ITEM #: IMPA1213D

KIT CONTENTS

- 1 Pair of Gloves
- 1 Towel
- 1 Gauze, 2" x 2"
- 1 Split Non-Adherent Gauze, 2" x 2"
- 3 Prefilled Saline Syringes
- 1 ChloraPrep™ 3mL
- 30/Case



Miscellaneous Kit

BIOHAZARD SPILL KIT

ITEM #: IM1300A

KIT CONTENTS

- 2 Pairs of Gloves
- 1 Chemo Waste Bag
- 1 Antiseptic Towelette
- 2 Absorbent Towels
- 1 Poly-Coated ISO Gown
- 20/Case



Custom Kits Kits

Can't find what you're looking for?

We offer custom kit solutions to ensure your needs are met in the most convenient and economical way possible.


For custom kit inquires please contact your Territory Manager.

Please visit our website at: integrated medsys.com to view our full line of available kits.



IMed Products Alcohol Prep Pads, a reliable solution for safe and effective skin preparation before injections or venipuncture. These premium-quality pads are saturated with 70% isopropyl alcohol, effectively combating bacteria and germs to reduce the risk of contamination. The non-woven prep pads are known for their exceptional softness and absorbency, ensuring enhanced performance. Each pad is designed for single use and is securely sealed in sterile pouches, ensuring the highest standards of hygiene.

- Sterile
- 2-ply
- Available in medium and large sizes
- Contains 70% isopropyl alcohol





Item #	Description	UOM
IM25092	Sterile Alcohol Prep Pads, Medium	20 Boxes/200; 4000/Case
IM25093	Sterile Alcohol Prep Pads, Large	40 Boxes/100; 4000/Case

The Cleo® 90 infusion system is a single-use disposable system to be used in a variety of infusion therapies (including insulin). This system makes insertions quick and convenient in the home. The enclosed needle in the cannula retracts automatically into the housing after insertion, reducing the risk of accidental needle-stick injuries. The tubing can be connected with a single motion using just one hand, and can be used with any pump that has a luer fitting.

PRODUCT FEATURES:

- Hypoallergenic dressing The small, transparent adhesive dressing breathes yet keeps water out – and stays in place with less tape. Its patented release tab makes removal easy
- Standard luer connection

Flexible 360° connection - Ability to rotate the connection of the tubing 360°, which
offers even more freedom and flexibility. The Flex Attach™ 360° site connection
provides on-site rotation offering greater flexibility

IM53081




Item #	Description	UOM
IM53080	Cleo® SubQ Infusion Cannula, 6mm	50/Box
IM53081	Cleo® SubQ Infusion Cannula, 9mm	50/Box

Item #	Description	UOM
IM53082	Cleo® SubQ Infusion Set Tubing, 24"	50/Box
IM53083	Cleo® SubQ Infusion Set Tubing, 31"	50/Box
IM53084	Cleo® SubQ Infusion Set Tubing, 42"	50/Box

IMed Products kits are convenient, economical, and have everything in one place, which prevents waste by providing exact quantities and supplies needed. This convenient Subcutaneous Dose Kit eliminates the need for component picking at the pharmacy.

- 2 Alcohol Prep Pads
- 2 Adhesive Bandages, 1x3"
- 1 SurGuard®3 Hypodermic Needle, 27g x 1/2"
- 1 SurGuard®3 Hypodermic Syringe with Safety Needle,1ml, 27g x 1/2"
- 1 SurGuard®3 Hypodermic Syringe with Safety Needle, 3ml, 25g x 5/8"








Item #	Description	UOM
IM46116	IMed Subcutaneous Dose Kit	100/Case

IMed Products[™] Empty Syringes are designed to provide accurate and effective dosing with an easy-to-read numbering system and transparent barrel. Minimal dead space allows for effective and smooth evacuations from the syringe barrel for increased patient comfort.

- Latex-free
- Sterile-packed in tamper-evident peel packs



Item #	Description	Size	UOM
IM54085	Empty Syringe with Luer Lock	1mL	100/Box
IM54086	Empty Syringe with Luer Lock	3mL	100/Box
IM54087	Empty Syringe with Luer Lock	5mL	100/Box
IM54088	Empty Syringe with Luer Lock	10mL	100/Box
IM54089	Empty Syringe with Luer Lock	20mL	50/Box
IM54090	Empty Syringe with Luer Lock	30mL	50/Box
IM54091	Empty Syringe with Luer Lock	60mL	50/Box

Enteral



IMed Products[™] Enteral Backpacks are durable, economical, and stylish for patients receiving nutrition via ambulatory enteral feeding pumps; ideal for patients on the go







Item #	Description	Size	UOM
IM58020	IMed Enteral Pump Backpack	500mL	1 Each
IM58009U	IMed Enteral Pump Backpack	1200mL	1 Each

IMed Products™ EnFit Syringes are designed to provide accurate and effective dosing with an easy-to-read numbering system and transparent barrel. Minimal dead space allows for effective and smooth evacuations from the syringe barrel for increased patient comfort.

- Latex-free
- Sterile-packed in tamper-evident peel packs



Item #	Description	Size	UOM
IM54092	EnFit Syringe	5mL	50/Box
IM54093	EnFit Syringe	10mL	50/Box
IM54094	EnFit Syringe	20mL	50/Box
IM54095	EnFit Syringe	35mL	50/Box
IM54096	EnFit Syringe	60mL	50/Box
IM54097	EnFit Syringe with Thumb Ring	60mL	50/Box

Cleanroom



IMed Products[™] 70% Isopropyl Alcohol (IPA) contains 70% by volume USP-grade Isopropanol (isopropyl alcohol) and 30% USP-purified water.

- Sterile
- Filled into cleaned containers, .1 micron filtered






Item #	Description	Size	UOM
IM27077	Sterile, 70% Isopropyl Alcohol	16oz	12/Case
IM27078	Sterile, 70% Isopropyl Alcohol	32oz	12/Case
IM27079	Sterile, 70% Isopropyl Alcohol	Gallon	4/Case
IM27094	Sterile, 70% Isopropyl Alcohol with Flip Top	16oz	12/Case

Respiratory



IMed Products™ Filters are single-use devices indicated for patient during the delivery of ventilator gasses. They provide filtration for reducing possible cross contamination between patient and equipment. This filter is for use in hospital, ICU, anesthesia, respiratory therapy, during transport, and with resuscitators. These filters can be used on adult patients and are indicated for use by qualified medical personnel only. These are designed for disposable use and should be changed every 24 hours.


- Latex-free
- DEHP-free



Item #	Description	UOM
IM55028	Viral & Bacterial Breathing Filter	50/Case
IM55029	Viral & Bacterial HEPA Filter	50/Case

IMed Products™ HME's offer a wide range of performance capabilities and configuration options to meet the unique needs of your patient's respiratory care.

- Latex-free
- DEHP-free



Item #	Description	UOM
IM55010	Straight In-Line HME	50/Case
IM55011	Pedi-Neo HME	50/Case
IM55012	HME Vent	50/Case
IM55030	HME with O2 & Suction Ports	50/Case



Customer Service: 800.755.3800 www.integratedmedsys.com/imed-products