

IMed Products™ Empty Syringes with Luer Lock Tips are designed for precise fluid and medication delivery in medical applications. Manufactured to meet the ISO 80369-7 standard, these syringes are available in sizes ranging from 1mL to 60mL. Constructed from high-quality, biocompatible materials, they offer durability, flexibility, and a secure connection for safe, reliable use.



## Features & Benefits

- Not Made with Natural Rubber Latex
- DEHP-Free
- Pyrogen-Free
- Non-Toxic
- Sterilized by Ethylene Oxide
- Individually-Packed
- Easy-to-Read Numbering
- Transparent Barrel
- Luer Lock Tip
- Meets ISO 80369-7

## Ordering Details

Item #	Description	Size	UOM
<b>IM54085*</b>	Empty Syringe with Luer Lock	1mL	100/Box
<b>IM54086</b>	Empty Syringe with Luer Lock	3mL	100/Box
<b>IM54087</b>	Empty Syringe with Luer Lock	5mL	100/Box
<b>IM54088</b>	Empty Syringe with Luer Lock	10mL	100/Box
<b>IM54089</b>	Empty Syringe with Luer Lock	20mL	100/Box
<b>IM54090*</b>	Empty Syringe with Luer Lock	30mL	50/Box
<b>IM54091</b>	Empty Syringe with Luer Lock	50/60mL	60/Box

\*Temporarily unavailable.

Proudly Distributed by:



Leading the healthcare industry towards a more efficient and effective future through our commitment to excellence in developing cutting-edge medical products and solutions.



To order, contact your IMS Territory Manager or Customer Service at 800.755.3800

The product details provided herein, including any images or supplements, reflect information from various sources regarding the product; this information has not been independently verified by Integrated Medical Systems ("IMS") or its affiliates. IMS is not liable for any inaccuracies or omissions in the product information, which is provided solely for your convenience.

The information presented here is subject to change without prior notice. Customers should review product packaging and labeling before use or purchase to ensure accuracy, especially when referencing information from IMS's online or printed materials. IMS disclaims any responsibility for the accuracy or completeness of the information reflected herein, regardless of whether any inaccuracies or omissions are due to errors or other factors.