

Microbial Resistance

**Bacterial spores
(*C.Diff*, Bacillus)**

**Use sporicidal disinfectants
once a week/2 weeks:**

- TexCide™ TX690 or
- TexTab™ TX6460 (sporicidal dilution)

Texwipe's Offering



**Mycobacterium
tuberculosis (TbC)
Non-lipid-coated viruses**

**Use bactericidal, virucidal
and fungicidal disinfectants
every day:**

- TexQ® Disinfectant or
- TexTab™ TX6460 (bactericidal dilution)



**Lipid-coated viruses
Vegetative bacteria
Fungi, fungal spores,
molds and yeast
HIV**

**Use cleaners and
sanitizers every day:**

- IPA (sterile, non-sterile)
- Denatured Ethanol
- Hydrogen Peroxide
- TexTab™ TX6460 (sanitizing dilution)



Basic microbes