

TABLE OF CONTENTS


03	AMBULATORY INFUSION PUMP
04	SPIKE SYSTEM
05	LUER LOCK SYSTEM
08	ACCESSORIES
11	ORDERING & CONTACT INFORMATION

AMBULATORY INFUSION PUMP NIMBUS™ II FLEX/PLUS

THE FUTURE OF AMBULATORY INFUSION THERAPY.

The Nimbus™ II Ambulatory Infusion System combines the accuracy of an electronic infusion pump with simple and intuitive programming.

This unique pump is engineered to be forward-thinking and capable of delivering advantages that meet the evolving needs of cancer care.




PRODUCT ID	PRODUCT NAME	QTY
IT1062	NIMBUS™ II FLEX PUMP	1
IT1078	NIMBUS™ II PLUS PUMP	1
INCLUDED WITH EITHER PUMP:		
SOL211	CLINICAL QUICK REFERENCE GUIDE	1
SOL212	PATIENT QUICK REFERENCE GUIDE	1
IT1069	NIMBUS BATTERY (EVE ER18505M 3.6V lithium battery)	1

ADMINISTRATION SETS & IV BAGS

SPIKE SYSTEM

HS-004

QTY: 20/CASE



(1) Luer Cap, (2) One Female Luer Lock Adaptor, (3) One Cassette, (4) One 1.2 Micron Filter, (5) One Slide Clamp, (6) One Backcheck Valve, (7) One Male Luer Lock Adaptor, (8) Luer Cap
 Priming Volume: 3mL (approx.)
 Length: 67" (170cm)

Our filters are made with low protein binding PES material.

SECURE® IV BAGS

QTY: 50/CASE





SOL202	SOL201	SOL203	SOL206	SOL207
50mL fill	100mL fill	250mL fill	500mL fill	1000mL fill

ADMINISTRATION SETS & IV BAGS

LUER LOCK SYSTEM

HS-008



QTY: 20/CASE



- (1) Luer Cap, (2) One Proximal Male Luer Lock Adaptor, (3) One Cassette, (4) One 0.2 Micron Filter, (5) One Slide Clamp, (6) One Backcheck Valve, (7) One Distal Male Luer Lock Adaptor, (8) Luer Cap
 Priming Volume: 2.5mL (approx.)
 Length: 53" (135cm)

HS-008-B

QTY: 20/CASE



- (1) Luer Cap, (2) One Proximal Male Luer Lock Adaptor, (3) One Cassette, (4) One 0.2 Micron Baby Filter, (5) One Slide Clamp, (6) One Backcheck Valve, (7) One Distal Male Luer Lock Adaptor, (8) Luer Cap
 Priming Volume: 2.5mL (approx.)
 Length: 53" (135cm)

Our filters are made with low protein binding **PES** material.

The Nimbus Administration Set is designed to administer fluids/medication from a container to a patient through a needle or catheter. The fluid path is sterile (EtO sterilized), non-pyrogenic, not made of DEHP, not made with natural latex. The set is single use only.

ADMINISTRATION SETS & IV BAGS

LUER LOCK SYSTEM

EVA BAGS

QTY: 50/CASE



SOL237

150mL fill



SOL238

250mL fill



SOL273

500mL fill



SOL247

105mL fill (for mini hardshell)

ADMINISTRATION SETS & IV BAGS

LUER LOCK SYSTEM

EVA BAGS 2 WAY NEEDLE-FREE VALVE QTY: 20/CASE



SOL260

150mL fill



SOL261

250mL fill



SOL262

500mL fill

CARRY POUCHES, HARDSHELL & MORE **ACCESSORIES**

CARRY POUCHES

QTY: EACH



SOL224

Fits Nimbus Pump with **50mL Spike System** OR **100mL LL System**



SOL247

Fits Nimbus Pump with **100mL Spike System** OR **250mL LL System**
with belt loop



while supplies last

SOL225

Fits Nimbus Pump with **100mL Spike System** OR **250mL LL System**



SOL226

Fits Nimbus Pump with **250mL Spike System** OR **500mL LL System**



SOL248

Fits Nimbus Pump with **500mL Spike System** OR **1L LL System**

CARRY POUCHES, HARDSHELLS & MORE ACCESSORIES

HARDSHELLS

QTY: 20/CASE



SOL228

For 100mL size
Max Bag fill is 110mL



SOL244

For 250mL size
Max Bag fill is 250mL



SOL254

For mini EVA bag
Max Bag fill is 105 mL

HARDSHELL POUCHES

QTY: EACH



SOL251

For 100mL hardshell with
window



SOL252

For 250mL hardshell with window



SOL265

For mini hardshell with window

CARRY POUCHES, HARDSHELLS & MORE **ACCESSORIES**

OTHERS



IT1069

Nimbus Battery

EVE ER18505M 3.6V lithium battery

QTY: EACH or 24/SLEEVE



SOL234

Chemo Spill Kit includes absorbant Bio-wipe Bag, Pr Chemo Gloves, Chemo Bag, Instructions

QTY: 5/PACK

INFORMATION

HOW TO ORDER & CONTACT US

OUR CUSTOMER SERVICE TEAM

Phone +1 (508) 650-2008
Fax +1 (508) 650-2006
Email flexorders@zynosolutions.com
Web www.zynosolutions.com

Operating days and hours: Monday - Friday, 8 AM - 5 PM

OUR TECH SUPPORT TEAM

Phone +1 (866) 395-1888
Email techsupport@zynosolutions.com
Web www.zynosolutions.com

Operating days and hours: 24/7

PREVENTATIVE MAINTENANCE

Customer is responsible for annual preventive maintenance procedures as required by the manufacturer. Zyno Solutions cannot assure safe use of the infusion pump if the customer fails to follow such procedures.

ADMINISTRATION SETS

Only manufacturer tubing is tested and calibrated for its products.



+1 (508) 650-2008
info@zynosolutions.com
177 Pine Street Natick MA 01760 USA
www.zynosolutions.com

ZYS-CT2F-002 Rev G