


# FreedomEdge®

Part of the Freedom Integrated Syringe Infusion System



## 123

### Easy to use & learn

No controls - just load, close and go!



### Safe constant pressure

Highly accurate constant pressure system, meeting ISO accuracy standards for safe drug delivery.



### DynEQ™ (Dynamic Equilibrium)

If the pressure in your sites increases, the system will immediately respond by automatically slowing down.



### Portable

No batteries or electricity needed. Comes with a travel pouch to promote patient independence.



### Quality assured

Made in the USA with the highest level of quality control for defect-free products. No maintenance needed.



### Fits BD® 20 & 30 ml syringes

Easy to handle and purge of air. Can be filled from 1 ml to 30 ml.



### Progress window

Easy to check infusion progress.



### Silent operation

Convenient for patients - no unnecessary alarms & noise.



### Ideal for all users

Maximum simplicity for all users - small hands, limited dexterity and limited strength.

**NOTE:** Use ONLY KORU Freedom System accessories with the FreedomEdge® Syringe Driver. Use of any other accessory may cause over delivery of fluids, unknown flow rates or additional site complications such as pain, swelling and redness.

## Related Freedom Integrated System Products:



### FREEDOM60® Syringe Infusion System

For larger infusion doses, BD® 50 ml syringe



### HiGH-Flo Subcutaneous Safety Needle Sets™

24G & 26G, 4 -14 mm lengths, 1-6 leg sets



### Precision Flow Rate Tubing™

For flow rates from F0.5 - F2400



### 24" Extension Set (FEXT)

For those who need longer tubing length



KORU Medical Systems | 845-469-2042 | korumedical.com | @korumedical

The FREEDOM60® and FreedomEdge® Syringe Infusion Systems, Precision Flow Rate Tubing™ and HiGH-Flo Subcutaneous Safety Needle Sets™ are trademarks of KORU Medical Systems and are compliant with Medical Device Directive 93/42/EEC. KORU Medical Systems is ISO 13485 certified. ©2020 KORU Medical Systems; All Rights Reserved. FE\_FL\_Benefits\_03-10007vC