

Product Specification Sheet

for

Secure[®] Fill Sets

 Gravity Fill Sets: Detachable Gravity Fill Legs with female screw connector For Secure[®] Legless bags with male screw connector This product is not made with natural rubber latex or DEHP Available in 1, 2 or 3 leg versions Fill legs have vented spike 							
Part #	Size (ml)	Description	Qty/Case	Case Price			
58794		Secure [®] 1 Leg Gravity Fill Set	50				
58782		Secure [®] 2 Leg Gravity Fill Set	50				
58783		Secure [®] 3 Leg Gravity Fill Set	50				

		Luer Fill Set:		
• H • 7	For Secure [®] Le This product is	Leg with female luer connector gless bags with male screw connector not made with natural rubber latex or DEHP ces requiring a female luer	<u>-</u>	C)#@
Part #	Size (ml)	Description	Qty/Case	Case Price
58784		Secure [®] 1 Leg Luer Adaptor (female)	50	

Manufactured by:

THE METRIX COMPANY 4400 Chavenelle Road Dubuque, IA 52002 800-752-3148

Distributed by: