


Product Specification Sheet for

Secure™ Dual Chamber TPN Bags Legless

Product Features:

- This product is not made with natural rubber latex or DEHP
- PVC container and tubing
- Stores fat emulsion separate from other TPN components for improved inspection and stability
- Legless Bags offer increased versatility
- Fill legs are available for most TPN and gravity filling applications
- 3000 ml (500ml/2500ml)
- 1500 ml (250ml/1250ml)



Part #	Size (ml)	Description	Qty/Case	Case Price
58713	3000	Secure™ Dual Chamber TPN Container - Legless	42	

Optional Accessories

Part #	Description	Qty/Case	Case Price
58781	1 Leg Compounder Fill Set	50	
58794	1 Leg Gravity Fill Set	50	
58782	2 Leg Gravity Fill Set	50	
58783	3 Leg Gravity Fill Set	50	
58784	1 Leg Luer Fill Set (female)	50	
58997	Non-Reopening Clamp	50/bag	

(All Secure™ EVA Legless Bags have a clamp attached, but some accounts require extra ones)

Manufactured by:

THE METRIX COMPANY
4400 Chavenelle Road
Dubuque, IA 52002
800-752-3148

Distributed by: