

IMed Products[®] Mini Transfer Pin is used to transfer fluids from source container vials. With a spike length of just 20 mm, it is perfect for accessing vials. This product also contains a 0.8 um vent filter to relieve vacuum pressure and provide easy transferring of fluids from the source container.

FEATURES & BENEFITS

- Not made with natural rubber latex
- DEHP-Free
- Sterile
- Non-pyrogenic
- 20 mm Spike with Cap
- Female Luer
- 0.8um Vent Filter
- Quantity: 100/Box
- UOM: Box



ORDERING DETAILS

ITEM #	DESCRIPTION	UOM
IM27127	Mini Transfer Pin	100/Box

Strengthened by Nationwide Distribution and Dedicated Customer Support

Proudly
Distributed by

IMS
INTEGRATED MEDICAL SYSTEMS



IMed[®]
Products

IMed Products delivers reliable, high-quality medical support solutions at value-driven pricing—designed to ensure consistent supply, trusted clinical performance, and operational ease for outpatient centers, clinics, hospitals, and pharmacies.

The product details provided herein, including any images or supplements, reflect information from various sources regarding the product; this information has not been independently verified by Integrated Medical Systems ("IMS") or its affiliates. IMS is not liable for any inaccuracies or omissions in the product information, which is provided solely for your convenience.

The information presented here is subject to change without prior notice. Customers should review product packaging and labeling before use or purchase to ensure accuracy, especially when referencing information from IMS's online or printed materials. IMS disclaims any responsibility for the accuracy or completeness of the information reflected herein, regardless of whether any inaccuracies or omissions are due to errors or other factors.

© 2026 IMed® Products. All rights reserved. 02/26 IMS25086 IMed Sell Sheets_Brand Refresh Mini Transfer Pin_01