

Get User-Friendly Options to Help Manage Hazardous Ambulatory Infusions

Enhance safety throughout the ambulatory pump process with CSTDs designed to help protect from beginning to end

Maximizing the safety of ambulatory IV oncology starts with minimizing your exposure to hazardous medications. We can help with a full line of extension sets with bonded closed system transfer devices (CSTDs) that feature intuitive, needlefree connections.


- › Intuitive closed system design requires fewer components and helps minimize exposure to hazardous drugs, ensuring USP <800> compliance
- › Bonded extension sets connect directly to the ambulatory pump cartridge to help improve workflow and reduce the risk of accidental disconnects and drug spill
- › Needlefree technology gives you a wide fluidpath to facilitate air removal and fluid transfer



ChemoLock™

Needlefree Closed System Transfer Device

Get a mechanically closed system with a membrane-to-membrane seal



ChemoClave™

Needlefree Closed System Transfer Device

Help keep your connections secure with familiar luer-based CSTDs




Improve Efficiency by Automating Air Removal and Drug Transfers


Manually removing air from ambulatory pump cassettes can be time consuming and result in repetitive stress injuries. With the Diana™ Automated Compounding Workflow System, you can easily and automatically remove air and transfer IV oncology drugs.

Get a Full Portfolio of Needlefree CSTD Bonded Kits


Get everything you need to prepare ambulatory cassettes, including intuitive needlefree components and convenience kits.




Item #	Length (in)	Case Qty
CL4120	61	50
Description		
Oncology Kit w/61" (155 cm) Smallbore Ext Set, ChemoLock™ w/Red Cap, Anti-Siphon Valve, ChemoLock™ Port, Y-Connector, Clamp, Rotating Luer; Double ChemoLock™ Syringe Transfer Set; ChemoLock™		



Item #	Length (in)	Case Qty
CL3960	60	50
Description		
Oncology Kit w/60" (152 cm) Smallbore Ext Set w/ ChemoLock™ Port, Clamp, Anti-Siphon Valve, Spiros™; ChemoLock™; Syringe Transfer Set w/MicroClave™, ChemoLock™ Port		



Item #	Length (in)	Case Qty
CL4151	45	50
Description		
45" (114 cm) Smallbore Ext Set w/Spiros™ w/Red Cap, Anti-Siphon Valve, Clamp, Y-Connector w/ChemoLock™ Port, Rotating LuerLock		



Item #	Length (in)	Case Qty
CH3606	60	50
Description		
Oncology Kit w/60" (152 cm) Smallbore Ext Set w/Spiros™, Anti-Siphon Valve, ChemoClave™, Clamp, Rotating Luer; Double ChemoClave™ Syringe Transfer Set; Spiros™ w/Red Cap		