

Pharmacy Dropper Bottle Series



- **PACKAGED STERILE DROPPER BOTTLES WITH VARIETY OF PACKAGING OPTIONS**
 - Single Pack, Dual Pack, Twelve Pack and Higher Options available
 - Custom Packaging configurations can be accommodated to optimize the efficiencies of any filling line set-up

- **VARIETY OF PACKAGING AND SIZING AVAILABLE:**
 - Natural colors available for the bottle portion, tips, and caps.

- **EXTENDED DESIGN OF DROPPING TIPS ALLOW FOR RELIABLE, REPEATABLE UNIFORM DROP SIZES**

- **FRICTION AND SNAP FIT TIPS PROVIDE A SECURE FIT TO THE BOTTLE**

- **VARIETY OF CUSTOM OPTION AVAILABLE IN SIZES RANGING FROM 1 ML TO 125 ML**

- **BOTTLE AND TIP (FLUID CONTACT SURFACES) COMPOSED OF LOW-DENSITY POLYETHYLENE (LDPE)**
 - Brittleness Temperature: < -75 ° C
 - Material meets standards set by the United States Pharmacopoeia (USP 22), the Biological Reactivity Tests for Class VI plastics and the requirements of USP 24 <661> for LDPE
 - Listed on the U.S. FDA Type III Drug Master File #1572
 - Listed on the Canadian Drug Master File #90147

- **PRODUCTS COME PRE-STERILIZED VIA ETO AND ARE GAMMA COMPATIBLE**

Pharmacy Dropper Bottle Series

NATURAL COLOR BOTTLES STOCKED AND READY TO SHIP TODAY:

Bottle Size	Bottle Color	Bottle Portion Diameter x Height (mm)	Dropper Tip Accuracy
3 mL	Natural	13 mm x 39 mm	44µL ±6µL
7 mL	Natural	19 mm x 51 mm	44µL ±6µL
10 mL	Natural	25mm x 45mm	44µL ±6µL
15 mL	Natural	25 mm x 66 mm	44µL ±6µL
30 mL	Natural	31 mm x 70 mm	44µL ±6µL

AVAILABLE ITEMS:

Catalog Number	Description	Number of Bottles Per Pack
AC211631-12	3 mL bottles (Natural), pre-assembled (12/pouch)	12 packs of 12 bottles each
AC211631-2	3 mL bottles (Natural), pre-assembled (2/pouch)	72 packs of 2 bottles each
AC211632-1	7 mL bottles (Natural), pre-assembled (1/pouch)	144 packs of 1 bottle each
AC211632-12	7 mL bottles (Natural), pre-assembled (12/pouch)	12 packs of 12 bottles each
AC211635-1	10 mL bottles (Natural), pre-assembled (1/pouch)	144 packs of 1 bottle each
AC211635-12	10 mL bottles (Natural), pre-assembled (12/pouch)	12 packs of 12 bottles each
AC211633-1	15 mL bottles (Natural), pre-assembled (1/pouch)	144 packs of 1 bottle each
AC211633-12	15 mL bottles (Natural), pre-assembled (12/pouch)	12 packs of 12 bottles each
AC211634-1	30 mL bottles (Natural), pre-assembled (1/pouch)	144 packs of 1 bottle each
AC211634-12	30 mL bottles (Natural), pre-assembled (12/pouch)	12 packs of 12 bottles each

CUSTOM OPTIONS (NOT LISTED) ARE AVAILABLE, WITH VARIETY OF PACKAGING, SIZING, AND COLORED CAPS, BOTTLES, AND TIPS. CONTACT MPS FOR ADDITIONAL INFORMATION.