

## TECHNICAL DATA SHEET

### Poly-Coated Impervious Isolation Gown: AAMI PB70 Level 2 Barrier Protection

#### Intended Use

Disposable Impervious Isolation Gown designed to protect health professionals and patients from cross infection.


#### Claim

Protective ultrasonically sealed gown that meets AAMI PB70 level II barrier performance for protective apparel and tested for chemotherapy drug permeation as per ASTM D6978.

Characteristic		Specification			
Material Composition		Polypropylene (PP) & Polyethylene (PE). Latex-Free			
Description		Gown made of non-woven material with full back closure. Features inclined sleeves, flat belt, elastic cuffs (elastic does not go all around the cuff) with thumb loops available, Over-The-Head (OTH) flexible neck closure made of two elastic strings.			
Part of Gown		Dimension (cm)			
		Small	Medium	Large	X-Large
A	Gown Length	100	105	112	129
B	Underarm Width	63	67	71	79
	Open Gown Width	126	134	142	158
C	Underarm Height	63	68	76	90
D	Sleeve Length (Excludes Thumb Loop)	Up: 65 Down: 50	Up: 72 Down: 56	Up: 76 Down: 59	Up: 76 Down: 60
E	Shoulder Height	32	36	38	40
F	Shoulder Width	10	10	10	13
G1	Belt Length (Front)	87	87	87	87
G2	Belt Length (Back)	93	93	93	93
	Belt Width	5	5	5	5
	Belt Position	20 cm under the arm, 7 cm towards the front	20 cm under the arm, 7 cm towards the front	20 cm under the arm, 7 cm towards the front	20 cm under the arm, 7 cm towards the front
H	Cuff Width (Stretched)	14	14	14	14
I-1	Trim Length	54	54	54	54

Part of Gown		Dimension (cm)			
		Small	Medium	Large	X-Large
I-2	Elastic Strings	9	9	9	9
I-3	Collar Height	10	10	10	10
I-4	Collar Width	19	19	22	22
	Tag Position	3,0 cm diameter, 2 cm from cuff edge	3,0 cm diameter, 2 cm from cuff edge	3,0 cm diameter, 2 cm from cuff edge	3,0 cm diameter, 2 cm from cuff edge
	Thumbloop Hole	23,5 (Left side of back only)	23,5 (Left side of back only)	23,5 (Left side of back only)	23,5 (Left side of back only)

Measurement Diagram



Refer to specific  
attribute for position  
of the belt

## Additional Information

Storage Condition	Country of Origin
Store under ambient temperature and dry condition	China

## Appendix 1 List of Item Codes and Attributes

## Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment
A69949	AMD-Ritmed Assure Wear VersaGown Impervious Isolation Gown	10 Gowns/Bag, 10 Bags/Case, 100 Gowns/Case	Small	Blue	Belt on the left side	Non-Sterile
A69950	AMD-Ritmed Assure Wear VersaGown Impervious Isolation Gown	10 Gowns/Bag, 10 Bags/Case, 100 Gowns/Case	Medium	Blue	Belt on the left side	Non-Sterile
A69951	AMD-Ritmed Assure Wear VersaGown Impervious Isolation Gown	10 Gowns/Bag, 10 Bags/Case, 100 Gowns/Case	Large	Yellow	Belt on the left side	Non-Sterile
A69954	AMD-Ritmed Assure Wear VersaGown Impervious Isolation Gown	10 Gowns/Bag, 10 Bags/Case, 100 Gowns/Case	Large	Blue	Belt on the left side	Non-Sterile
A69955	AMD-Ritmed Assure Wear VersaGown Impervious Isolation Gown	10 Gowns/Bag, 10 Bags/Case, 100 Gowns/Case	X-Large	Blue	Belt on the left side	Non-Sterile