



Sterile Cleanroom Wipers

Our sterile cleanroom wipers are sterile gamma irradiated to a 10⁻⁶ sterility assurance level. They are constructed of an extremely low linting 68gsm 55/45 cellulose blend. They are appropriate for critical areas up to and including ISO 5 cleanroom environments. The Cleanroom Wipers feature maximum absorbencies of 20-25 mL, and are available in 9" x 9" and 12" x 12" sizes, and available in different quantities.

Product # TCBWIP09

300 wipers per pack, 12 packs per case, a total of 3600 wipers per case

Product # TCBWIP09SP-20

20 wipers per bag, 15 bags per bundle, 12 bundles per case, a total of 3600 wipers per case

Product # TCBWIP12

150 wipers per pack, 8 packs per case, a total of 1200 wipers per case

Product # TCBWIP12SP-20

20 wipers per bag, 8 bags per bundle, 8 bundles per case, a total of 1280 wipers per case



Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBWIP09
Product # TCBWIP09SP-20
Product # TCBWIP12
Product # TCBWIP12SP-20