






PureWay Mailback Solutions

PureWay Sharps Mailback Solutions




The PureWay Mailback Service provides a complete turnkey collection and disposal process for simplified management of used needles, syringes and other small quantities of medical waste. The PureWay Mailback System is comprised of a custom designed, FDA approved mailable sharps container, a postage prepaid return-by-mail shipping box, manifest and instructions for use.

When the drug therapy is concluded or when the sharps container is full the patient or nurse places the sharps container back into the provided return packaging box. The return packaging and the enclosed sharps container are then mailed to a PureWay designated treatment facility for disposal. Upon receipt, PureWay documents the weight and destruction of each system. Hard copy proof of disposal is available to the generator along with being available online at purewayportal.com.



Item Number	Item Description	Case Ordering
	<p>1.4 Quart Sharps Disposal Mailback System Item: 80014</p> <p>1.4 Quart container will hold approximately 50-70 1cc syringes. includes One - 1.4 Quart sharps container, postage prepaid return packaging box, disposal of the waste along with documentation.</p>	<p>Order SKU# 80014CS</p> <p>Quantity: 12 Kits to shipping case.</p>
	<p>1.2 Gallon Sharps Disposal Mailback System Item: 40001</p> <p>1.2 Gallon sharps container will hold approximately 250-270 1cc syringes. Includes One - 1.2 gallon sharps container, postage prepaid return packaging box, disposal of the waste along with documentation.</p>	<p>Order SKU# 40001CS</p> <p>Quantity: 6 Kits to shipping case.</p>
	<p>2 Gallon Sharps Disposal Mailback System Item: 40002</p> <p>2 Gallon sharps container will hold approximately 530-570 1cc syringes. Includes One - 2 gallon sharps container, postage prepaid return packaging box, disposal of the waste along with documentation.</p>	<p>Order SKU# 40002CS</p> <p>Quantity: 4 Kits to shipping case.</p>
	<p>3 Gallon Sharps Disposal Mailback System Item: 40003</p> <p>3 Gallon sharps container will hold approximately 840-860 1cc syringes. Includes One - 3 gallon sharps container, postage prepaid return packaging box, disposal of the waste along with documentation. Accommodates 100 Mil cassettes for disposal.</p>	<p>Order SKU# 40003CS</p> <p>Quantity: 3 Kits to shipping case.</p>

Item Number	Item Description	Case Ordering
	<p>18 Gallon Trace Chemo Disposal Mailback System Item: 40018TC</p> <p>18 Gallon Trace Chemo Mailback disposal system. Includes one-18 gallon puncture-resistant container, postage prepaid return packaging box, disposal of the waste along with documentation.</p>	<p>Order SKU# 40018TC</p> <p>Quantity: 1 Kit to shipping case.</p>

Multi-Pack Configurations

	<p>1.4 Quart Sharps Disposal Mailback System - 6 Pack Item: 40016</p> <p>Includes six 1.4 Quart sharps container and one postage prepaid return-by-mail shipping box. Once all 6 containers are filled, send to PureWay for disposal. Each 1.4 Quart container will hold approximately 50-70 1cc syringes.</p>	<p>Order SKU# 40016</p> <p>Quantity: 1 Kit to shipping case.</p>
	<p>1.2 Gallon Sharps Disposal Mailback System Twin Pack Item: 40011</p> <p>Includes two 1.2 Gallon sharps containers and one postage prepaid return-by-mail shipping box. Once both containers are filled, send to PureWay for disposal. Each 1.2 Gallon container will hold approximately 250-270 1cc syringes.</p>	<p>Order SKU# 40011</p> <p>Quantity: 1 Kit to shipping case.</p>
	<p>2 Gallon Sharps Disposal Mailback System Twin Pack Item: 40022</p> <p>Includes two 2 Gallon sharps containers and one postage prepaid return-by-mail shipping box. Once both containers are filled, send to PureWay for disposal. Each 2 Gallon sharps container will hold approximately 530-570 1cc syringes.</p>	<p>Order SKU# 40022</p> <p>Quantity: 1 Kit to shipping case.</p>

Mailback Pail And Tub Solutions

	<p>5 Gallon Medical Waste Mailback Pail System Item: 40005</p> <p>Use for red bag waste or small existing sharps container that may be present. Includes One - 5 gallon mailback pail, postage prepaid return packaging box, disposal of the waste along with documentation of disposal. Not for loose sharps.</p>	<p>Order SKU# 40005</p> <p>Quantity: 1 Kit</p>
	<p>8 Gallon Medical Waste Mailback Tub System Item: 40008</p> <p>Use for red bag waste or small existing sharps container that may be present. Includes One - 8 gallon mailback tub, postage prepaid return packaging box, disposal of the waste along with documentation of disposal. Not for loose sharps.</p>	<p>Order SKU# 40008</p> <p>Quantity: 1 Kit</p>

PureWay mailback system advantages:

- Our containers accommodate 20-25% more medical waste than competitors.
- Smaller case pack is better for your budget, warehousing storage needs.
- Manifests preprinted for your convenience.
- Freight included in pricing to branch

Contact PureWay at **877-765-3030** to learn more about our mailback solutions for your small quantity medical waste, nursing or patient generated waste needs.