


Respiratory Care Product Guide



Solutions across the continuum of care.



We help meet every need—
from hospital to home

Medline + Hudson RCI



Your partner for world-class respiratory solutions and service

Optimize your operations with our comprehensive portfolio. Driven by best practices, our products deliver the uncompromising quality and performance required to care for the complex needs of your patients across the continuum of care.

We control the entire supply chain, which enables us to reliably deliver the products you need—when and where you need them.

30+

Respiratory
Specialists



50+

distribution centers
across North America



1,300+

MedTrans trucks
nationwide



30+


million

square feet of
warehousing



20+

North American
manufacturing sites



**1-day ship to
99% of the US**

with an average 60-day
inventory on hand for each
customer



200+

clinicians
on staff



We help you elevate respiratory care

Medline offers the depth and breadth of products to meet your needs with innovative solutions to enhance patient experience, improve outcomes and increase efficiency. Our solutions and products are clinically supported by our Respiratory Specialists and delivered by our unparalleled supply chain.

Humidification—Active

Neptune® Heated Humidifier	5
Concha™ Inhalation Water	6
Mounting Brackets	7
Heated Wire Breathing Circuits	8–9

Humidification—Passive

Filters	11
Bypass HMEs	12
HMEs & Filtered HMEs	13–14
Non-Heated Breathing Circuits	15–18

Non-Invasive Interfaces

Non-Invasive Ventilation Masks	20
Circuits & Accessories	21
CPAP Masks	22–24
Go-PAP™ Systems	25
O ₂ -Max™ Systems	25–26
Babi Plus® Infant Nasal CPAP	27

Heated High Flow

Comfort Flo® Plus Nasal Cannula	29
High Flow Humidification Systems	30
MaxVenturi High Flow	31
MicroMax™ High Flow	32
MicroMax™ Low Flow	33
Accessories	34

Resuscitation

Manual Resuscitators	36–39
Colorimetric CO ₂ Detectors	40
Resuscitation Masks & Accessories	41
Hyperinflation Systems	42

Airway Management

AG Cuffill®	44
Closed Suction	45–47
Open Suction	48
Addipak® Sterile Solutions	49
Trach Care & Cleaning Trays	50–52
Tracheostomy Tube Holders	53
SubG™ Endotracheal Tubes	54
Mucus/Specimen Traps	55

O₂ Delivery

OxyMask™	57–58
O ₂ Masks	58
Oxygen and Aerosol Masks	58–60
Face Tents	61
Oxygen Cannulas	62–65
Oxygen Supply Tubing	66
CO ₂ Cannulas & Masks	67–69
Bubble Humidifiers	70
Large Volume Nebulizers	71
Corrugate Tubing	72
Oxygen Kits	73

Medication Delivery

Medication Nebulizers	75–77
Aerosol Masks	78
MDI Chambers/Spacers	79–80
Nebulizer Compressors	81–82

Bronchial Hygiene

Incentive Spirometers	84
Peak Flow Meters	85
PEP Devices	86

Adapters and Connectors

General	88–89
Elbows	90
Wyes	91
MDI Adapters	91–92
NeoNatal Adapters	92
ConchaTherm Adapters	93
Drainage Systems	93
One-Way Valves	94
Temperature Monitoring	94
Water Traps	95
Tees	95–96
Disposable Nose Clips	96
Disposable Mouthpiece	96
Miscellaneous	97

Minor Capital

Fingertip Pulse Oximeters	99
Suction Regulators	100–101
Air & Oxygen Flowmeters	102–103
Oxygen Regulators	104
Oxygen Analyzers	105
Suction Aspirators	106

Humidification— Active

Neptune Heated Humidifier	5
Concha Inhalation Water	6
Mounting Brackets	7
Heated Wire Breathing Circuits	8-9



Hudson RCI Neptune® Heated Humidifier

The Neptune Heated Humidifier provides humidification levels for patients—from neonates to adults—with both invasive and non-invasive ventilation, and features an Auto Settings mode, low water notification, and expanded temperature range to help clinicians improve the quality of patient care.

- Auto Settings mode allows for pre-set and lock desired settings while still allowing for mode selection
- Low water notification feature lets users know when water needs to be replaced
- Temperature range for the Invasive Mode is 32°C–40°C and Non-Invasive Modes 28°C–37°C
- Heated high flow cannula therapy flow rates ranging from 1–60 liters per minute (LPM)



Item No.	Description	Pkg.
HUD42500	Neptune Heated Humidifier with ConchaSmart Technology	1 ea
HUD872DTPCS	Disposable Temperature Probe	10/cs
HUD39590	Right-angle Dual Temperature Probe	1 ea
HUD39599	Straight Dual Temperature Probe	1 ea
HUD38681	Water Reservoir Bracket	1 ea
HUD42599	Diagnostic Probe Kit	1 ea



Humidification Essentials

A suite of Hudson RCI® Humidification Essentials designed to promote patient safety and clinician efficiency, including ventilator circuits that support AARC clinical practice guidelines.

Hudson RCI Concha® Reservoirs

- The Hudson RCI Concha Reservoir is a self-contained plastic bottle with sterile water for inhalation
- For use exclusively with the Neptune Heated Humidifier
- Features unique 1,650 mL sterile water packaging and connections
- Water Source: USA

Item No.	Description	Pkg.
HUD38150	Concha Reservoir Water, 1,650 mL	6/cs
HUD38520	Conchapak Humidification System with 1,650 mL Concha Reservoir and Low-Compliance ConchaColumn Column with 1 Vertical and 1 Right-Angle Port	6/cs
HUD38530	Conchapak Humidification System with 1,650 mL Concha Reservoir and Low-Compliance ConchaColumn Column with 2 Vertical Ports	6/cs
HUD38560	Conchapak Humidification System with 1,650 mL Concha Reservoir and Standard ConchaColumn Column with 1 Vertical and 1 Right-Angle Port	6/cs
HUD38570	Conchapak Humidification System with 1,650 mL Concha Reservoir and Standard ConchaColumn Column with 2 Vertical Ports	6/cs



HUD38150

Hudson RCI ConchaSmart® and Concha Columns

- ConchaSmart® low water notification lets users know when water needs to be replaced

Item No.	Description	Pkg.
HUD38210	ConchaSmart Column	10/cs
HUD38270	Concha Standard Column with Vertical Ports	48/cs
HUD38230	Concha Low Compliance Column	48/cs



HUD38210

Hudson RCI Mounting Brackets

Mounting Brackets provide multiple mounting capabilities for Hudson RCI® Neptune® Heated Humidifiers.

Item No.	Description	Pkg.
HUD38610	Universal Pole Mount	1 ea
HUD38673	Rail Mounting Bracket	1 ea
HUD38675	Universal Ventilator Bracket; for use with Siemens®, Puritan Bennett®, and Respironics™ Esprit™ Ventilators	1 ea
HUD38682	Ventilator Bracket; for use with Puritan Bennett® 840 Ventilators	1 ea
HUD38688	Ventilator Bracket; for use with MAQUET® SERVO-i® Ventilators	1 ea
HUD38681	Concha® Reservoir Bracket	1 ea
HUD38692	Dräger® Bracket Adapter with Stop	1 ea
HUD38693	Universal Pole Bracket	1 ea
HUD38694	Dräger® Water Reservoir Bracket Adapter	1 ea



HUD38675



HUD38681



HUD38610



HUD38693




HUD38682

Heated Wire Breathing Circuits

- Hudson RCI® Breathing Circuits are available in a variety of configurations
- Proprietary stranded wire technology enhances circuit safety and durability
- Designed to provide humidification, while minimizing condensation in inspiratory and expiratory limbs¹

ConchaSmart® and ISO-Gard® Circuits

Hudson RCI ConchaSmart® and ISO-Gard® are validated for a 30-day useful life (adult circuits only). These circuits are designed to reduce circuit breaks, confirm humidification at the patient wye, and streamline clinician efficiency.¹ The ISO-Gard Circuit helps clinicians address the unique humidification needs of each patient, while minimizing circuit breaks and managing condensation.



*Only applicable for adult ISO-Gard and ConchaSmart Circuits

1. Restrepo RD, Walsh BK. AARC CPG: Humidification during invasive and noninvasive mechanical ventilation. *Respiratory Care*. 2012; 57 (5): 782-788.
2. Branson, Richard D. Conditioning Inspired Gases: The Search for Relevant Physiological End Points. *Respiratory Care*. 2009; 54 (4).

ISO-Gard Circuits[†]

Item No.	Description	Circuit Configurations								Length	Pkg.
		Limb	H-IL	H-EL	ILD	ELD	RTP	PMT			
HUD88001KIT	Infant Circuit**	Single	•	•			18"	•	70"	20/cs	
HUD88009KIT	Infant Circuit**	Dual	•	•	•	•		•	65"	20/cs	
HUD88013KIT	Infant Circuit**	Dual	•				18"	•	48"	20/cs	
HUD88014KIT	Infant Circuit**	Dual	•		•	•	18"	•	48"	20/cs	
HUD88015KIT	Infant Circuit**	Dual	•	•	•	•	18"	•	48"	20/cs	
HUD88017KIT	Infant Circuit**	Dual	•	•			18"	•	48"	20/cs	
HUD88023KIT	Pediatric Circuit**	Dual	•					•	66"	10/cs	
HUD88024KIT	Pediatric Circuit**	Dual	•	•	•	•		•	66"	10/cs	
HUD88025KIT	Pediatric Circuit**	Dual	•	•				•	66"	10/cs	
HUD88028KIT	Pediatric Circuit**	Dual	•		•	•		•	66"	10/cs	
HUD88019KIT	Adult Circuit	Single	•		•				70"	10/cs	
HUD88033KIT	Adult Circuit	Dual	•			•			70"	10/cs	
HUD88034KIT	Adult Circuit	Dual	•		•	•			70"	10/cs	
HUD88036KIT	Adult Circuit	Dual	•	•	•	•			70"	10/cs	
HUD88037KIT	Adult Circuit	Dual	•	•		•			70"	10/cs	
HUD88098KIT	Adult Circuit	Single	•		•			•	70"	10/cs	
HUD800	ISO-Gard Suction Wand	n/a								50/cs	
HUD810	ISO-Gard Drain—Neonate	n/a								25/cs	
HUD820	ISO-Gard Drain—Pediatric	n/a								25/cs	
HUD830	ISO-Gard Drain—Adult	n/a								25/cs	




HUD88019KIT

*All ISO-Gard Circuits and ConchaSmart Circuit Kits are packaged with ConchaSmart Columns.
[†]Heated-Wire Circuits are designed specifically for use with ConchaSmart and ConchaTherm (21 volt Neptune Heated Humidifiers).

Circuit Matrix Keys:

- EL-Trap = Expiratory Limb Water Trap
- H-IL = Heated Inspiratory Limb
- H-EL = Heated Expiratory Limb
- H-LIMB = Humidifier Limb
- ILD = Inspiratory Limb Drain, ELD = Expiratory Limb Drain
- IL-TRAP = Inspiratory Limb Water Trap
- RTP = Remote Temperature Monitoring Port
- ** = Pressure Monitoring Tubing



HUD87035KIT

ConchaSmart Circuits with Column Kits[†]

Item No.	Column	Description	Size	Circuit Configurations								Pkg.	
				Limb	H-IL	H-EL	RTP	60" L	48" L	PMT	H-Limb		
HUD78007KIT	ConchaSmart Column	Infant	10 mm	Dual	•	•	18"			•	•	•	20/cs
HUD78008KIT	Low Compliance	Infant	10 mm	Dual	•	•				•	•	•	20/cs
HUD78009KIT	Low Compliance	Infant	10 mm	Dual	•	•				•	•	•	20/cs
HUD78018KIT	ConchaSmart Column	Infant	10 mm	Dual	•	•	18"			•	•	•	20/cs
HUD78001KIT	ConchaSmart Column	Infant	10 mm	Single	•		18"			•	•	•	20/cs
HUD78017KIT	ConchaSmart Column	Infant	10 mm	Dual	•	•	12"			•	•	•	20/cs
HUD78024KIT	ConchaSmart Column	Pediatric	15 mm	Dual	•	•				•	•	•	10/cs
HUD87035KIT	ConchaSmart Column	Adult	22 mm	Dual	•	•				•	•	•	10/cs
HUD87098KIT	ConchaSmart Circuit	Adult	22 mm	Single	•					•			10/cs
HUD87098KITF	ConchaSmart Column	Adult	22 mm	Single	•					•			10/cs
HUD87019KIT	ConchaSmart Column	Adult	22 mm	Single	•					•	•		10/cs
HUD78096	N/A	Adult	22 mm	Single	•								15/cs

[†]Heated-Wire Circuits are designed specifically for use with ConchaSmart and ConchaTherm (21 volt) Neptune Heated Humidifiers.



Humidification— Passive

Filters	11
Bypass HME	12
HMEs & Filtered HMEs	13–14
Non-Heated Breathing Circuits	15–18

Filters

Bacterial/Viral Filters

Medline's Hydrophobic Filters feature translucent housings that reveal built-up secretions or other obstructions, for respiratory and anesthesia applications.

Performance Specifications

	DYNJAABV1	DYNJAABV100	DYNJAAHF
Description	Bacterial/Viral Filter	Bacterial/Viral Filter	HEPA Filter
Tidal Volume (mL)	≥250	≥250	≥200
BFE (%)	99.999	99.999	99.99999
VFE (%)	99.99	99.99	99.9999
Volume (mL)	34	34	63
Weight (g)	15.1	15.1	40
Flow Resistance @ 30 LPM	0.71	0.71	1.3
Flow Resistance @ 60 LPM	1.02	1.02	2.89
CO₂ Port	No	No	Yes
Filter Type	Electrostatic	Electrostatic	HEPA
Pkg.	50/cs	100/cs	40/cs



DYNJAABV1
DYNJAABV100



DYNJAAHF

Used for both
Respiratory
and Anesthesia
applications

Bypass Heat and Moisture Exchanger

The AARC recommends that ventilator circuits not be routinely changed and that passive humidifiers can be used for up to 48 hours.¹ It's widely understood that the opening of a respiratory circuit increases the risk of hospital-acquired infection. Therefore, delivering aerosolized medication creates a challenge in maintaining a closed circuit.

The Medline Bypass Heat and Moisture Exchanger supports AARC clinical recommendations by eliminating the need to open the circuit during aerosol medication delivery for up to 72 hours.

Features

- Easy and intuitive to use
- Clear marking of “HME” and “Aerosol” modes
- Use for up to 72 hours

Performance Specifications

- Tidal volume 250–1500 mL
- Moisture output 32.09 mg H₂O/L @ 500 mL TV at 72 hours
- Resistance to flow 0.6 cm H₂O @ 30 L/min flow
- Dead space 29.2 mL
- Weight 39.3 g

Item No.	Description	Pkg.
DYNJAABHME	15 mm ID/22 mm OD x 22 mm ID	25/cs



1. AARC Evidence-Based Clinical Practice Guidelines, Care of the Ventilator Circuit And Its Relation to Ventilator-Associated Pneumonia. Respiratory Care, September, 2003, Volume 48, No 9.

Heat and Moisture Exchangers and Filtered Heat and Moisture Exchangers

Our broad portfolio of HMEs and HMEFs meet your patient's humidification needs throughout the care continuum.

- Multiple tidal volume ranges and sizes to meet your clinical needs
- Non-filtered, electrostatic-filtered, and HEPA-filtered options available
- Angled options available to reduce the need for additional adapters and dead space
- Translucent housing for viewing of possible obstructions



	DYNJAAHME36	DYNJAAHME10	DYNJAAHME12	DYNJAAHME14	DYNJAPHME20
Description	HME, Adult	HME, Compact	HME, Compact	HME, Compact	HME, Pediatric
Type	Straight	Straight	Straight	Angled	Straight
Connection	22OD/15ID x 22ID/15OD	22OD/15ID x 22ID/15OD	22OD/15ID x 22ID/15OD	22OD/15ID x 22ID/15OD	22OD/15ID x 22ID/15OD
Vt (mL)	200-1500	120-1500	120-1500	120-1500	90-250
Vd (mL)	90	32.1	32.1	34.5	31.3
Weight (g)	22	14.7	14.7	16.8	13.3
CO₂ Port	Yes	Yes	No	Yes	Yes
Resistance (cm H₂O)	0.71 @ 30 LPM	0.6 @ 30 LPM	0.6 @ 30 LPM	0.5 @ 30 LPM	0.4 @ 30 LPM
Moisture Output (Mg H₂O/L)*	30 @ 500 mL	33.7 @ 500 mL	33.7 @ 500 mL	33.2 @ 500 mL	34.0 @ 500 mL
Pkg.	50/cs	50/cs	50/cs	50/cs	50/cs
HME Media	Foam	Foam	Foam	Foam	Foam



	DYNJAPHME38	DYNJAPHME34	DYNJAAHME30	DYNJAAHME40	DYNJAAHMEFLX
Description	HME, Infant	HME, Neonatal	HME, Trach	HME, Trach with Suction and Oxygen Ports	HME, Flex Tube
Type	Straight	Straight	Trach	Trach	Straight
Connection	N/A	15ID x 15OD	15 mm ID	15 mm ID	15mm ID x 22mm ID
Vt (mL)	75-500	15 - 50	> 25	>50	250-1500
Vd (mL)	11.5	3.1	14	13.4	90
Weight (g)	13	4.1	3.5	7.1	15.5
CO₂ Port	Yes	No	No	No	No
Resistance (cm H₂O)	2 @ 15 LPM	1.05 @ 10 LPM	0.51 @ 30 LPM	0.42 @ 30 LPM	1.4 @ 60 LPM
Moisture Output (Mg H₂O/L)*	28.5 @ 250 mL	24.8 @ 50 mL	25.9 @ 500 mL	26 @ 500 mL	30 @ 1 L
Pkg.	50/cs	50/cs	50/cs	50/cs	25/cs
HME Media	Paper	Paper	Paper	Paper	Foam


 Filtered HMEs



	DYNJAAHME1	DYNJAAHME5	DYNJAAHME9	DYNJAAHME11
Description	Filtered HME, Large	Filtered HME, Large	Filtered HME, Compact	Filtered HME, Compact
Type	Straight	Straight	Straight	Straight
Connection	22 M / 15 F x 22 F / 15 M	22 M / 15 F x 22 F / 15 M	22 M / 15 F x 22 F / 15 M	22 M / 15 F x 22 F / 15 M
Vt (mL)	150–1500	150–1500	120–1500	120–1500
Vd (mL)	40.7	40.7	31.3	31.3
Weight (g)	18.5	18.5	15	15
CO₂ Port	Yes	No	Yes	No
Resistance (cm H₂O)	0.8 @ 30 LPM	0.8 @ 30 LPM	2.1 @ 30 LPM	2.1 @ 30 LPM
Moisture Output (Mg H₂O/L)*	34.1 @ 500 mL	34.1 @ 500 mL	33.6 @ 500 mL	33.6 @ 500 mL
Filter Media	Electrostatic	Electrostatic	Electrostatic	Electrostatic
BFE(%)	99.999	99.999	99.999	99.999
VFE(%)	99.99	99.99	99.99	99.99
Pkg.	50/cs	50/cs	50/cs	50/cs
HME Media	Foam	Foam	Foam	Foam




	DYNJAAHME13	DYNJAPHME17	DYNJAAHMEHEPA	DYNJAAHME25
Description	Filtered HME, Compact	Filtered HME, Pediatric	Filtered HME, HEPA	Filtered HME, Flat, Adult
Type	Angled	Straight	Straight	Straight
Connection	22 M / 15 F x 22 F / 15 M	22 M / 15 F x 22 F / 15 M	22 M / 15 F x 22 F / 15 M	22 ID x 22 OD / 15 ID
Vt (mL)	120–1500	90–1500	200–1500	200–1500
Vd (mL)	33.7	20.5	52	62
Weight (g)	17.2	13.6	40	37.2
CO₂ Port	Yes	Yes	Yes	Yes
Resistance (cm H₂O)	1.9 @ 30 LPM	2.4 @ 30 LPM	1.5 @ 30 LPM	<1 @ 30 LPM
Moisture Output (Mg H₂O/L)*	33.4 @ 500 mL	34.0 @ 250 mL	30 @ 500 mL	32.5 @ 500 mL
Filter Media	Electrostatic	Electrostatic	HEPA	Electrostatic
BFE(%)	99.999	99.999	99.9999	>99.99%
VFE(%)	99.99	99.99	99.9999	>99.99%
Pkg.	50/cs	50/cs	50/cs	50/cs
HME Media	Foam	Foam	Foam	Foam

Non-Heated Breathing Circuits


Hudson RCI® Adult Conventional (non-heated) ventilator circuits are available in a variety of different configurations and styles.




HUD1607




HUD1617 *Includes two water traps




HUD1609



HUD1619



HUD1613



HUD69040

Adult Non-Heated Conventional Circuits

Circuit Configurations


Item No.	Limb	Tubing	IL-TRAP	EL-TRAP	H-Limb	E-Valve	Length	PMT	RTP	Notes
HUD1607	Dual	22 mm					72"			HME only
HUD1608	Dual	22 mm			•		72"			
HUD1609	Dual	22 mm	•	•	•		72"			
HUD1613	Dual	22 mm					72"			
HUD1614	Dual	22 mm	•	•			48"			
HUD1615	Dual	22 mm	•	•			60"			
HUD1617	Dual	22 mm	•	•			72"	•		Circuit with 1/8" and 1/4" ID pressure
HUD1619	Dual	22 mm					72"	•		
HUD1627	Dual	22 mm				•	72"	•		
HUD69030	Dual	22 mm					72"	•		Made of Kraton
HUD69040	Dual	22 mm					72"			Made of Kraton
DYNJ300	Dual	22 mm					72"			

Circuit Matrix Key:

H-IL = Heated Inspiratory Limb, H-EL = Heated Expiratory Limb, IL-TRAP = Inspiratory Limb Water Trap, EL-Trap = Expiratory Limb Water Trap, RTP = Remote Temperature Monitoring Port, PMT = Pressure Monitoring Tubing, H-LIMB = Humidifier Limb

Adult Non-Heated Single Limb NIV Circuits


Item No.	Description	Pkg.
HUD1698	Single Limb, 6 ft.	10/ea
HUD1698F	Single Limb, 6 ft., with Bacterial Viral Filter	10/ea
DYNJ500EP	Single Limb, Smooth Bore Circuit with Exhalation Port, 6 ft.	10/ea
DYNJ500EPF	Single Limb, Smooth Bore Circuit with Exhalation Port and Filter, 6 ft.	10/ea



DYNJ500EPF

Adult Portable & IPPB Ventilator Circuits

Item No.	Description	Pkg.
HUD1696	Adult Portable Ventilator Circuit, Includes Main Circuit Hose, Exhalation Valve Manifold, Aerosol Hose, Patient Elbow Connector, Proximal Airway Pressure Line, Exhalation Valve Line, and Humidifier Limb	10/cs
HUD1560	Universal IPPB Circuit, Complete with High-Output Nebulizer with 6 MI Jar Capacity, Mouthpiece, 6" and 39" Corr-A-Flex® Tubing, Multiple Adapter Kit, and Exhalation Valve And Nebulizer Drivelines	50/cs
HUD1667	Ventilator Monitoring Adapter Unit, Includes 6" tube, Tee, Nipple Adapter, 15 mm/22 mm Adapter, 2 One-Way Valves	50/cs
HUD1668	Trach Tee Oxygenator, Includes Oxygen Nipple Adapter, Aerosol Tee, and 6" Corr-A-Flex Tubing	50/cs



HUD1696

MRI-Conditional Transport Ventilator Circuits

- Ready-to-use transport ventilator circuits include expandable tubing and a MRI-conditional PEEP Valve
- Tubing is expandable up to 120" in length
- Patients can be scanned with the device in operation under the following conditions: static magnetic field up to 3.0 T; maximum spatial gradient field less than or equal to 30 T/m


Item No.	Description	Pkg.
DYNJAA300	Transport Ventilator Circuit with Bacterial Filter	20/cs
DYNJAA300NF	Transport Ventilator Circuit	20/cs




DYNJAA300

Pediatric Non-Heated Conventional Circuits

- Pediatric ventilator circuits are available in a variety of configurations and styles
- Circuits use 15 mm tubing, include two pressure lines, and incorporate standard 15 mm and 22 mm connectors
- Ported wyes allow pressure and temperature monitoring and include tethered caps
- All pediatric, heated wire and conventional circuits are 60" long



HUD78021




HUD78022

Circuit Configurations


Item No.	Limb	Tubing	IL-TRAP	EL-TRAP	H-Limb	E-Valve	Length	PMT	RTP	Notes	Pkg.
HUD78021	Dual	15 mm			•		60"	•			20/cs
HUD78022	Dual	15 mm	•	•	•		60"	•		Two separate water traps	20/cs
DYNJ200	Dual	15 mm					60"				20/cs

Neonatal Non-Heated Conventional Circuits


- Circuits use 10 mm tubing, include two pressure lines, and incorporate standard 15 mm and 22 mm connectors
- Ported wyes allow pressure and temperature monitoring and include tethered caps
- All neonatal circuits, heated-wire and conventional, are 48" long, unless noted
- All neonatal circuits include a humidifier limb
- Neonatal conventional and heated-wire circuits are available with and without water traps



HUD78010



HUD78011



HUD78012

Circuit Configurations

Item No.	Limb	Tubing	IL-TRAP	EL-TRAP	H-Limb	E-Valve	Length	PMT	RTP	Notes	Pkg.
HUD78010	Dual	10 mm	•	•	•		48"	•			20/cs
HUD78011	Dual	10 mm	•	•	•		48"	•	•	Two separate water traps	20/cs
HUD78012	Dual	10 mm	•	•	•		48"	•	•	Two separate water traps	20/cs



Non-Invasive Interfaces

Non-Invasive Ventilation Masks	20
Circuits & Accessories	21
CPAP Masks	22-24
Go-PAP™ Systems	25
O2-Max™ Systems	25-26
Babi Plus Infant Nasal CPAP	27

NIV Masks and Circuits

Flex™ Non-Invasive Ventilation (NIV) Adult Full Face Masks with FlexiPlane™ Technology

- Flex masks offer patient comfort and reduce the risk of skin breakdown
- Dual swivel, flexible elbow with silicone mask shell allows the mask to “flex” which aids to reduce pressure points.
- Large, adjustable silicone forehead pad designed to reduce pressure points and aid in stability
- Soft and flexible cushion reduces the need to overtighten the mask
- Comfortable headgear with clips for quick release

Item No.	Description	Pkg.
HCSFFVAS	Small, Vented, Anti-Asphyxia Valve (AAV) and Headgear	10/cs
HCSFFVAM	Medium, Vented, Anti-Asphyxia Valve (AAV) and Headgear	10/cs
HCSFFVAL	Large, Vented, Anti-Asphyxia Valve (AAV) and Headgear	10/cs
HCSFFNVAS	Small, Non-Vented, Anti-Asphyxia Valve (AAV) and Headgear	10/cs
HCSFFNVAM	Medium, Non-Vented, Anti-Asphyxia Valve (AAV) and Headgear	10/cs
HCSFFNVAL	Large, Non-Vented, Anti-Asphyxia Valve (AAV) and Headgear	10/cs
HCSFFNVS	Small, Non-Vented, with Headgear	10/cs
HCSFFNVM	Medium, Non-Vented, with Headgear	10/cs
HCSFFNVL	Large, Non-Vented, with Headgear	10/cs



HCSFFVAM

Non-Invasive Ventilation (NIV) Adult Nasal Masks

- Designed to improve patient comfort
- Vented and non-vented options
- Silicone forehead pad helps reduce pressure points and aids in stability
- Comfortable headgear with clips for quick release

Item No.	Description	Pkg.
HCSNVAS	Small, Vented, Anti-Asphyxia Valve (AAV) and Headgear	10/cs
HCSNVAM	Medium, Vented, Anti-Asphyxia Valve (AAV) and Headgear	10/cs
HCSNVAL	Large, Vented, Anti-Asphyxia Valve (AAV) and Headgear	10/cs
HCSNNVAS	Small, Non-Vented, Anti-Asphyxia Valve (AAV) and Headgear	10/cs
HCSNNVAM	Medium, Non-Vented, Anti-Asphyxia Valve (AAV) and Headgear	10/cs
HCSNNVAL	Large, Non-Vented, Anti-Asphyxia Valve (AAV) and Headgear	10/cs



HCSNVAL

Flex™ Non-Invasive Ventilation (NIV) adult mask accessories

- Replacement 4-point, breathable headgear for full face and nasal NIV masks
- Breathable material stretches to offer comfort and aid in mask cushion seal
- Includes clips for quick release




Item No.	Description	Pkg.
HCSFFHG	Headgear with Clips for Full Face NIV Masks	100/cs
HCSNHG	Headgear with Clips for Nasal NIV Masks	100/cs


Non-Invasive Ventilation (NIV) Circuits

- Flexible non-invasive circuits are compatible with vented and non-vented masks with anti-asphyxia valve

Item No.	Description	Pkg.
DYNJ500EP	Single Limb, Smooth Bore Circuit with Exhalation Port, 6 ft.	10/ea
DYNJ500EPF	Single Limb, Smooth Bore Circuit with Exhalation Port and Filter, 6 ft.	10/ea
DYNJ400EPF	Single Limb, Smooth Bore Circuit with Exhalation Port and Filter, 6 ft. without Pressure Line	10/ea
HUD1698	Universal Vent Circuit, Nonheated, NIV	10/ea
HUD1698F	Adult Single-Limb NPPV Breathing Circuit with Bacterial Viral Filter, 72"	10/ea
HUD87098KIT	ConchaSmart Circuit, Adult NIV	10/ea
HUD87098KITF	ConchaSmart Circuit, Adult NIV with Filter	10/ea
HUD88019KIT	ISO-Gard 70" Single-Limb Adult Circuit with Heated Inspiratory Limb and Limb Drain	10/ea
HUD88098KIT	ISO-Gard Single-Limb Adult Breathing Circuit with Heated Inspiratory Limb and Limb Drain	10/ea



DYNJ500EPF



HUD87098KITF

CPAP Masks & Accessories

AIRgel® CPAP Masks

- Amazing comfort and customized fit technology with patented AIRgel replaceable cushioning that minimizes discomfort and skin irritation while effectively sealing airflow
- Easily detached and reattached for easy cleaning, saving time and money
- Proprietary AIRgel cushioning minimizes discomfort and skin irritation while effectively sealing airflow
- The flexible and formable mask provides a custom fit

Aura Nasal Masks

- The shapeable shell of each mask can be easily molded to fit individual needs
- Flexible spacebar is designed to provide superior stability without touching the skin
- Active headgear connector moves with the patient to help maintain an effective seal

Item No.	Description	Mask Size	Pkg.
CPAP50176	Adjustable Gel Cushion	Small	1 ea
CPAP50177	Adjustable Gel Cushion	Medium	1 ea
CPAP50178	Adjustable Gel Cushion	Large	1 ea

Mojo® Full-Face Gel Masks

- Full-face design in four sizes
- Adjustable tension wheel minimizes discomfort while ergonomic cushions adapt to unique facial features for a comfortable fit
- Magnetic headgear clips hold firmly and disconnects with one hand easy on-and-off

Item No.	Description	Headgear Size	Mask Size	Pkg.
CPAP50835	Vented	One Size	Small	1 ea
CPAP50836	Vented	One Size	Medium	1 ea
CPAP50837	Vented	One Size	Large	1 ea
CPAP50838	Vented	One Size	X-Large	1 ea
CPAP55045	Headgear Only	One Size	One Size	1 ea



CPAP50177



CPAP50836

CPAP50838



Headgear included

CPAP50177

AIRgel® CPAP Masks

iQ Blue® Nasal Gel Mask

- One-size nasal mask design fits a wide variety of patients
- Clearest field of vision possible with no forehead pad required—patients can keep their glasses on while reading or watching TV
- Easy on-and-off thanks to quick release, breathable headgear

Item No.	Mask Style	Headgear Size	Mask Size	Pkg.
CPAP50655	Vented	Stable Fit	One Size	1 ea
CPAP55044	Headgear Only	Stable Fit	One Size	1 ea



CPAP50655

MiniMe® 2 Pediatric Gel Masks

- Nasal mask in four sizes for children ages 2 to 12
- Headgear with 5 adjustment points
- Soft, breathable headgear for all-night comfort
- Available in vented and non-vented designs
- Adjustable cushion with memory wire for a true custom fit

Item No.	Description	Headgear Size	Mask Size	Pkg.
CPAP60216	Vented	X-Small/Small	X-Small	10/cs
CPAP602172	Vented	X-Small/Small	Small	10/cs
CPAP602182	Vented	Medium/Large	Medium	10/cs
CPAP60221	Vented	Medium/Large	Large	10/cs
CPAP60256	Non-vented	X-Small/Small	X-Small	10/cs
CPAP602572	Non-vented	X-Small/Small	Small	10/cs
CPAP602582	Non-vented	Medium/Large	Medium	10/cs
CPAP60261	Non-vented	Medium/Large	Large	10/cs
CPAP55055	Headgear	X-Small/Small	N/A	10/cs
CPAP55056	Headgear	Medium/Large	N/A	10/cs



CPAP60217

Ascend® Masks

- Patented AIRgel replaceable cushion is easily detached and reattached for easy cleaning and replacement
- Leak ports are built-in, so a separate exhalation port is not needed
- Minimizes pressure points and skin irritation
- Ergonomic cushion technology delivers an effective seal for the most comfortable fit; flexible spacebar with unique bendable design holds a seal without contacting skin
- 360° rotating swivel for optimal tube positioning and easier mobility



CPAP50815

CPAP Full-Face Masks

Item No.	Description	Pkg.
CPAP50825	Starter Kit with Small, Medium, Large Cushions and Headgear	1 ea
CPAP50813	Full-Face Mask, Small	1 ea
CPAP50814	Full-Face Mask, Medium	1 ea
CPAP50815	Full-Face Mask, Large	1 ea
CPAP50207	Replacement Cushion, Small	1 ea
CPAP50208	Replacement Cushion, Medium	1 ea
CPAP50209	Replacement Cushion, Large	1 ea
CPAP55051	CPAP Headgear Replacement, One Size	1 ea

CPAP Nasal Masks

Item No.	Description	Pkg.
CPAP50174	Starter Kit with Small, Medium, Large Cushions and Headgear	1 ea
CPAP50187	Nasal Mask, Small	1 ea
CPAP50188	Nasal Mask, Medium	1 ea
CPAP50189	Nasal Mask, Large	1 ea
CPAP50201	Replacement Cushion, Small	1 ea
CPAP50202	Replacement Cushion, Medium	1 ea
CPAP50203	Replacement Cushion, Large	1 ea
CPAP55051	CPAP Headgear Replacement, One Size	1 ea

Silicone CPAP Masks

- A consistent seal without excessive tightness for prescribed therapy with enough comfort to help improve compliance
- OmniClip and silicone forehead pad for a custom fit at the bridge of the nose for patients of all head shapes and sizes
- SilentVent elbow flushes out exhaled air quietly, without an additional inline exhalation port
- Headgear attaches with convenient clips, includes a handy wash bag to keep it clean and fresh
- Cushions, forehead pads, elbows, cones and clips are individually replaceable

Item No.	Size	Description	Pkg.
CPAP9400	Small	Full Face Mask with Head Strap	1 ea
CPAP9401	Medium	• SilentVent Anti-Asphyxia Elbow, 22 mm Male	
CPAP9402	Large	• OmniClip with Silicone Forehead Pad	
CPAP9403	X-Large		
CPAP9404	Small	Full Face Mask, No Head Strap	1 ea
CPAP9405	Medium	• SilentVent Anti-Asphyxia Elbow, 22 mm Male	
CPAP9406	Large	• OmniClip with Silicone Forehead Pad	
CPAP9407	X-Large		
CPAP9408	Small	Nasal Mask with Head Strap	1 ea
CPAP9409	Medium	• SilentVent Elbow, 22 mm Male	
CPAP9410	Large	• Silicone Forehead Pad	
CPAP9411	X-Large	• OmniClip with Silicone Forehead Pad	
CPAP9412	Small	Nasal Mask, No Head Gear Strap	1 ea
CPAP9413	Medium	• SilentVent Elbow, 22 mm Male	
CPAP9414	Large	• Silicone Forehead Pad	
CPAP9415	X-Large	• Omniclip with Silicone Forehead Pad	

CPAP Tubing

Item No.	Description	Pkg.
CPAPTUB10	Standard CPAP/BiPAP Tubing, 6"	10/cs




CPAPTUB10



CPAP9408

CPAP9400



CPAP9408

Go-PAP™ CPAP Systems

- Offers emergency CPAP with integrated nebulization using the flow from an oxygen flowmeter or cylinder
- Connects to barbed connector and can generate up to 90 LPM flow
- Includes 3 options for PEEP of 5.0, 7.5 and 10.0 cm H₂O
- Integrated nebulization without interrupting the delivered flow to the patient



Item No.	Description	Pkg.
CPAP3134602	With BiTrac™ ED Mask and Standard Headgear, Adult Medium	10/cs
CPAP3134602N	With BiTrac™ ED Mask, Standard Headgear and Nebulizer Canister, Adult Medium	10/cs
CPAP3134602NC	With BiTrac™ ED Mask, Standard Headgear, Nebulizer Canister and Neb-Connect, Adult Medium	10/cs
CPAP3134603	With BiTrac ED Mask and Standard Headgear, Adult Large	10/cs
CPAP3134603N	With BiTrac ED Mask, Standard Headgear and Nebulizer Canister, Adult Large	10/cs
CPAP3134603NC	With BiTrac ED Mask, Standard Headgear, Nebulizer Canister and Neb-Connect, Adult Large	10/cs

O₂-Max™ CPAP Systems

- Offers emergency CPAP with integrated nebulization using the flow from an oxygen flowmeter or cylinder
- Connects to DISS port and can generate up to 140 LPM flow
- Includes 3 options for PEEP of 5.0, 7.5 and 10.0 cm H₂O
- Integrated nebulization without interrupting the delivered flow to the patient



O₂-MAX CPAP System with 3-Set CPAP (5.0, 7.5 and 10.0 cm H₂O)

Item No.	Description	Pkg.
CPAP3137555X	With BiTrac ED Mask, Adult Large	10/cs
CPAP3137555XN	With BiTrac ED Mask and Nebulizer Canister, Adult Large	10/cs
CPAP3137555X1	With BiTrac ED Mask and Ohmeda Quik-Connect™, Adult Large	10/cs
CPAP3137555N1	With BiTrac ED Mask, Ohmeda Quik-Connect and Nebulizer Canister, Adult Large	10/cs
CPAP3137555X2	With BiTrac ED Mask and Chemtron Quik-Connect, Adult Large	10/cs
CPAP3137556X	With BiTrac ED Mask, Adult Medium	10/cs
CPAP3137556XN	With BiTrac ED Mask and Nebulizer Canister, Adult Medium	10/cs
CPAP3137556X1	With BiTrac ED Mask, Ohmeda Quik-Connect and Nebulizer Canister, Adult Medium	10/cs
CPAP31375561	With BiTrac ED Mask and Ohmeda Quik-Connect, Adult Medium	10/cs
CPAP3137556X2	With BiTrac ED Mask and Chemtron Quik-Connect, Adult Medium	10/cs

O₂-MAX CPAP System with 5-Set CPAP (5.0, 7.5, 10.0, 12.5 and 15 cm H₂O)

Item No.	Description	Pkg.
CPAP3137554X	With BiTrac ED Mask, Adult Large	10/cs
CPAP3137554X1	With BiTrac ED Mask, Ohmeda Quik-Connect and Nebulizer Canister, Adult Large	10/cs
CPAP3137554XN	With BiTrac ED Mask and Nebulizer Canister, Adult Large	10/cs
CPAP3137557X	With BiTrac ED Mask, Adult Medium	10/cs
CPAP3137557X1	With BiTrac ED Mask, Ohmeda Quik-Connect and Nebulizer Canister, Adult Medium	10/cs
CPAP3137557XN	With BiTrac ED Mask and Nebulizer Canister, Adult Medium	10/cs

O₂-MaxCPAP systems continues on the next page

O₂-Max™ CPAP Systems (continued)

O₂-MAX CPAP Systems

(With fixed 10 cm H₂O CPAP)

Item No.	Description	Pkg.
CPAP3137535X	With BiTrac ED Mask, Adult Large	10/cs
CPAP3137538X	With BiTrac ED Mask, Adult Medium	10/cs
CPAP3137538XN	With BiTrac ED Mask and Nebulizer Canister, Adult Medium	10/cs
CPAP31375381	With BiTrac ED Mask and Ohmeda Quik-Connect, Adult Medium	10/cs
CPAP3137538X1	With BiTrac ED Mask, Ohmeda Quik-Connect and Nebulizer Canister, Adult Medium	10/cs

O₂-MAX Trio CPAP Systems

(Can be set to 30%, 60% or 90% Fio₂)

Item No.	Description	Pkg.
CPAP3138007X	3-SET PEEP System with BiTrac ED Mask with Head Strap, Adult Medium	10/cs
CPAP3138017X	3-SET PEEP System with BiTrac ED Mask with Head Strap, Adult Large	10/cs
CPAP3138007X1	3-SET PEEP System with BiTrac ED Mask with Head Strap Ohmeda Quik-Connect and Nebulizer Canister, Adult Medium	10/cs
CPAP3138107X	3-SET PEEP System with BiTrac ED Mask with Full Cap Headgear, Adult Medium	10/cs
CPAP3138017X1	3-SET PEEP System with BiTrac ED Mask with Head Strap, Ohmeda Quik-Connect and Nebulizer Canister, Adult Large	10/cs
CPAP3138117X	3-SET PEEP System with BiTrac ED Mask with Full Cap Headgear, Adult Large	10/cs
CPAP3138117XN	3-SET PEEP System with BiTrac ED Mask with Full Cap Headgear, Small Volume Nebulizer, Adult Large	10/cs
CPAP31381171	3-SET PEEP System with BiTrac ED Mask with Full Cap Headgear, Adult Large	10/cs
CPAP3138008X	5-SET PEEP System with BiTrac ED Mask with Head Strap, Adult Medium	10/cs
CPAP3138108X	5-SET PEEP System with BiTrac ED Mask with Full Cap Headgear, Adult Medium	10/cs
CPAP3138018X	5-SET PEEP System with BiTrac ED Mask with Head Strap, Adult Large	10/cs
CPAP3138118X	5-SET PEEP System with BiTrac ED Mask with Full Cap Headgear, Adult Large	10/cs
CPAP3138003X	10 cm H ₂ O fixed PEEP System with BiTrac ED Mask with Head Strap, Adult Medium	10/cs
CPAP3138013X	10 cm H ₂ O fixed PEEP System with BiTrac ED Mask with Head Strap, Adult Large	10/cs

O₂-Max Accessories

Item No.	Description	Pkg.
CPAP3137028	Small Mask for Go-PAP™ and O ₂ -Max CPAP Systems	10/cs
CPAP3137029	Medium Mask for Go-PAP and O ₂ -Max CPAP Systems	10/cs
CPAP3137030	Large Mask for Go-PAP and O ₂ -Max CPAP Systems	10/cs
CPAP313111111	Quik-Connect™ Ohmeda DISS	10/cs
CPAP3138000	O ₂ -Max Trio with 6" Tubing and DISS	10/cs

Infant Nasal CPAP

Hudson RCI Babi.Plus® Bubble Pap Valve

- 360° swivel inlet connector accepts 15 mm OD and 22 mm ID expiratory circuit connectors
- Positive engagement is designed to reduce the risk of inadvertent change
- Water level adjustment port allows water to be added or removed without disconnecting the expiratory circuit
- Gas discharge ports direct gas away from clinician
- Rotating pressure adjustment lever

Item No.	Description	Pkg.
HUD2690	Babi.Plus Bubble PAP Valve 10 cm H ₂ O	12/cs
HUD2691	Babi.Plus Pressure Relief Manifold	12/cs
HUD2692	Babi.Plus Single Pole Clamp	1 ea
HUD2693	Babi.Plus Dual Pole Clamp	1 ea



HUD2690

Infant Nasal Prong CPAP Cannulas

- Designed to reduce trauma associated with the delivery of infant nasal CPAP
- Soft, anatomically curved prongs enhance fit
- Luer fitting on expiratory connector allows proximal airway pressure monitoring
- 6 prong sizes allow greater choice for appropriate sizing for each infant
- Available in a system or in a set

Item No.	Description	Infant Weight Range	Suggested Cannula Size	French Size	Pkg.
HUD1683	System	<700 g	0	10	10/cs
HUD1690	Set	<700 g	0	10	10/cs
HUD1685	System	700–1,250 g	1	12	10/cs
HUD1691	Set	700–1,250 g	1	12	10/cs
HUD1686	System	1,250–2000 g	2	13	10/cs
HUD1692	Set	1,250–2000 g	2	13	10/cs
HUD1687	System	2,000–3,000 g	3	15	10/cs
HUD1693	Set	2,000–3,000 g	3	15	10/cs
HUD1688	System	>3,000 g	4	16	10/cs
HUD1694	Set	>3,000 g	4	16	10/cs
HUD1689	System	1–2 years	5	20	10/cs
HUD1695	Set	1–2 years	5	20	10/cs



HUD1687



HUD1690

CPAP Systems

Each system includes the following:

- Cannula
- Inspiratory elbow connector
- Expiratory elbow connector
- Knit cap
- Two 6" hook and loop fastener sections
- Two 4 ft., 10 mm ID corrugated tubing sections
- 4 ft. pressure monitoring line
- 22 to 10 mm humidifier adapter

CPAP Sets

Each system includes the following:

- Cannula
- Inspiratory elbow connector
- Expiratory elbow connector
- Knit cap
- Two 6" hook and loop fastener sections
- Two 10 to 7.5 mm adapters



Heated High Flow

Comfort Flo® Plus Nasal Cannula	29
High Flow Humidification Systems	30
MaxVenturi High Flow	31
MicroMax™ High Flow	32
MicroMax™ Low Flow	33
Accessories	34

High Flow Interfaces

Comfort Flo® Plus High Flow Nasal Cannulas

- Capable of delivering heated and humidified gas between 1–60 LPM
- Chin strap designed to encourage closed-mouth treatments, which may increase PEP¹
- Wide lanyard is designed to help distribute the weight of the circuit
- Adjustable bifurcated head strap offers a secure fit during patient movement

Item No.	Description	Pkg.
HUD241213	Comfort Flo Plus, Small	20/cs
HUD241212	Comfort Flo Plus, Medium	20/cs
HUD241211	Comfort Flo Plus, Large	20/cs
HUD241210	Comfort Flo Plus Chin Strap	50/cs
HUD241101	Nasal Cannula—Adult (Flows up to 60 LPM)	25/cs
HUD241102	Nasal Cannula—Pediatric (Flows up to 20 LPM)	25/cs
HUD241103	Nasal Cannula—Infant (Flows up to 8 LPM)	25/cs
HUD241104	Nasal Cannula—Premature (Flows up to 8 LPM)	25/cs



HUD241101



Large Bore Nares
Soft nasal prongs are designed to enhance patient comfort

Tubing
Flexible to accommodate head movement, rigid to resist kinking

Bifurcated Head Strap
Designed to encourage a secure fit during patient movement

Chin Strap (Sold Separately)
Can assist patients to keep mouth closed, which may increase PEP¹

HUD241211

Trach-Flex Plus Tracheostomy Adapter

- The Hudson RCI Trach-Flex Plus Tracheostomy Adapter facilitates oxygen and humidification delivery for your tracheostomy patients
- Features a transparent cap deflector that helps protect the clinician
- Adjustable bifurcated head strap offers a secure fit during patient movement

Item No.	Description	Pkg.
HUD241501	Tracheostomy Adapter	20/cs




HUD241501

1. Hirst KR, Patel A, Vines DL. Evaluation of Bronchial Pressures and Tidal Volume Using Three Different Adult High Flow Nasal Cannula (HFNC) Devices. AARC 2011 Open Forum # 13 Presentation. i. Cited in support of PEP effect only

High Flow Humidification Systems

Comfort Flo® Humidification Systems

- Delivers high flow nasal cannula therapy (HFNCT) for enhanced patient comfort.¹
- 1–60 LPM flow rates are available with the Neptune Heated Humidifier with Comfort Flo Column
- Provides continuous heated high flow therapy to a variety of patients ranging from neonatal to adult
- The Neptune Heated Humidifier is able to provide inspired gases at a temperature range of 28°–37° C to maximize patient comfort
- Temperature gradient selection for humidification control



HUD2415

Item No.	Description	Pkg.
HUD2410	Comfort Flo Humidification System*	20/cs
HUD2414	Comfort Flo Humidification System with Remote Temperature Port†	20/cs
HUD2415	Corrugated Comfort Flo Humidification System*	20/cs
HUD2416	Corrugated Comfort Flo Humidification System with Remote Temperature Port†	20/cs
HUD2417	Comfort Flo 15 mm Oxygen Tubing Adapter	50/cs

*System includes circuit, Comfort Flo column, supply tubing, 5-PSI system pressure relief valve and tubing adapter. Water and nasal cannula sold separately.

†System includes circuit with remote port, Comfort Flo column, supply tubing, 5-PSI system pressure relief valve and tubing adapter. Water and nasal cannula sold separately.

1. Restrepo RD, Walsh BK. AARC CPG: Humidification during invasive and noninvasive mechanical ventilation. Respiratory Care. 2012; 57 (5): 782-788.

MaxVenturi® Blender Setup

- Only requires piped oxygen
- Uses the venturi principle to mix ambient air and oxygen without the need for piped medical grade air
- 15' Single Coiled O2 Hose (HUDR127P35)
- Overall flow range of 0-60 LPM*

Item No.	Description	Pkg.
HUDR211P03036	MaxVenturi® Blender	1 ea
HUDR127P35	15' Single Coiled O2 Hose	1 ea
HUDR211P38	MaxVenturi® Muffler	1 ea



Heated High Flow



HUDR211P03036



HUDR127P35



HUDR211P38

MicroMax™ High Flow (DISS Connection)

- Requires oxygen and medical air connection
- Features delivery of accurate FiO2 mixtures from two outlet ports and has a primary outlet flow range of 2-120 LPM*
- 10' Dual Blender Hose (HUDR129P02)

Item No.	Description	Pkg.
HUDR203P14	MicroMax High Flow Blender	1 ea



Blender Buddy for MicroMax™ 0-70 LPM

- Robust design provides quality and reliability
- Manifold for use with the MicroMax™ Blender that allows for the attachment of a flow meter and oxygen analyzer
- Oxygen analyzer port on the top of the manifold allows gas to be analyzed in a clean and dry environment away from the patient to reduce the risk of cross-contamination
- Includes one high flow (0-70 LPM)* flow meter and one MaxO2+A Analyzer

Item No.	Description	Pkg.
HUDR219P20002	Blender Buddy for MicroMax™ 0-70 LPM	1 ea
HUDR219P20003	Blender Buddy for MicroMax™ 0-70 LPM with MaxO2+A Analyzer	1 ea
HUDR100P26	Blender Pole Mount	1 ea
HUDR129P01	Blender Hose, Dual-Coil, Diss, 10'	1 ea



HUDR219P20003



HUDR129P01



HUDR100P26

MicroMax™ Low Flow

- Requires piped oxygen and medical air connection
- Features delivery of accurate FiO2 mixtures from two outlet ports and has a primary outlet flow range of 0–30 LPM*



Item No.	Description	Pkg.
HUDR203P13	MicroMax Low Flow Blender	1 ea

Blender Buddy MicroMax™ 0–3.5 LPM and 0–15 LPM with MaxO2+A

- Robust design provides quality and reliability
- Manifold for use with the MicroMax™ Blender that allows for the attachment of a flow meter and oxygen analyzer
- Oxygen analyzer port on the top of the manifold allows gas to be analyzed in a clean and dry environment away from the patient to reduce the risk of cross-contamination
- Includes two low flow (0-3 LPM and 0-15 LPM)* flow meters and one MaxO2+A Analyzer



HUDR219P20004

Item No.	Description	Pkg.
HUDR219P20004	Blender Buddy MicroMax™ 0-3.5 LPM and 0-15 LPM with MaxO2+A	1 ea
HUDR100P26	Blender Pole Mount	1 ea
HUDR129P01	Blender Hose, Dual-Coil, Diss, 10'	1 ea



HUDR100P26



HUDR129P01

Accessories

MaxO2+A Oxygen Analyzer

- Designed with simplicity and portability in mind
- Ideal for checking oxygen percentage accuracy in concentrators

Item No.	Description	Pkg.
HUDR217P62	MaxO2+A Oxygen Analyzer	1 ea



SmartStack Pole Stand

- Specially designed 24" pole base provides extra stand support and durability
- Patented design stacks neatly with other bases of its kind (22" and 24"), allowing for extra storage space and easy transportation of multiple stands all at once

Item No.	Description	Pkg.
HUDR100P49001	SmartStack Pole Stand	1 ea



HUDR100P49001

IV Pole Basket (1")

- Fits any IV pole with a 1" diameter and includes a robust, solid aluminum mounting bracket that securely attaches the basket to the pole
- Protected by a durable, powder-coated finish

Item No.	Description	Pkg.
HUDR105P17	IV Pole Basket, 1"	1 ea



Dual Tank Holder

- Attaches to most IV pole stands

Item No.	Description	Pkg.
HUDR105P04001	Dual Tank Holder	1 ea





Resuscitation

Manual Resuscitators	36–39
Colorimetric CO2 Detectors	40
Resuscitation Masks & Accessories	41
Hyperinflation Systems	42

Manual Resuscitators

Easy-to-use manual resuscitators are available in infant, pediatric and adult sizes, and include a color-coded drawstring bag to indicate size.

A Better Bag

- Delivers up to 99% FiO₂
- Maneuverable—the body of the bag swivels 360° for multiple clinician hand positions, to help reduce clinician fatigue
- Translucent—lightweight bag body and components offer visualization of vomitus or occlusions in the system
- Textured—for an enhanced, non-slip grip
- Repositionable—patient port swivels 360°, allowing the clinician freedom to move around the patient without placing excessive torque or stress on the tracheal tube or mask
- Manometer allows for a clear visual reading of the centimeters of water pressure being delivered to the patient

An Improved Mask

- Premium material offers a soft cushion and a firm seal

	Adult (Mask Size #6)	Small Adult (Mask Size #6)	Pediatric (Mask Size #3)	Infant (Mask Size #2)
Bag Volume	1,900 mL	1,500 mL	1,000 mL	300 mL
Stroke Volume	950–960 mL	640–660 mL	600–620 mL	140–170 mL
Dead Space with Mask	10 mL	10 mL	10 mL	10 mL
Inspiratory Resistance	<5 cm H ₂ O	<5 cm H ₂ O	<5 cm H ₂ O	<5 cm H ₂ O
Exhalation Resistance	<5 cm H ₂ O 50 LPM flow	<5 cm H ₂ O 50 LPM flow	<5 cm H ₂ O 50 LPM flow	<5 cm H ₂ O 50 LPM flow
Attainable Delivery Pressure	≥80 cm H ₂ O	≥80 cm H ₂ O	≥80 cm H ₂ O	≥80 cm H ₂ O
Ventilatory Frequency	≥60 cm bpm (no load)	≥60 cm bpm (no load)	≥60 cm bpm (no load)	≥90 cm bpm (no load)
Body Mass Use Range	≥40 kg (88 lbs)	≥40 kg (88 lbs)	≤40 kg (88 lbs)	≤10 kg (22 lbs)
Minimal Tidal Volume	600 mL	600 mL	500 mL	150 mL



Also available with Manometer.

Manual Resuscitator Ordering Guide

Adult Resuscitation Bags

Item No.	Mask Size	Reservoir Style		Tubing Length		Pressure Relief Valve	CO ₂ Detector	CO ₂ FLOCAP	Peep Valve	Exhalation Filter	Manometer	Pkg.
		Bag	40" Tubing	7-ft.	14-ft.							
CPRM1116	6	•		•								6/cs
CPRM11163	4, 5, 6	•		•								6/cs
CPRM1116B	6	•		•								6/cs
CPRM1116C	6	•		•			•					6/cs
CPRM1116C2	5, 6	•		•			•					6/cs
CPRM1116F	6	•		•						•		6/cs
CPRM1116FC	6	•		•			•			•		6/cs
CPRM1116FM	6	•		•						•	•	6/cs
CPRM1116FMT	6	•			•					•	•	6/cs
CPRM1116FP	6	•		•					•	•		6/cs
CPRM1116FPC	6	•		•			•		•	•		6/cs
CPRM1116FPM	6	•		•					•	•	•	6/cs
CPRM1116L	6	•		•								6/cs
CPRM1116M	6	•		•							•	6/cs
CPRM1116MC	6	•		•			•				•	6/cs
CPRM1116P	6	•		•					•			6/cs
CPRM1116PC	6	•		•			•		•			6/cs
CPRM1116PCM	6	•		•			•		•		•	6/cs
CPRM1116PD5	6	•		•					•			6/cs
CPRM1116PM	6	•		•					•		•	6/cs
CPRM1116PMT	6	•			•				•		•	6/cs
CPRM1116PW	6	•		•					•			6/cs
CPRM1116T	6	•			•							6/cs
CPRM1116WM	6	•		•							•	6/cs
CPRM1126	6		•	•								6/cs
CPRM1126C	6		•	•			•					6/cs
CPRM1126F	6		•	•						•		6/cs
CPRM1126FC	6		•	•			•			•		6/cs
CPRM1126FL	6		•	•						•		6/cs
CPRM1126FM	6		•	•						•	•	6/cs
CPRM1126FP	6		•	•					•	•		6/cs
CPRM1126FPC	6		•	•			•		•	•		6/cs
CPRM1126FPM	6		•	•					•	•	•	6/cs
CPRM1126FPWM	6		•	•					•	•	•	6/cs
CPRM1126FPWT	6		•		•					•		6/cs
CPRM1126M	6		•	•							•	6/cs
CPRM1126P	6		•	•					•			6/cs
CPRM1126PC	6		•	•			•		•			6/cs

Small Adult Resuscitation Bags

Item No.	Mask Size	Reservoir Style		Tubing Length		Pressure Relief Valve	CO ₂ Detector	CO ₂ FLOCAP	Peep Valve	Exhalation Filter	Manometer	Pkg.
		Bag	40" Tubing	7-ft.	14-ft.							
CPRM4416	6	•		•								6/cs
CPRM4416M	6	•		•							•	6/cs
CPRM4416P	6	•		•					•			6/cs
CPRM4426	6		•	•								6/cs
CPRM4426M	6		•	•							•	6/cs

Manual Resuscitator Ordering Guide

Pediatric Resuscitation Bags

Item No.	Mask Size	Reservoir Style		Tubing Length		Pressure Relief Valve	CO ₂ Detector	CO ₂ FLOCAP	Peep Valve	Exhalation Filter	Manometer	Pkg.
		Bag	40" Tubing	7-ft.	14-ft.							
CPRM2216	3	●		●		40 cm H ₂ O						6/cs
CPRM2216C	3	●		●		40 cm H ₂ O	●					6/cs
CPRM2216F	3	●		●		40 cm H ₂ O				●		6/cs
CPRM2216FC	3	●		●		40 cm H ₂ O	●			●		6/cs
CPRM2216FM	3	●		●		40 cm H ₂ O				●	●	6/cs
CPRM2216FM3	2, 3, 4	●		●		40 cm H ₂ O				●	●	6/cs
CPRM2216FP	3	●		●		40 cm H ₂ O			●	●		6/cs
CPRM2216FPC	3	●		●		40 cm H ₂ O	●		●	●		6/cs
CPRM2216FPM	3	●		●		40 cm H ₂ O			●	●	●	6/cs
CPRM2216FPM3	2, 3, 4	●		●		40 cm H ₂ O			●	●	●	6/cs
CPRM2216M	3	●		●		40 cm H ₂ O					●	6/cs
CPRM2216M3	2, 3, 4	●		●		40 cm H ₂ O					●	6/cs
CPRM2216P	3	●		●		40 cm H ₂ O			●			6/cs
CPRM2216P3M3	2, 3, 4	●		●		40 cm H ₂ O			●		●	6/cs
CPRM2216PC	6	●		●		40 cm H ₂ O	●		●			6/cs
CPRM2216PCM	3	●		●		40 cm H ₂ O			●		●	6/cs
CPRM2216PM	3	●		●		40 cm H ₂ O			●		●	6/cs
CPRM2216PM6	3	●		●		40 cm H ₂ O			●		●	6/cs
CPRM2216PMT	3	●			●	40 cm H ₂ O			●		●	6/cs
CPRM22162	2, 3	●		●		40 cm H ₂ O						6/cs
CPRM22163	1, 2, 3	●		●		40 cm H ₂ O						6/cs
CPRM2226	3		●	●		40 cm H ₂ O						6/cs
CPRM2226B	3		●	●		40 cm H ₂ O						6/cs
CPRM2226C	3		●	●		40 cm H ₂ O	●					6/cs
CPRM2226F	3		●	●		40 cm H ₂ O				●		6/cs
CPRM2226F2	3, 4		●	●		40 cm H ₂ O				●		6/cs
CPRM2226FC	3, 4		●	●		40 cm H ₂ O	●			●		6/cs
CPRM2226FM	3		●	●		40 cm H ₂ O				●	●	6/cs
CPRM2226FM3	2, 3, 4		●	●		40 cm H ₂ O				●	●	6/cs
CPRM2226FP	3		●	●		40 cm H ₂ O			●	●		6/cs
CPRM2226FPC	3		●	●		40 cm H ₂ O	●		●	●	●	6/cs
CPRM2226FPM	3		●	●		40 cm H ₂ O			●	●	●	6/cs
CPRM2226FPWM	3		●	●		40 cm H ₂ O			●	●	●	6/cs
CPRM2226M	3		●	●		40 cm H ₂ O					●	6/cs
CPRM2226P	3		●	●		40 cm H ₂ O			●			6/cs
CPRM2226PC	3		●	●		40 cm H ₂ O	●		●			6/cs
CPRM2226PW	3		●	●		40 cm H ₂ O			●			6/cs

If you do not see your preferred configuration, Medline can customize a resuscitator bag to your preference.

Manual Resuscitator Ordering Guide

Infant Resuscitation Bags


Item No.	Mask Size	Reservoir Style		Tubing Length		Pressure Relief Valve	CO ₂ Detector	CO ₂ FLOCAP	Peep Valve	Exhalation Filter	Manometer	Pkg.
		Bag	40" Tubing	7-ft.	14-ft.							
CPRM3312	2	•		•		40 cm H ₂ O						12/cs
CPRM3312C	1, 2	•		•		40 cm H ₂ O	•					12/cs
CPRM3312F	1, 2, 3	•		•		40 cm H ₂ O				•		12/cs
CPRM3312FC	2	•		•		40 cm H ₂ O				•		12/cs
CPRM3312FM	2	•		•		40 cm H ₂ O				•	•	12/cs
CPRM3312FM123	2	•		•		40 cm H ₂ O				•	•	12/cs
CPRM3312FP	2	•		•		40 cm H ₂ O			•	•		12/cs
CPRM3312FPC	1, 2, 3	•		•		40 cm H ₂ O	•		•	•		12/cs
CPRM3312FPM	2	•		•		40 cm H ₂ O			•	•	•	12/cs
CPRM3312FPM3	2	•		•		40 cm H ₂ O			•	•	•	12/cs
CPRM3312M	2	•		•		40 cm H ₂ O					•	12/cs
CPRM3312M3	0, 1, 2	•		•		40 cm H ₂ O					•	12/cs
CPRM3312P	2	•		•		40 cm H ₂ O			•			12/cs
CPRM3312P3M3	1, 2, 3	•		•		40 cm H ₂ O			•		•	12/cs
CPRM3312PC	2	•		•		40 cm H ₂ O	•		•			12/cs
CPRM3312PM	0, 1, 2	•		•		40 cm H ₂ O			•		•	12/cs
CPRM3312PM1	2	•		•		40 cm H ₂ O			•		•	12/cs
CPRM3312PM2	2	•		•		40 cm H ₂ O			•		•	12/cs
CPRM3312PM3	2	•		•		40 cm H ₂ O			•		•	12/cs
CPRM3312PMC2	2	•		•		40 cm H ₂ O	•		•		•	12/cs
CPRM3312PMTD1	1, 2, 3	•			•	40 cm H ₂ O			•		•	12/cs
CPRM3312PMX	1, 2	•		•		40 cm H ₂ O			•		•	12/cs
CPRM3312PWM	1	•		•		40 cm H ₂ O			•		•	12/cs
CPRM33122		•		•		40 cm H ₂ O						12/cs
CPRM33123	2	•		•		40 cm H ₂ O						12/cs
CPRM3322	2		•	•		40 cm H ₂ O						12/cs
CPRM3322C	2		•	•		40 cm H ₂ O	•					12/cs
CPRM3322F	2		•	•		40 cm H ₂ O				•		12/cs
CPRM3322FC	2		•	•		40 cm H ₂ O	•			•	•	12/cs
CPRM3322FM	2		•	•		40 cm H ₂ O				•		12/cs
CPRM3322FP	2		•	•		40 cm H ₂ O			•	•		12/cs
CPRM3322FPC	2		•	•		40 cm H ₂ O	•		•	•		12/cs
CPRM3322M	2		•	•		40 cm H ₂ O					•	12/cs
CPRM3322MC	2		•	•		40 cm H ₂ O	•				•	12/cs
CPRM3322MG	2		•	•		40 cm H ₂ O					•	12/cs
CPRM3322P	2		•	•		40 cm H ₂ O			•			12/cs
CPRM3322PC	2		•	•		40 cm H ₂ O	•		•			12/cs
CPRM3322PM	2		•	•		40 cm H ₂ O			•		•	12/cs
CPRM3412P	0, 2	•		•		40 cm H ₂ O			•			12/cs

If you do not see your preferred configuration, Medline can customize a resuscitator bag to your preference.

Colorimetric CO₂ Detectors

- A simple way to indicate the presence of carbon dioxide during a patient's breathing cycle.
- Quick response of changes in color between Blue and Yellow make it easy to visualize and determine whether exhaled CO₂ is present
- Long lasting—can be used for up to 24 hours
- Both CO₂ detectors have a small internal dead space.
Item No. CPRM2400 is 16 mL and CPRM2401 is 3 mL
- Supports patients weighing 1 kg or greater

Item No.	Description	Pkg.
CPRM2400	Colorimetric CO ₂ Detector for Patients ≥10 kg (≥22-lbs.)	10/cs
CPRM2401	Colorimetric Pediatric CO ₂ Detector for Patients 1–15 kg (2.2–33-lbs.)	10/cs



Resuscitation Masks and Accessories

Resuscitation Masks

- An ergonomic cushion design for an optimal seal and clear cone to visualize patient breathing
- Removable hook ring for hand-held or head-strap application (hook rings included on sizes 3–6)
- No valve prevents mask deflation
- Available in all sizes, including preemie

Item No.	Description	Pkg.
CPRMASK11R	Neonatal, Size 1	40/cs
CPRMASK12R	Infant, Size 2	40/cs
CPRMASK13R	Child, Size 3	40/cs
CPRMASK14R	Small Adult, Size 4	40/cs
CPRMASK15R	Medium Adult, Size 5	20/cs
CPRMASK16R	Large Adult, Size 6	20/cs

Medline offers a wide variety of masks which are not shown here. Please contact us for more information.



PEEP Valves and Filters

- Work with all Medline resuscitation devices
- Customize to patient need

Item No.	Description	Pkg.
CPRMP19	5-20 cm H2O Adjustable PEEP Valve, 19 mm ID and 22 mm OD Connections	20/cs
CPRMP30	5-20 cm H2O Adjustable PEEP Valve, 19 mm ID, 22 mm OD, and 30 mm ID Adapter	20/cs
CPRMFILTER	Expiratory Filter	50/cs



CPRMP19



CPRMFILTER

Mouth-to-Mask Resuscitator with Valve

- Clear plastic design for easy placement with a highly effective, one-way filtering system
- Soft inflated cuff for better fit
- Standard 15mm fittings are compatible with any resuscitator
- Easy-to-understand illustrations are printed on hard plastic carrying case
- Standard mask size

Item No.	Description	Pkg.
HCS64185	Resuscitator Mask with One-Way Valve	12/cs



Medline offers a wide variety of masks which are not shown here. Please contact us for more information.



HCS64185

Hyperinflation Systems

- Help prevent or treat lung collapse
- Bag design uses a “dipped” process to enhance tactile sensitivity and improve control
- Uniform inflation and the highest level of performance
- 100% compliant with ISO 5362 standards
- Variety of sizes to more accurately fit every patient



Hyperinflation Systems

7-ft. (2.13 m) of O₂ tubing included*

Item No.	Breathing Bag	O ₂ Connector	Mask	Manometer	Pkg.
DYNJAAHY005	0.5 L	Standard	N/A		20/cs
DYNJAAHY600	0.5 L with Tail-end Bleed Valve	Universal	N/A		50/cs
DYNJAAHY30	0.5 L	Standard	1 Neonatal		20/cs
DYNJAAHY29	0.5 L	Universal	2 Infant (Pre-inflated)		20/cs
DYNJAAHY01	1 L	Standard	N/A		20/cs
DYNJAAHY601	1 L with Tail-end Bleed Valve	Universal	N/A		50/cs
DYNJAAHY02	2 L	Standard	N/A		20/cs
DYNJAAHY03	3 L	Standard	N/A		20/cs
DYNJAAHY03M5	3 L	Universal	5 Adult		20/cs
DYNJAAHY5005	0.5 L	Standard	N/A	●	50/cs
DYNJAAHY5053P	0.5 L	Universal	N/A	●	50/cs
DYNJAAHY32	0.5 L	Standard	1 Neonatal	●	50/cs
DYNJAAHY33	0.5 L	Standard	2 Infant	●	50/cs
DYNJAAHY35	0.5 L	Universal	2 Infant	●	50/cs
DYNJAAHY34	0.5 L	Universal	2 Infant (Pre-inflated)	●	50/cs
DYNJAAHY501	1 L	Universal	N/A	●	50/cs
DYNJAAHY502	2 L	Universal	N/A	●	50/cs
DYNJAAHY503	3 L	Universal	N/A	●	50/cs

*Custom options are available



DYNJAAHY601



DYNJAAHY5005



Airway Management

AG Cuffill	44
Closed Suction	45-47
Open Suction	48
Addipak Sterile Solutions	49
Trach Care & Cleaning Trays	50-52
Tracheostomy Tube Holders	53
SubG™ Endotracheal Tubes	54
Mucus/Specimen Traps	55

AG CUFFILL® Cuff Pressure Manometer

Medline
Exclusive

Precise airway cuff pressure is critical to minimizing complications associated with intubation. Studies show that conventional methods caregivers use for cuff pressure readings are unreliable.¹ In fact, 52% of intubated patients have been found to be out of therapeutic cuff pressure range, putting them at risk for VAE or tracheal damage.²

AG CUFFILL can help. The pocket-sized and easy-to-use AG CUFFILL device is an accurate solution for both controlling cuff volume and measuring cuff pressure for cuffed ET tubes, trach tubes and LMAs.

Simple: Easily measure and adjust cuff pressure and volume of cuffed airways

Accurate: Digital volume and pressure readout, +/-2 (cmH2O) accuracy with simple-one-handed calibration

Safe: Designed for single-patient use to help reduce the chances of cross-contamination risk

Accessible: Pocket-sized, light-weight and stays at the patient's bedside

Cost Effective: Priced to make regular cuff pressure checks and adjustments a new standard of care

Item No.	Description	Pkg.
HCSCUFF0041	AG Cuffill Device, 10 cc	10/bx
HCSCUFF0041H	AG Cuffill Device, 10 cc	1 ea



References: 1. Mueen Ullah Khan, et. al., Measurement of endotracheal tube cuff pressure: Instrumental versus conventional method. Saudi J Anaesth. Oct-Dec 2016;10(4):428-431. Accessed Sept 27, 2020. <https://pubmed.ncbi.nlm.nih.gov/27833487/> 2. Sole ML, Su X, Talbert S, Penoyer DA, Kalita S, Jimenez E, Ludy JE, Bennett M. Evaluation of an intervention to maintain endotracheal tube cuff pressure within therapeutic range. Am J Crit Care. 2011 Mar; 20(2):109-17. PMID: 21362715, <https://www.ncbi.nlm.nih.gov/pubmed/21362715>

ClearPro™ Closed Suction Catheters

The power to set—and follow—best practice.

ClearPro provides the performance you expect without a time-based change-out recommendation on the label. This gives you the freedom to set your own change-out policy to meet both clinical and regulatory recommendations—and drive down costs.

Plus, it includes advanced features to help provide better care:

- Isolated cleaning chamber for rapid irrigation
- Pressure wiper seal to help maintain PEEP

Additional details:


- Kitted with all necessary adapters and connectors
- Does not contain DEHP
- Available in French sizes 10–16 with catheter depth markings and color-coded by ISO standards

Endotracheal


Item No.	Description	Pkg.	Thumb Valve Color
DYNCPDS10	Double Swivel Elbow, 10 Fr	20/cs	Black
DYNCPDS12	Double Swivel Elbow, 12 Fr	20/cs	White
DYNCPDS14	Double Swivel Elbow, 14 Fr	20/cs	Green
DYNCPDS16	Double Swivel Elbow, 16 Fr	20/cs	Orange
DYNCPTP10	T-Piece Adapter, 10 Fr	20/cs	Black
DYNCPTP12	T-Piece Adapter, 12 Fr	20/cs	White
DYNCPTP14	T-Piece Adapter, 14 Fr	20/cs	Green
DYNCPTP16	T-Piece Adapter, 16 Fr	20/cs	Orange

Tracheostomy


Item No.	Description	Pkg.	Thumb Valve Color
DYNCPDS10T	Double Swivel Elbow, 10 Fr	20/cs	Black
DYNCPDS12T	Double Swivel Elbow, 12 Fr	20/cs	White
DYNCPDS14T	Double Swivel Elbow, 14 Fr	20/cs	Green
DYNCPDS16T	Double Swivel Elbow, 16 Fr	20/cs	Orange
DYNCPTP10T	T-Piece Adapter, 10 Fr	20/cs	Black
DYNCPTP12T	T-Piece Adapter, 12 Fr	20/cs	White
DYNCPTP14T	T-Piece Adapter, 14 Fr	20/cs	Green
DYNCPTP16T	T-Piece Adapter, 16 Fr	20/cs	Orange



Colored Thumb Valves



Accessories



Closed Suction Catheters

Adult

- Labeled for facilities to determine changeout policy
- Soft tip helps reduce damage to mucous membranes
- Transparent connector for easy observation
- High-quality, translucent EVA sleeve helps reduce cross-contamination and facilitates operation when advancing the catheter

Endotracheal

Item No.	Description	Pkg.	Thumb Valve Color
DYNCSDS12	Double Swivel Elbow, 12 Fr	10/cs	White
DYNCSDS14	Double Swivel Elbow, 14 Fr	10/cs	Green
DYNCSTP12	T-Piece Adapter, 12 Fr	10/cs	White
DYNCSTP14	T-Piece Adapter, 14 Fr	10/cs	Green

Tracheostomy

Item No.	Description	Pkg.	Thumb Valve Color
DYNCSDS12T	Double Swivel Elbow, 12 Fr	10/cs	White
DYNCSDS14T	Double Swivel Elbow, 14 Fr	10/cs	Green
DYNCSTP12T	T-Piece Adapter, 12 Fr	10/cs	White
DYNCSTP14T	T-Piece Adapter, 14 Fr	10/cs	Green



DYNCSDS12T

Also available with T-Piece adapter

Pediatric

- Colored depth markings for clear viewing in low-light settings
- Available with an elbow or Y-style modular patient adapter
- Three adapter sizes to ensure compatibility
- Labeled for facilities to determine change out policy
- Compatible with all modes of ventilation

Elbow Adapter

Item No.	Description	Pkg.	Thumb Valve Color
DYNCPE5	Elbow Adapter, 5 Fr	10/cs	Black
DYNCPE6	Elbow Adapter, 6 Fr	10/cs	Light Green
DYNCPE7	Elbow Adapter, 7 Fr	10/cs	Purple
DYNCPE8	Elbow Adapter, 8 Fr	10/cs	Light Blue

Y-Connector

Item No.	Description	Pkg.	Thumb Valve Color
DYNCPY5	Y-Connector, 5 Fr	10/cs	Black
DYNCPY6	Y-Connector, 6 Fr	10/cs	Light Green
DYNCPY7	Y-Connector, 7 Fr	10/cs	Purple
DYNCPY8	Y-Connector, 8 Fr	10/cs	Light Blue



DYNCYP6

Elbow Adapter and Y-Connector products come with three adapter sizes in each kit

Open Suction Catheters

- Our D419 series catheters offers a full range of straight packaging
- Available in a complete range of continuous adult and pediatric sizes
- Easy-to-read depth markings
- DeLee tip is flat, soft and smooth to minimize airway trauma for pediatric patients
- Whistle tip is angled, soft and smooth to maximize suctioning capability for adult patients, while minimizing trauma potential

Single Sterile-Straight Pack

Item No.	French Size	Pkg.
DYND41906	6 Fr	100/cs
DYND41908	8 Fr	100/cs
DYND41900	10 Fr	100/cs
DYND41901	12 Fr	100/cs
DYND41902	14 Fr	100/cs
DYND41903	16 Fr	100/cs
DYND41904	18 Fr	100/cs




DYND41902




DYND41904

Coiled

Item No.	French Size	Pkg.
DYND41906C	6 Fr	100/cs
DYND41908C	8 Fr	100/cs
DYND41900C	10 Fr	100/cs
DYND41901C	12 Fr	100/cs
DYND41902C	14 Fr	100/cs
DYND41903C	16 Fr	100/cs
DYND41904C	18 Fr	100/cs



DYND40992



DYND41962



DYND41902C

Poly-Cath

Item No.	French Size	Pkg.
DYND40990	10 Fr	100/cs
DYND40992	14 Fr	100/cs

Continuous Valve

Item No.	French Size	Pkg.
DYND41960	10 Fr	100/cs
DYND41962	14 Fr	100/cs
DYND41964	18 Fr	100/cs


Sleeved Suction Catheters

- D407 series open-suction catheters are contained in a sleeve to help avoid cross-contamination
- Also includes a pop-up cup

Item No.	French Size	Pkg.
DYND40708F	8 Fr	50/cs
DYND40700F	10 Fr	50/cs
DYND40701F	12 Fr	50/cs
DYND40702F	14 Fr	50/cs



DYND40708F



Addipak[®] Unit Dose Solutions

- Sterile solutions
- Tamper-proof, color-coded vials for easy solution identification
- Inverted mL graduations and precise opening for accurate and controlled dispensing of solution
- 15 mL vials packaged 48 to a box, 3 boxes to a case
- 3 mL and 5 mL vials packaged 100 to a box, 10 boxes to a case
- Not for use in parenteral applications

Item No.	Size	Sterile Water (Blue)	0.45% Saline (Green)	0.9% Saline (Red)	Pkg.
HUD20031	3 mL	•			1000/cs
HUD20034	3 mL		•		1000/cs
HUDRHU39U	3 mL			•	1000/cs
HUD20051	5 mL	•			1000/cs
HUDRHU59U	5 mL			•	1000/cs
HUDRHU69U	15 mL			•	144/cs



Sterile Water

0.45% Saline

0.9% Saline

Tracheostomy Care and Cleaning Trays

- In standard and customized configurations to meet diverse clinical needs
- Available in economy, basic and full-compartment kits
- May be customized

Did you know?
 We can build a custom kit to your exact specifications. Medline is a leading manufacturer of custom procedure kits.

	DYND40582	DYND40622	DYND40610	DYND40620	DYND40589	DYND40615	DYND40580	DYND40650	DYND40511	DYND40594	CC3T3030A	CC3T4691A	CC4682A	CC4681A
Tray Compartments	3	1	1	1	3	3	3	3	3	3	1	3	3	3
AloeTouch Gloves	2		2	2	2	2	2	2	2		2			2
Drape	1		1	1	1	1			1	1		1	1	1
Pipe Cleaners		2	2							3	3	3	3	3
Non-Woven Gauze 4" x 4"	4	2	2	2	2	2	2	1	3	4	1	3	3	3
Cotton Tip Applicator	2	2	5	5	2	5	2	2	4	2	2	2	3	3
Pliable Cleaning Brush	1	1	1		1	1	1		1	1	1	1	1	1
Twill Tape, 34"	1	1	1	1	1	1	1	1			1	1	1	1
6-Ply Trach Dressing	1		1	1	1	1	1	1	1	1	1	1	1	1
Hydrogen Peroxide					2	2		1		1				
Saline 100 mL					1					1	1			
Water 100 mL								1						
Suction Catheter 14 fr								1		1				
Forceps													1	1
Pkg.	20/cs	30/cs	20/cs	20/cs	20/cs	20/cs	20/cs	20/cs	20/cs	20/cs	30/cs	20/cs	20/cs	20/cs



DYND40589

DYND40582

Kits and Trays

- Irrigation solutions can be poured into the tray or the pop-up cup
- All kits and trays contain Medline Aloetouch® 3G gloves
- Catheters have easy-to-read depth markings
- DeLee tip is flat, soft and smooth to help minimize airway trauma for pediatric patients
- Whistle tip is angled, soft and smooth to maximize suctioning capability for adult patients, while helping minimizing trauma potential

Medline manufactures over 125,000 sterile procedure kits every day.

Suction Catheter Mini Tray with Gloves

Item No.	French Size	Pkg.
DYND40986	6 Fr	100/cs
DYND40988	8 Fr	100/cs
DYND40980	10 Fr	100/cs
DYND40981	12 Fr	100/cs
DYND40982	14 Fr	100/cs
DYND40983	16 Fr	100/cs
DYND40984	18 Fr	100/cs



Suction Catheter Pack with Gloves and Pop-up Cup

Item No.	French Size	Pkg.
DYND40976	6 Fr	50/cs
DYND40978	8 Fr	50/cs
DYND40970	10 Fr	50/cs
DYND40971	12 Fr	50/cs
DYND40972	14 Fr	50/cs
DYND40973	16 Fr	50/cs
DYND40974	18 Fr	50/cs



Suction Catheter Mini Tray with Gloves, Poly-Cath

Item No.	French Size	Pkg.
DYND48980A	10 Fr	100/cs
DYND48982	14 Fr	100/cs



Suction Catheter Trays with Sterile Solutions

- Include everything needed for airway clearance, including a choice of irrigation fluids
- Include Medline Aloetouch® 3G gloves
- Catheters have easy-to-read depth markings
- DeLee tip is flat, soft and smooth to minimize airway trauma for pediatric patients
- Whistle tip is angled, soft and smooth to maximize suctioning capability for adult patients, while helping to minimize trauma potential

Water

Item No.	Description	Pkg.
DYND41476	Suction Catheter Tray with Sterile Water, 100 mL; 6 Fr	36/cs
DYND41478	Suction Catheter Tray with Sterile Water, 100 mL; 8 Fr	36/cs
DYND41470	Suction Catheter Tray with Sterile Water, 100 mL, 10 Fr	36/cs
DYND41471	Suction Catheter Tray with Sterile Water, 100 mL; 12 Fr	36/cs
DYND41472	Suction Catheter Tray with Sterile Water, 100 mL; 14 Fr	36/cs
DYND41474	Suction Catheter Tray with Sterile Water, 100 mL; 18 Fr	36/cs

Saline

Item No.	Description	Pkg.
DYND41478SA	Suction Catheter Tray with Sterile Saline, 100 mL; 8 Fr	36/cs
DYND41470SA	Suction Catheter Tray with Sterile Saline, 100 mL; 10 Fr	36/cs
DYND41442	Suction Catheter Tray with Sterile Saline, 100 mL; 14 Fr	36/cs
DYND40442	Suction Catheter Tray with Saline and Drape; 14 Fr	36/cs

Sterile Solutions

- Convenient, ready-to-use containers reduce the contamination risks and potential waste associated with bulk containers of sterile water and saline
- Ideal for device irrigation, suction, dressing moistening and wound debridement
- Cup size provides the proper amount for suctioning

Water

Item No.	Description	Pkg.
RDI30295	Sterile Water, 100 mL	48/cs
DYND40570	Sterile Water, 110 mL	48/cs
PCS1550	Sterile Water, 250 mL	24/cs
DYND500MLW	Sterile Water, 500 mL	18/cs
DYND1000MLW	Sterile Water, 1,000 mL	6/cs

Saline

Item No.	Description	Pkg.
RDI30296	Sterile Saline, 0.9%, 100 mL	48/cs
DYND40540	Sterile Saline, 0.9%, 110 mL	48/cs
PCS1650	Sterile Saline, 0.9%, 250 mL	24/cs
DYND500MLS	Sterile Saline, 0.9%, 500 mL	18/cs
DYND1000MLS	Sterile Saline, 0.9%, 1,000 mL	6/cs



Tracheostomy Tube Holders

Soft, comfortable straps have adjustable hook-and-loop closures to keep tracheostomy tubes securely in place, ensuring uninterrupted air flow.

Premium Tracheostomy Tube Holders

- Soft, no stretch cotton ensures patient safety and comfort while reducing skin irritation
- Choose from a wide range of sizes

Two-piece

Item No.	Description	Pkg.
HCS540	Adult, Fits up to 20" Neck	10/bx
HCS541	Pediatric, Fits up to 18" Neck	10/bx
HCS542	Neonatal/Infant, Fits up to 9" Neck	10/bx

One-piece

Item No.	Description	Pkg.
HCS597LA	Large Adult, Fits 13–19" Neck	10/bx
HCS597M	Small Adult, Fits 9–14" Neck	10/bx
HCS597S	Neonatal/Pediatric, Fits 7–9" Neck	10/bx

Standard Medline Tracheostomy Tube Holders

- Available in two sizes
- Fits a broad range of patients

Two-piece

Item No.	Description	Pkg. ^a
ZPP501	Adult, Fits up to 18.25" Neck	20/bx
ZPP501P	Pediatric, Fits up to 13.5" Neck	20/bx

One-piece

Item No.	Description	Pkg.
ZPP301	Adult, Fits 18"–20" Neck	20/bx
ZPP301P	Pediatric, Fits 11.5"–13.5" Neck	20/bx
ZPP401	Anti-disconnect	20/bx



HCS540



HCS597M



HCS597LA



ZPP301



ZPP501


SubG™ Endotracheal Tubes

- Help reduce the risk of nosocomial infections due to aspiration of oral and upper airway secretors
- Polyurethane cuff helps to prevent channel formation
- High volume, low pressure barrel-shaped cuff for optimal sealing
- Kink-resistant tubing with easy-to-read insertion markings
- Subglottic suction port for removing secretions above the cuff; suction line cap keeps contaminants out of lumen

Item No.	Description	Pkg.
DYNJ18860	SubG ET Tube, 6 mm	10/bx
DYNJ18865	SubG ET Tube, 6.5 mm	10/bx
DYNJ18870	SubG ET Tube, 7 mm	10/bx
DYNJ18875	SubG ET Tube, 7.5 mm	10/bx
DYNJ18880	SubG ET Tube, 8 mm	10/bx
DYNJ18885	SubG ET Tube, 8.5 mm	10/bx

Advantages of the SubG ET Tube

- Subglottic suction line and evacuation port for removal of pooled secretions above the cuff
- Barrel-shaped cuff made of polyurethane to aid in the reduction of microaspirations
- Thermosensitive material, a bullnose tip and easy to read marking for ease in insertion



Mucus/Specimen Traps


- More secure specimen collection
- Sterile plastic vials, graduated in 1 cc increments
- Includes threaded cap, tubing not made from natural rubber latex and patient ID label

Sterile DeLee Mucus Traps

Item No.	Description	Pkg.
DYND44108	Trap with Catheter and Cap, Capacity 20 cc; 8 Fr	50/cs
DYND44110	Trap with Catheter and Cap, capacity 20 cc; 10 Fr	50/cs

Sterile Mucus Specimen Traps

Item No.	Description	Pkg.
DYND44140	Trap with Transport Cap, Capacity; 40 cc	50/cs
DYND44180	Trap with Transport Cap, Capacity; 80 cc	50/cs



A young girl with dark hair is lying in a hospital bed, wearing a nasal cannula. She is looking towards the camera with a neutral expression. The background shows white hospital linens and a person's hand resting on the bed near her. A dark blue overlay covers the bottom left portion of the image, containing the title and table of contents.

O₂ Delivery

Southmedic OxyMask™	57-58
O ₂ Masks	58
Oxygen and Aerosol Masks	58-60
Face Tents	61
Oxygen Cannulas	62-65
Oxygen Supply Tubing	66
CO ₂ Cannulas & Masks	67-69
Bubble Humidifiers	70
Large Volume Nebulizers	71
Corrugate Tubing	72
Oxygen Kits	73

OxyMask™

Everyone breathes easier with OxyMask

Southmedic OxyMask is a unique open-style oxygen mask with a proprietary design that allows one mask to provide a broad range of oxygen concentrations—while increasing patient comfort, reducing safety risks for patients, and helping to improve overall clinical and operational efficiencies.



Patient Experience



- Reduce risk of CO₂ re-breathing
- Easier communication
- Ability to sip fluids to take oral medications
- Fewer interruptions in care

Caregiver Efficiency



- Reduce risk of errors
- Only one product to learn
- Standardization of care
- Time saving

Operational Efficiency



- Reduce SKUs
- Reduce O₂ consumption
- Save storage space
- Reduce waste

OxyMask™

- Open mask design
- Delivers oxygen flow 1–15+ Lpm at 24%–90% FiO₂

Item No.	Size	Universal Tubing Length	Pkg.
OM11258ML	Adult	7-ft.	25/bx
OP11258ML	Adult Plus	7-ft.	25/bx
OK11258ML	Pediatric	7-ft.	25/bx
OT10258ML*	Infant	7-ft.	25/bx
OM112514ML	Adult	14-ft.	25/bx
OP112514ML	Adult Plus	14-ft.	25/bx
OAPLUS11258**	Adult Plus OxyArm	7-ft.	25/bx

*Infant Mask delivers from ¼ L/min- 6+ at 22%-65% FiO₂

**OxyArm™ delivers 24–50% FiO₂

OxyMask™ Non-PVC

Item No.	Size	Universal Tubing Length	Pkg.
OM70258ML	Adult	7-ft.	25/bx
OP70258ML	Adult Plus	7-ft.	25/bx
OK70258ML	Pediatric	7-ft.	25/bx
OT70258ML*	Infant	7-ft.	25/bx

*Infant Mask delivers from ¼ L/min- 6+ at 22%-65% FiO₂



Visit
[medline.com/
 capabilities/
 respiratory/
 oxy-mask/
 for more info](https://www.medline.com/capabilities/respiratory/oxy-mask/)

OxyMulti-Mask™

- Open mask design delivers oxygen flows from 1–15 liters per minute at a range of 23–83% FiO₂
- Delivers oxygen through standard small bore tube using 22mm adapter or can be used with aerosol nebulizers

Item No.	Description	Universal Tube Length	Pkg.
OHH10258ML	Adult Mask with 22 mm Adapter	7 ft.	25/bx
OPHH10258ML	Adult Plus Mask with 22 mm Adapter	7 ft.	25/bx
OKHH10258ML*	Pediatric Mask with 22 mm Adapter	7 ft.	25/bx
OAT1025ML	Adult Trach Mask	7 ft.	25/bx
OAT1025CSML	Adult Trach Mask with Comfort Neck Strap	7 ft.	25/bx

*Pediatric mask delivers 22–93% FiO₂.



OHH10258ML

Oxygen and Aerosol Masks

The Medline—Hudson RCI portfolio offers a complete range of oxygen and aerosol masks with both standard and universal connections to meet your clinical needs.

- Most oxygen masks are elongated to fit under chin unless otherwise noted
- A secure fit with elastic straps
- Convenient choice of universal or standard connectors
- All oxygen masks come packaged with 7' oxygen tubing unless noted otherwise

Hudson RCI Adult O₂ Masks

Item No.	Description	Connector	Short Style	Reservoir Bag	Safety Vent	Pkg.
HUDRHO41U	Medium Concentration	Standard				50/cs
HUD1041CC	Medium Concentration	ClearConnect				50/cs
HUD1930	Medium Concentration	Universal				50/cs
HUD1049	Medium Concentration	No Tubing				50/cs
HUD1040	Medium Concentration	Standard	●			50/cs
HUD1048	Medium Concentration	No Tubing	●			50/cs
HUD1009	High Concentration	Standard	●			50/cs
HUD1007	High Concentration	Standard		●		50/cs
HUD1059	Non-Rebreather	Standard		●	●	50/cs
HUD1059CC	Non-Rebreather	ClearConnect		●	●	50/cs
HUD1069	Non-Rebreather	No Tubing		●		50/cs
HUD1935	Non-Rebreather	Universal		●		50/cs
HUD1060	Non-Rebreather	Standard		●		50/cs
HUD1936	Non-Rebreather	Universal		●		50/cs
HUD1061	3-in-1 Mask	Standard		●		50/cs
HUD1088	Multi-Vent Venturi Mask	Standard				50/cs
HUD1940	Multi-Vent Venturi Mask	Universal				50/cs
HUD1098	Select-A-Vent Mask	Standard				50/cs
HUD1943	Select-A-Vent Mask	Universal				50/cs



HUDRHO41U



HUD1059
Non-Rebreather with Safety Vent



HUD1060

*3-in-1 Masks come packaged with two 1-way valves, reservoir bag, 7' of oxygen tubing, oxygen connector, and oxygen mask



Medline Adult O2 Masks

Item No.	Description	Connector	Reservoir Bag	Safety Vent	Pkg.
HCS4600B	Medium Concentration	Standard			50/cs
HCSU4600B	Medium Concentration	Universal			50/cs
HCS4640B	Non-Rebreather	Standard	●	●	50/cs
HCSU4640B	Non-Rebreather	Universal	●	●	50/cs
HCS120E	Venturi Mask with 6 Jet Dilluters	Standard			50/cs
HCSU120E	Venturi Mask with 6 Jet Dilluters	Universal			50/cs
HCS121E	Venturi Mask with Single Jet Dilluter	Standard			50/cs
HCSU121E	Venturi Mask with Single Jet Dilluter	Universal			50/cs
HCS123E	Venturi Mask with Single Jet Dilluter and Corrugate Tubing	Standard			50/cs
HCS4670B	Non-Rebreather	Standard	●		50/cs
HCSU4670B	Non-Rebreather	Universal	●		50/cs
HCS4680B*	3-in-1 Mask	Standard	●	●	50/cs
HCSU4680B*	3-in-1 Mask	Universal	●	●	50/cs

Hudson RCI Pediatric O2 Masks

Item No.	Description	Connector	Short Style	Reservoir Bag	Safety Vent	Pkg.
HUD1042	Medium Concentration	Standard				50/cs
HUD1050	Medium Concentration	No Tubing				50/cs
HUD1011	High Concentration	Standard		●		50/cs
HUD1058	Non-Rebreather	Standard		●	●	50/cs
HUD1089	Multi-Vent Venturi Mask	Standard				50/cs
HUD1940	Multi-Vent Venturi Mask	Universal				50/cs
HUD1096	Select-A-Vent Mask	Standard				50/cs
HUD1035	Medium Concentration	Standard	●			50/cs
HUD1000	High Concentration	Standard	●	●		50/cs

Medline Pediatric O2 Masks

Item No.	Description	Connector	Reservoir Bag	Safety Vent	Pkg.
HCS4601B	Medium concentration	Standard		●	50/cs
HCSU4601B	Medium concentration	Universal		●	50/cs
HCS4642B	Non-rebreather	Standard	●	●	50/cs
HCSU4642B	Non-rebreather	Universal	●	●	50/cs
HCS4675B	Non-rebreather	Standard	●		50/cs
HCSU4675B	Non-rebreather	Universal	●		50/cs
HCS4685B*	3-In-1 mask	Standard	●	●	50/cs
HCSU4685B*	3-In-1 mask	Universal	●	●	50/cs

Hudson RCI Infant O₂ Masks

Item No.	Description	Connector	Short Style	Reservoir Bag	Safety Vent	Pkg.
RSH396218	Medium Concentration	Standard	●			50/cs
RSH395499	Medium Concentration	No Tubing	●			50/cs
RSH395407	High Concentration	Standard	●	●		50/cs
RSH395406	High Concentration	No Tubing	●	●		50/cs
RSH395497	Non-Rebreather	Standard	●	●		50/cs
RSH395498	Non-Rebreather	Standard	●	●		50/cs

Aerosol Masks

- Anatomical form and a soft, light and comfortable vinyl material
- Elastic straps
- Standard 22 mm connector

Adult—Short Style

Item No.	Description	Pkg.
HUD1084	Short Style Aerosol Mask	50/cs

Adult—Elongated

Item No.	Description	Pkg.
HCS4630B	Aerosol Mask with Straight Connector	50/cs
HCS4630	Aerosol Mask with Swivel Connector	50/cs
HUD1083	Aerosol Mask with Straight Connector	50/cs
HUD1385	Aerosol Mask with 6 ft. Corr-A-Flex® Corrugated Tubing	50/cs

Pediatric—Short Style

Item No.	Description	Pkg.
HUD1080	Pediatric Aerosol Mask	50/cs

Pediatric—Elongated

Item No.	Description	Pkg.
HCS4631B	Aerosol Mask with Swivel Connector	50/cs
HUD1085	Aerosol Mask with Straight Connector	50/cs

Infant—Short Style

Item No.	Description	Pkg.
RSH395826	Infant Aerosol Mask	50/cs

SideStream Pediatric Aerosol Masks

Item No.	Description	Pkg.
1044382ML	SideStream Tucker the Turtle Pediatric Aerosol Mask	50/cs
HS880ML	SideStream Sami the Seal Pediatric Aerosol Mask	50/cs



HCS4630B



HCS4631B



HCS4630



1044382ML



HS880ML

OxyMask™ Aerosol

- Open bottom design with right angle swivel adapter
- Optimal aerosol delivery, even when patient is laying down or sitting

Item No.	Description	Pkg.
OMN1025ML	Adult Mask with 22 mm Adapter	25/bx
OTN1025ML	Infant Mask with 22 mm Adapter	25/bx
OKN1025ML	Pediatric Mask with 22 mm Adapter	25/bx



Trach Masks and Face Tents

Adult Trach Mask

Item No.	Description	Pkg.
HCS245E	Trach Mask	50/cs
HUD1075	Trach Mask	50/cs

Adult Tent Masks

Item No.	Description	Pkg.
HCS4632B	Face Tent with 22 mm Connector	50/cs
HCS4632BCT	Face Tent with 6 ft. Corrugate Tubing	25/cs
HUD1095	Face Tent Mask	50/cs
HUD1388	Face Tent Mask with 6 ft. Corrugated Tubing	50/cs
HUD1389	Face Tent with 6" Clear Reservoir Tube, 7 ft. Oxygen Tubing and O ₂ Adapter	25/cs

Pediatric Trach Mask

Item No.	Description	Pkg.
HCS4621B	Trach Mask	50/cs
HUD1076	Trach Mask	50/cs



Oxygen Cannulas

Medline offers a variety of nasal cannulas in various lengths, softness and adapters, ensuring you have the right product to meet your clinical or financial goals.

Medline SuperSoft Oxygen Cannulas

- Medline's softest nasal cannula available
- Even short term use of an oxygen cannula can lead to the development of pressure ulcers on the ears and cheeks
- The Medline SuperSoft nasal cannula is made of an advanced material that is soft, flexible and non-tacky. Also features curved, non-flared nasal prongs. You will notice a difference right away
- Promote patient compliance with added comfort—patients are less apprehensive about wearing a cannula that feels soft and comfortable
- Available with a standard or universal connector. Universal connector can be attached directly to compatible flowmeters



HCSS4514S

Item No.	Size	Tubing Length	Star Lumen Tubing	Standard or Universal Connector	Pkg.
HCS4510S	Adult	None			50/cs
HCSS4504S	Adult	4 ft.	●	Standard	50/cs
HCSU4504S	Adult	4 ft.	●	Universal	50/cs
HCSS4514S	Adult	7 ft.	●	Standard	50/cs
HCSU4514S	Adult	7 ft.	●	Universal	50/cs
HCSS4516S	Adult	14 ft.	●	Standard	50/cs
HCSU4516S	Adult	25 ft.	●	Universal	50/cs
HCSS4515S	Adult	25 ft.	●	Standard	50/cs
HCSU4515S	Adult	7 ft.	●	Universal	50/cs
HCSS4518S	Pediatric	7 ft.	●	Standard	50/cs
HCSU4518S	Pediatric	14 ft.	●	Universal	50/cs
HCSS4517S	Infant	7 ft.	●	Standard	50/cs
HCSU4517S	Infant	7 ft.	●	Universal	50/cs
HCSS4520S	Neonatal	50 ft.	●	Standard	50/cs




Crush-resistant tubing
This multi-channel interior feature helps to maintain a strong airflow

Hudson RCI Softech® Nasal Cannulas

- Notch between nasal prongs designed to provide flexibility and a more customized fit for patients (excluding the Neonatal Softech Plus Nasal Cannula, HUD1873)
- Soft durometer lariat tubing to increase patient comfort
- Available with a standard or universal connector. Universal connector can be attached directly to compatible flowmeters

Item No.	Size	Tubing Length	Star Lumen Tubing	Standard or Universal Connector	Pkg.
HUD1870	Adult	7 ft.	●	Standard	50/cs
HUD1874	Adult	14 ft.	●	Universal	50/cs
HUD1877	Adult	14 ft.	●	Standard	50/cs
HUD1871	Pediatric	7 ft.	●	Standard	50/cs
HUD1876	Pediatric	7 ft.	●	Universal	50/cs
HUD1878	Pediatric	14 ft.	●	Standard	50/cs
HUD1872	Infant	7 ft.	●	Standard	50/cs
HUD1873	Neonatal	7 ft.	●	Standard	50/cs



HUD1874

Hudson RCI Softech® and Medline Soft Touch Oxygen Cannulas

- Designed to be softer and offer more comfort than a standard nasal cannula
- The Hudson RCI® Softech Nasal Cannula and Medline Soft Touch cannulas offer a standard over-the-ear design with exceptionally soft nasal prongs
- Extremely lightweight and flexible
- Over-the-ear design designed to optimize fit and stability while maximizing patient comfort
- Available with a standard or universal connector. Universal connector can be attached directly to compatible flowmeters



HUD1820



HUD1818

Hudson RCI Softech

Item No.	Size	Tubing Length	Star Lumen Tubing	Standard Tubing	Standard or Universal Connector	Pkg.
HUD1821	Adult	None				50/cs
HUD1827	Adult	4 ft.	•		Standard	50/cs
HUD1820	Adult	7 ft.	•		Standard	50/cs
HUD1822	Adult	7 ft.	•		Universal	50/cs
HUD1819	Adult	7 ft.		•	Standard	50/cs
HUD1816	Adult	10 ft.	•		Standard	50/cs
HUD1817	Adult	14 ft.	•		Standard	50/cs
HUD1818	Adult	14 ft.	•		Universal	50/cs
HUD1823	Adult	25 ft.	•		Standard	25/cs
HUD1824	Adult	25 ft.	•		Universal	25/cs
HUD1831	Adult	50 ft.	•		Standard	25/cs
HUD1825	Adult	50 ft.	•		Universal	25/cs
HUD1837	Pediatric	None				50/cs
HUD1838	Pediatric	7 ft.	•		Standard	50/cs
HUD1839	Pediatric	14 ft.	•		Standard	50/cs
HUD1828	Infant	7 ft.	•		Standard	50/cs

Medline Soft Touch

Item No.	Size	Tubing Length	Star Lumen Tubing	Standard Tubing	Standard or Universal Connector	Pkg.
HCS4510	Adult	None				50/cs
HCS4504B	Adult	4 ft.	•		Standard	50/cs
HCSU4504B	Adult	4 ft.	•		Universal	50/cs
HCS4514	Adult	7 ft.	•		Standard	50/cs
HCSU4514	Adult	7 ft.	•		Universal	50/cs
HCS4516B	Adult	14 ft.	•		Standard	50/cs
HCSU4516B	Adult	14 ft.	•		Universal	50/cs
HCS4515B	Adult	25 ft.	•		Standard	25/cs
HCSU4515B	Adult	25 ft.	•		Universal	25/cs
HCS4518	Pediatric	7 ft.	•		Standard	50/cs
HCSU4518	Pediatric	7 ft.	•		Universal	50/cs
HCS4517	Infant	7 ft.	•		Standard	50/cs
HCSU4517	Infant	7 ft.	•		Universal	50/cs

Over-The-Ear and Standard Oxygen Cannulas

- The Hudson RCI® Over-the-ear Cannula and Medline standard cannulas are cost-effective and are offered in a wide range of adult configurations
- Flexible, contoured lip tab provides a high level of stability
- Available with a standard or universal connector. Universal connector can be attached directly to compatible flowmeters
- Codes offered with standard (non-flared) or flared nasal prongs



HUD1103

Hudson Standard

Item No.	Tubing Length	Star Lumen Tubing	Standard Tubing	Standard or Universal Connector	Clear Connect	Flared Nasal Prong	Pkg.
HUD1108	None			Standard		●	50/cs
HUD1109	None			Standard			50/cs
HUD1103	7 ft.	●		Standard			50/cs
HUD1104	7 ft.			Standard		●	50/cs
HUD1920	7 ft.		●	Universal			50/cs
HUD1103CC	7 ft.	●			●		50/cs
HUD1810	14 ft.	●		Standard			50/cs
HUD1810CC	14 ft.	●			●		50/cs
HUD1812	25 ft.	●		Standard			25/cs

Medline Standard

Item No.	Tubing Length	Star Lumen Tubing	Standard Tubing	Standard or Universal Connector	Clear Connect	Flared Nasal Prong	Pkg.
HCS4511B	7 ft.	●		Standard			50/cs
HCSU4511B	7 ft.	●		Universal			50/cs
HCS4512B	14 ft.	●		Standard			50/cs
HCSU4512B	14 ft.	●		Universal			50/cs
HCS4513B	25 ft.	●		Standard			50/cs
HCSU4513B	25 ft.	●		Universal			50/cs

Medline SuperSoft High Flow (1–15 LPM) Oxygen Cannulas

- Adult liter flow from 1–15 LPM
- For use with empty bubble humidifiers

Item No.	Tubing Length	Connector	Pkg.
HCS4514HFG	7 ft.	Standard	25/cs
HCS4516HFG	14 ft.	Standard	25/cs
HCS4515HFG	25 ft.	Standard	10/cs



HCS4514HFG

Hudson RCI Comfort Fit™ High Flow (1–15 LPM) Oxygen Cannulas


- Adult liter flow from 1–15 LPM
- For use with empty bubble humidifiers

Item No.	Tubing Length	Connector	Pkg.
HUD1866CS	7 ft.	Standard	50/cs
HUD1867CS	14 ft.	Standard	50/cs

Medline Soft Touch Oxygen Cannulas with Foam Ear Covers

- Pre-attached foam ear covers to help reduce skin irritation and increase patient comfort
- Curved, non-flared nasal prongs

Item No.	Tubing Length	Connector	Pkg.
HCS4514FM	7 ft.	Standard	25/cs
HCS4516FM	14 ft.	Standard	25/cs
HCS4515FM	25 ft.	Standard	25/cs



Foam ear covers

HCS4514FM

Hudson RCI Comfort Fit™ Oxygen Cannula with Foam Ear Covers

Item No.	Tubing Length	Connector	Pkg.
HUD1865	7 ft.	Standard	50/cs

Foam Ear Covers for Oxygen Cannulas

Item No.	Description	Pkg.
HCS0440	EZ Wrap Cushions for Cannulas	50/cs
HUD1868PRCS	Comfort Fit Pre-Slit Ear Cushion	25/cs

Oxygen Supply Tubing

- Selected products with distinctively designed Star Lumen® tubing features a channeled inner lumen to resist kinks and occlusions
- Selected products in green color to provide greater visual contrast to assist in patient safety
- For standard and high flow applications

Hudson RCI Clear Tubing

Item No.	Tubing Length	Smooth Bore Tubing	Standard or Universal Connector	Clear Connect	Pkg.
HUD1115	7 ft.		Standard		50/cs
HUD1925	7 ft.		Universal		50/cs
HUD1115CC	7 ft.			•	50/cs
HUD1991	7 ft.	•	Standard		50/cs
HUD1118	14 ft.		Standard		50/cs
HUD1119	25 ft.		Standard		25/cs
HUD1927	25 ft.		Universal		25/cs
HUD1120CS	50 ft.		Standard		10/cs

Medline Clear Tubing

Item No.	Tubing Length	Smooth Bore Tubing	Standard or Universal Connector	Clear Connect	Pkg.
HCS4507	7 ft.		Standard		50/cs
HCSU4507	7 ft.		Universal		50/cs
HCS4524	14 ft.		Standard		50/cs
HCSU4524	14 ft.		Universal		50/cs
HCS4525	25 ft.		Standard		25/cs
HCSU4525	25 ft.		Universal		25/cs
HCS455010	50 ft.		Standard		10/cs
HCSU455010	50 ft.		Universal		10/cs

Hudson RCI Green Tubing

Item No.	Tubing Length	Standard or Universal Connector	Pkg.
HUD1980	14 ft.	Standard	50/cs
HUD1981	25 ft.	Standard	25/cs
HUD1984	50 ft.	Standard	10/cs

Medline Green Tubing

Item No.	Tubing Length	Standard or Universal Connector	Pkg.
HCS4507G	7 ft.	Standard	50/cs
HCSU4507G	7 ft.	Universal	50/cs
HCS4524G	14 ft.	Standard	50/cs
HCSU4524G	14 ft.	Universal	50/cs
HCS4525G	25 ft.	Standard	25/cs
HCSU4525G	25 ft.	Universal	25/cs
HCS4550GB	50 ft.	Standard	10/cs
HCSU4550GB	50 ft.	Universal	10/cs

ClearConnect Secure Tubing Connector




Universal Connector

Standard Connector



HUD1115CC



HCS4507G

CO₂ Cannulas & Masks

Medline SuperSoft CO₂ Cannulas


- Features a “Pull N’ Peel” design to resist tangling
- Provides accurate, uninterrupted CO₂ sampling and oxygen delivery
- Universal O₂ connector is designed to fit any flowmeter without the need for adapters
- Male, female, or reflective luer connections available to meet caregiver needs



O₂ Delivery


CO₂ Nasal Sampling Cannulas with Male Luer

Item No.	Description	Gas Sampling Line Length	Pkg.
HCS4560S	Adult, 7 ft.	4" Pigtail	25/cs
HCS4562S	Pediatric, 7 ft.	7 ft.	25/cs
HCS4568BS	Adult, 7 ft.	7 ft.	100/cs
HCS4568S	Adult, 7 ft.	7 ft.	25/cs
HCS4571S	Adult, 7 ft.	7 ft.	25/cs
HCS4565S	Adult, 10 ft.	10 ft.	25/cs
HCS4565BS	Adult, 10 ft.	10 ft.	100/cs
HCS4566BS	Adult, 14 ft.	14 ft.	100/cs
HCS4566S	Adult, 14 ft.	14 ft.	25/cs



CO₂ Nasal Sampling Cannulas with Female Luer

Item No.	Description	Gas Sampling Line Length	Pkg.
HCS4561S	Adult, 7 ft.	4" Pigtail	25/cs
HCS4563S	Pediatric, 7 ft.	7 ft.	25/cs
HCS4568FBS	Adult, 7 ft.	7 ft.	100/cs
HCS4568FS	Adult, 7 ft.	7 ft.	25/cs
HCS4572S	Adult, 7 ft.	7 ft.	25/cs
HCS4565FS	Adult, 10 ft.	10 ft.	25/cs



CO₂ Oral Sampling Cannulas with Male Luer

Item No.	Description	Gas Sampling Line Length	Pkg.
HCS4573S	Adult	7 ft.	25/cs



Oral scoop sampling for mouth-breathing patients



4-port nasal sampling allows for dual nares CO₂ sampling and O₂ delivery

CO₂ Oral Sampling Cannulas with Female Luer

Item No.	Description	Gas Sampling Line Length	Pkg.
HCS4574S	Adult	7 ft.	25/cs

Medline MicroFilter CO₂ Sampling

- Compatible with reflective-style capnography monitors
- Intended for short-term applications
- Proven to measure end-tidal CO₂ within 2 mmHg of CapnoLine® Plus
- Four-port design available with or without oral sampling scoop

CO₂ Sampling Cannulas

Item No.	Sampling Type	Size	Universal Tubing Length	Microfilter Line	Pkg.
HCS4607MO	Oral/ nasal	Adult	7 ft.	7 ft.	25/cs
HCS4607MN	Nasal	Adult	7 ft.	7 ft.	25/cs
HCS4707MN	Nasal	Pediatric	7 ft.	7 ft.	25/cs
HCS4614MO	Oral/ nasal	Adult	14 ft.	14 ft.	25/cs
HCS4614MN	Nasal	Adult	14 ft.	14 ft.	25/cs
HCS4714MN	Nasal	Pediatric	14 ft.	14 ft.	25/cs



ET Tube Adapters

- In-line sampling line with 22 mm ID x 15 mm ID adapter
- Intended for short-term applications for adult and pediatric

Item No.	MicroFilter Line	Pkg.
HCS4607ET	7 ft.	25/cs



Accessories

- In-line sampling line intended for short-term application
- For adult and pediatric

Item No.	MicroFilter Line	Pkg.
HCS4607GS	7 ft.	25/ cs



OxyMask™ ETCO₂

- Open design mask that delivers oxygen flows from 1–15 liters per minute
- Range of 24–65% FiO₂ with end-tidal CO₂ monitoring

With Female Pigtail Adapter

Item No.	Size	Universal Tube Length	Female Pigtail Adapter Length	Pkg.
OM21258ML	Adult	7 ft.	8"	25/bx
OP21258ML	Adult Plus	7 ft.	8"	25/bx
OK21258ML	Pediatric	7 ft.	8"	25/bx
OM212514ML	Adult	14 ft.	15"	25/bx



OM21258SLMML

With Male Gas Sampling Line

Item No.	Size	Universal Tubing Length	Male Gas Sampling Line Length	Pkg.
OM21258SLMML	Adult	7 ft.	8 ft.	25/bx
OP21258SLMML	Adult Plus	7 ft.	8 ft.	25/bx
OK21258SLMML	Pediatric	7 ft.	8 ft.	25/bx
OM212514SLM	Adult	14 ft.	15 ft.	25/bx

CO₂ Masks

Southmedic Capnoxygen® Masks

- Allow sampling of exhaled carbon dioxide from the mouth and nose while at the same time administering oxygen
- Soft mask provides form-fitting comfort on the patient's face

With Female Pigtail Adapters

Item No.	Size	Universal Tubing Length	Female Pigtail Adapter Length	Pkg.
CO201ML	Adult	7 ft.	8"	25/bx
CO203ML	Adult Plus	7 ft.	8"	25/bx
CO204ML	Pediatric	7 ft.	8"	25/bx



CO201SLMML

With Male Gas Sampling Lines

Item No.	Size	Universal Tubing Length	Male Gas Sampling Line Length	Pkg.
CO201SLMML	Adult	7 ft.	8 ft.	25/bx
CO203SLMML	Adult Plus	7 ft.	8 ft.	25/bx

Capnoxygen® is a registered trademark of Southmedic INC.

Bubble Humidifiers

Medline-Hudson RCI offers an expanded line of bubble humidifiers, with all Hudson models manufactured in the US. Choose from a wide variety of sizes and configurations.

Hudson RCI AquaPak Prefilled Bubble Humidifer Kits

Designed to add humidity to inspired gas, designed to help control drying and irritation of the respiratory mucosa.

Item No.	Size	Sterile Water	Adapter	Pkg.
HUDRHP340U	340 mL	●	●	20/cs
HUD00540	540 mL	●	●	12/cs
HUD00640	640 mL	●	●	10/cs
HUD00301	340 mL	●	●	20/cs
HUD00501	540 mL	●	●	12/cs
HUD00601	640 mL	●	●	10/cs

Medline Prefilled Bubble Humidifier Kits

Item No.	Size	Sterile Water	Adapter	Pkg.
HCS00350	350 mL	●	●	20/cs
HCS00550	550 mL	●	●	12/cs

Empty Bubble Humidifiers

Item No.	Description	Pkg.
WST0480	Empty Bubble Humidifier, 6 PSI	50/cs
VTB6900	350 mL, Hi-Flow 6 PSI	10/cs

Humidifier Connections for Oxygen Concentrators

Item No.	Description	Pkg.
HCS676	15" Humidifier Connector	50/cs
HCS1790	21" Humidifier Connector	50/cs



HCS00350

HCS00550



HUD00640

HUD00540

HUDRHP340U



VTB6900



HCS1790

Large Volume Nebulizers



Hudson RCI Aquapak® Prefilled Nebulizers

- Select from 2 precision nebulizer adapters for delivery of a consistent aerosol therapy
- Prefilled sterile reservoirs come in 3 convenient sizes: 440 mL, 760 mL, and 1,070 mL
- Available with both sterile prefilled reservoir and adapter conveniently, packaged together
- Reservoirs may be used with other Hudson RCI Humidification and Nebulization Products, making the system truly versatile
- Heater option available separately for heated aerosol therapy

Hudson RCI Aquapak Prefilled Reservoir and Nebulizer Adapter Configurations

Item No.	Size	Sterile Water	Sterile Saline	HUD03128 (028) Adapter	HUD3133 (033) Adapter	Pkg.
HUD00400	440 mL	•				20/cs
HUD00405	440 mL		0.45%			20/cs
HUD04428	440 mL	•		•		20/cs
HUD04433	440 mL	•			•	20/cs
HUD04435	440 mL		0.45%	•		20/cs
HUD04470	440 mL		0.45%		•	20/cs
HUD03700	760 mL	•				10/cs
HUD3735	760 mL		0.45%			10/cs
HUD3709	760 mL		0.90%			10/cs
HUDRHP728U	760 mL	•		•		10/cs
HUD03733	760 mL	•			•	10/cs
HUD3735	760 mL		0.45%	•		10/cs
HUD3739	760 mL		0.90%	•		10/cs
HUD3770	760 mL		0.45%		•	10/cs
HUD04000	1070 mL	•				10/cs
HUD04045	1070 mL		0.45%			10/cs
HUD04128	1070 mL	•		•		10/cs
HUD04133	1070 mL	•			•	10/cs
HUD04135	1070 mL		0.45%	•		10/cs
HUD03128	Adapter only			•		10/cs
HUD3133	Adapter only				•	10/cs



AquaTherm III Heater

- Prefilled nebulizer system designed to provide heated aerosol
- Adjustable aerosol temperature helps meet varying patient needs
- Illuminated power switch for easy on and off

Item No.	Description	Pkg.
HUD05014	Heater for Prefilled Nebulizer	1 ea

Hudson RCI Smooth-Flo® Corr-A-Tube® II Tubing

- Disposable, lightweight tubing available in popular lengths, featuring a corrugated exterior with a smooth interior
- No internal corrugates for condensate to collect
- Smooth interior designed to offer quiet passage of airflow
- 22 mm ID adapts to most respiratory therapy applications

Item No.	Length	Length	Pkg.
HUD1410	6"	22 mm	50/cs
HUD1418	6 ft.	22 mm	50/cs



Hudson RCI Corr-A-Flex Tubing

- 22 mm disposable, flexible, and lightweight corrugated tubing available in popular pre-cut lengths and with cuttable cuffs every 6 in.
- Tapered cuffs securely attach to most 22 mm respiratory therapy applications

Item No.	Length	Cuttable Cuff	Pkg.
HUD1510	6"		50/cs
HUD1518	6 ft.	•	50/cs
HUD1680	100 ft.	•	1 ea



Oxygen Kits



HCSO2KIT8

Medline can customize an oxygen set up kit to your needs. Contact your Medline Representative.

Item No.	Description	Pkg.
HCS4514HFKG	7 ft. High-Flow Cannula, High-Flow Empty Bubble Humidifier	7/cs
HCS4516HFKG	14 ft. High-Flow Cannula, High-Flow Empty Bubble Humidifier	7/cs
HCSO2KIT6	7 ft. Soft Touch Adult Cannula (2), 4 ft. Soft Touch Adult Cannula, 50 ft. Green Oxygen Tubing, Swivel Oxygen Connector, Nut and Nipple Green Adapter, 6 psi Empty Bubble Humidifier, No Smoking Sign	10/cs
HCSO2KIT7	7 ft. Soft Touch Adult Cannula (2), 4 ft. Soft Touch Adult Cannula, 25 ft. Green Oxygen Tubing, Swivel Oxygen Connector, Nut and Nipple Green Adapter, 6 psi Empty Bubble Humidifier, No Smoking Sign	18/cs
HCSO2KIT8	7 ft. Soft Touch Adult Cannula (4), 50 ft. Green Oxygen Tubing, Swivel Oxygen Connector, Nut and Nipple Green Adapter, Water Trap, 6 psi Empty Bubble Humidifier, Humidifier Connector, Cylinder Wrench, No Smoking Sign	20/cs
HCSO2KIT9	7 ft. SuperSoft Adult Cannula (6), 50 ft. Green Oxygen Tubing, Swivel Oxygen Connector, Nut and Nipple Green Adapter, No Smoking Sign	15/cs

Medication Delivery


Medication Nebulizers	75-77
Aerosol Masks	78
MDI Chambers/Spacers	79-80
Nebulizer Compressors	81-82



Medication Nebulizers

Hudson RCI Micro Mist® Nebulizers

- Consistent performance at angles up to 90 degrees allowing flexible medication delivery
- Jet stays in place unless intentionally removed
- Easy-seal, threaded, cap and 6 cc capacity anti-spill jar to prevent leakage
- Recommended flow rate: 8 LPM; nebulization rate: >20 mL/min
- Available with standard, universal, or ClearConnect™ Secure Tubing Connector



HUDRHS883U

Item No.	Nebulizer	Tee	Mouthpiece	Reservoir Tube	7' Star Lumen Tubing	Connector	Mask	Pkg.
HUD1880	•							50/cs
HUD1881	•	•	•					50/cs
HUD1882	•	•	•		•	Standard		50/cs
HUDRHS883U	•	•	•	•	•	Standard		50/cs
HUD1883CC	•	•	•	•	•	ClearConnect		50/cs
HUD1884	•	•	•	•	•	Universal		50/cs
HUDRHS885U	•				•	Standard	Adult, Elongated (HUD1083)	50/cs
HUD1886CS	•				•	Standard	Pediatric (HUD1085)	50/cs
HUD1945	•				•	Universal	Adult	50/cs
HUD1946	•				•	Universal	Pediatric	50/cs

Hudson RCI Up-Draft.II® Opti-Neb® Nebulizers

- Nebulizes at angles up to 45° to accommodate special treatment requirements
- Produces fine, dense particles for effective therapy, with easy-seal, threaded cap and 8 cc capacity jar
- Anti-spill design prevents loss of medication in any position
- Available with standard, universal, or ClearConnect™ Secure Tubing Connectors
- Recommended flow rate: 5 to 9 LPM



HUDRHS705U

Item No.	Nebulizer	Tee	Mouthpiece	Reservoir Tube	7' Star Lumen Tubing	Connector	Mask	Pkg.
HUDRHS705U	•				•	Standard	Adult, Elongated (HUD1083)	50/cs
HUD1706	•				•	Standard	Adult (HUD1084)	50/cs
HUD1707	•				•	Standard	Pediatric (HUD1085)	50/cs
HUD1730	•							50/cs
HUD1731	•	•	•					50/cs
HUD1732	•	•	•		•	Standard		50/cs
HUD1734	•	•	•	•	•	Standard		50/cs
HUD1737	•	•	•	•	•	Universal		50/cs
HUD1734CC	•	•	•	•	•	ClearConnect		50/cs

SideStream® Nebulizers

- Offers treatment times up to twice as fast compared to conventional nebulizers with an improved respirable fraction
- Diamond Jet provides a high respirable fraction and consistent delivery
- Four minute average treatment time

Disposable

Item No.	Description	Pkg.
HS800ML	SideStream Nebulizer Kit	50/cs
1028219ML	SideStream Nebulizer (neb cup only, no additional components)	50/cs
1028125ML	SideStream Custom Nebulizer Kit (2 nebulizer kits + 2 replacement cups)	25/cs
1036160ML	SideStream Nebulizer Kit with Adult Mask	50/cs
1041906ML	SideStream Nebulizer Kit with Pediatric Mask	50/cs
1025529ML	SideStream Pediatric Aerosol Mask	50/cs
1025531ML	SideStream Adult Aerosol Mask	50/cs



HS800ML



HS860050ML

Reusable

Item No.	Description	Pkg.
HS860050ML	SideStream Reusable Nebulizer Kit	50/cs
1072286ML	SideStream Reusable Nebulizer Kit with Tucker the Turtle Mask	10/cs

Hudson RCI ISO-Neb™ Filtered Nebulizer System

- Anti-spill Up-Draft II® Opti-Neb® Nebulizer produces fine, dense particles for site-specific therapy
- High-output, 8 cc capacity nebulizer operates vertically or at a 45° angle
- Recommended flow rate: 6 to 8 LPM

Item No.	Description	Pkg.
HUD1755	ISO-Neb Filtered Nebulizer System	20/cs



HUD1755

Hudson RCI Neb-U-Mask® Systems

- The system includes a wye design featuring a nebulizer connection and MDI port, a high-concentration mask, and a 750 mL gas reservoir bag
- The nebulizer connection features a valved port, which maintains a closed system when the nebulizer is not in use
- Packaged with a Micro Mist® Nebulizer and color-coded gas supply tubing to facilitate ease of use
- Nebulizes at angles up to 90° and features a 360° swivel connector, allowing the patient to find a comfortable position while using the system

Item No.	Description	Pkg.
HUD1895	Adult Nebulizer and MDI Ports	20/cs
HUD1896	Nebulizer with Pediatric-Sized Mask and MDI Ports	20/cs



HUD1895

Medline Aeromist BAN

- Delivers the performance you require from a high-efficiency nebulizer, with additional features to improve operational efficiency
- Consistent drug delivery with minimal waste; demonstrated clinical efficacy in lab testing; easily switch between breath-actuated and continuous mode
- Fast treatment time improves patient satisfaction and saves caregiver time
- Limited drug loss to the environment reduces caregiver exposure to aerosol
- Connects directly with a face mask, no additional pieces needed



Standard Breath-actuated Nebulizers

Item No.	Mouthpiece	Mask Type	Tubing Length	Tubing Connector	Pkg.
HCS4433BAN	●		7 ft.	Standard	50/cs
HCSU4433BAN	●		7 ft.	Universal	50/cs
HCS4435BAN		Adult	7 ft.	Standard	50/cs
HCSU4435BAN		Adult	7 ft.	Universal	50/cs
HCS4436BAN		Pediatric	7 ft.	Standard	50/cs
HCSU4436BAN		Pediatric	7 ft.	Universal	50/cs

Filtered Breath-actuated Nebulizer

Item No.	Mouthpiece	Mask Type	Tubing Length	Tubing Connector	Pkg.
HCS4433FBAN	●		7 ft.	Standard	50/cs



HCS4433BAN/
HCSU4433BAN



HCS4433FBAN

Aerosol Masks

- Extra patient comfort comes from an anatomical form and a soft, light and comfortable vinyl material
- All aerosol masks feature an elongated fit for optimal patient fit
- A secure fit with elastic straps
- Standard 22 mm connector

Adult

Item No.	Description	Pkg.
HCS4630B	Aerosol Mask with Straight Connector	50/cs
HCS4630	Aerosol Mask with Swivel Connector	50/cs
HUD1083	Aerosol Mask with Straight Connector	50/cs

Pediatric

Item No.	Description	Pkg.
HCS4631B	Aerosol Mask with Swivel Connector	50/cs
HUD1085	Aerosol Mask with Straight Connector	50/cs

SideStream

Item No.	Description	Pkg.
1044382ML	SideStream Tucker the Turtle Pediatric Aerosol Mask*	50/cs
HS880ML	SideStream Sami the Seal Pediatric Aerosol Mask**	50/cs

*only SideStream disposable and reusable

**SideStream and SideStream Plus



HUD1083



HS880ML



1044382ML

MDI Chambers/Spacers

OptiChamber® Diamond Dual Valved Holding Chambers

- Small, intuitive design enhances medication delivery and compliance for all patients
- Detachable LiteTouch mask has a comfortable, soft-seal interface, designed to create an optimal seal with minimal force
- Anti-static valved holding chamber is smaller than most conventional chambers which enhance medication delivery—even under low pediatric pressures and flows, the valves open freely
- Stepped mouthpiece helps patients transition from pediatric facemasks to mouthpiece; allows attachment to a 22 mm connector
- Mouthpiece and adapter can be removed from the chamber for easy cleaning



Item No.	Description	Pkg.
1077478ML	OptiChamber Diamond	10/cs
HCS1143035	OptiChamber Diamond	50/cs
1079823ML	OptiChamber Diamond with Small LiteTouch Mask	10/cs
1079826ML	OptiChamber Diamond with Medium LiteTouch Mask	10/cs
1079827ML	OptiChamber Diamond with Large LiteTouch Mask	10/cs



1077478ML

LiteTouch Masks

- Provide greater comfort and a more secure seal for increased treatment effectiveness
- Mask quickly and comfortably conforms to facial contours, providing a soft, secure fit for increased compliance and treatment effectiveness
- Flexible seal makes treatment more comfortable and less technique dependent
- Ensures a proper fit with as little as one pound of force

Item No.	Description	Pkg.
1082714ML	LiteTouch Mask, Small	10/cs
1082715ML	LiteTouch Mask, Medium	10/cs
1082716ML	LiteTouch Mask, Large	10/cs



InspiraChamber® Valved Holding Chambers

- Anti-static chamber allows the medication particles to be suspended, giving the patient time to take a slow, deep breath¹
- One-way valve and exhalation port designed to direct exhaled medication away from patient's face
- Flexible MDI adapter end cap accommodates most commonly used pMDIs
- Available with two types of masks in a variety of sizes to meet patient needs
 - InspiraMask®, available in small, medium, and large sizes, is designed with soft, lightweight material to provide an effective seal with a gentle touch
 - SootherMask®, available in small and medium sizes, is designed with a unique, pacifier-ready slot, which may help to keep children asleep during administration for enhanced therapeutic delivery¹



HUD100150

Item No.	Description	Pkg.
HUD100110	InspiraChamber	10/cs
HUD100150	InspiraChamber	50/cs
HUD100110SS	InspiraChamber with Small SootherMask	10/cs
HUD100110SM	InspiraChamber with Medium SootherMask	10/cs
HUD100110S	InspiraChamber with Small InspiraMask	10/cs
HUD100110M	InspiraChamber with Medium InspiraMask	10/cs
HUD100110L	InspiraChamber with Large InspiraMask	10/cs



HUD100110SS

HUD100110S

1. Amirav I, Luder AS, Halamish A, Marzuk C, Daitzchman M, Newhouse MT. Computerized dead-space volume measurement of face masks applied to simulated faces. *Respir Care*. 2015 Sep;60(9):1247-51. doi: 10.4187/respcare.03813. Epub 2015 May 5. Dr. Newhouse is the inventor of the InspiraChamber® and SootherMask®.

Nebulizer Compressors



Aeromist Compact™ Nebulizer Compressors

- Our smallest nebulizer compressor in standard white model or color options
- Kits include a disposable nebulizer medicine cup, tee mouthpiece, 6" corrugate reservoir, 7 ft.(2.13 m) of tubing, and 5 replacement filters
- Durable hard plastic shell
- Five-year warranty
- Replacement Kits: HCS4483 (Disposable) and HCS4483R (Reusable)

Item No.	Description	Pkg.
HCS70004	White Compressor and Nebulizer Kit	8/cs
HCS70004RD	White Compressor with Reusable and Disposable Neb Kit	8/cs
HCS70004G	Neon Green Compressor and Nebulizer Kit	8/cs
HCS70004BL	Electric Blue Compressor and Nebulizer Kit	8/cs
HCS70004P	Neon Pink Compressor and Nebulizer Kit	8/cs
HCS70004C	Charcoal Compressor and Nebulizer Kit	8/cs
HCS70004ADLT	White Compressor and Nebulizer Kit with Adult Mask	8/cs
HCS70004DP	Sample Pack, 2 of Each Color	8/cs
HCS70004F	Replacement Filters	5/pk
HCS70004CAP	Replacement Cap for HCS70004	8/pk



HCS70004



Components included with each compressor

Specifications

- Method of nebulization: dry piston compressor
- Electrical requirements: 120 VAC/60 Hz
- Maximum jet nebulizer capacity: 10 mL
- Sound noise level: ≤55 dBA
- Operation air pressure: 9–16 PSI, Max air pressure: 35–50 PSI
- Operation air flow rate: 5–8 LPM



Scan here to view our nebulizer compressor in-service video

Aeromist Buddies™ Nebulizer Compressors

- Kid-friendly nebulizers help put young patients at ease, while providing safe and effective aerosol treatments
- Kit includes Disposable Nebulizer Medicine Cup, 7 ft.(2.13 m) of tubing, 6" Reservoir Tube, T-Mouthpiece, Pediatric Mask, and 5 Replacement Filters
- Durable hard plastic shell
- Five-year warranty

Item No.	Description	Pkg.
HCSMONKEY	Monkey Compressor and Nebulizer Kit	8/cs
HCSMONKEYCAP	Replacement Filter Cap for HCSMONKEY	8/pk
HCSFISHNEB*	Fish Compressor and Nebulizer Kit	8/cs
HCSPUPPY*	Puppy Compressor and Nebulizer Kit	8/cs
HCSPUPPYCAP	Replacement Filter Cap for HCSPUPPY	8/pk
HCSPUPPYRD*	Puppy Compressor with Disposable and Reusable Nebulizer Kit	8/pk
HCSTURTLE	Turtle Compressor and Nebulizer Kit	8/cs
HCSTURTLECAP	Replacement Filter Cap for HCSTURTLE	8/pk
HCS70004F	Replacement Filter	5/pk

*Comes with Grab-N-Go backpack for convenient storage.

All Compressors include a complementary Nebulizer Kit

Specifications

- Dry piston compressor nebulization
- Electrical requirements: 120 VAC/60 Hz
- Maximum nebulizer capacity: 10 mL
- Sound noise level: ≤60 dBA
- Operation air pressure: 10–17 PSI, Max air pressure: 35–50 PSI
- Operation air flow rate: 4–7.5 LPM



HCSTURTLE

HCSPUPPY

HCSFISHNEB

Kid-friendly
Nebulizers

Replacement Nebulizer Kits for Compressors

- We offer a wide selection of reusable and disposable nebulizer kits for medication delivery
- Ideal particle size of less than five microns, increasing the amount of medication delivered deep within lungs
- Standard fittings accommodate most clinical applications
- High respirable dose, leading to more efficient treatments and less wasted medication

Item No.	Adapter	Mouthpiece	Reservoir Tube	Tubing Length	Tubing Connector	Mask Type	Disposable/ Reusable	Pkg.
HCS4482	22 mm OD/ID	●		7 ft.	Standard		Disposable	50/cs
HCS4483	22mm OD/ID	●	●	7 ft.	Standard		Disposable	50/cs
HCS4485				7 ft.	Standard	Adult	Disposable	50/cs
HCS4486				7 ft.	Standard	Pediatric	Disposable	50/cs
HCS4482R	22 mm OD/ID	●		7 ft.	Standard		Reusable	50/cs
HCS4483R	22mm OD/ID	●	●	7 ft.	Standard		Reusable	50/cs
HCS4482RR	45° angled	●		7 ft.	Standard		Reusable	50/cs



Bronchial Hygiene

Incentive Spirometers	84
Peak Flow Meters	85
PEP Devices	86

Incentive Spirometers

Hudson RCI Voldyne® 4000 Spirometer

- Volume measurement up to 4,000 mL accommodates large inspired volumes
- Good-Better-Best flow window
- Compact, ergonomic design
- Built-in handle

Item No.	Description	Pkg.
HUD8884719010	Incentive Spirometer, 4,000 mL	12/cs



HUD8884719010

Hudson RCI Voldyne® 2500 Spirometer

- Volume measurement up to 2,500 mL accommodates smaller inspired volumes
- Good-Better-Best flow window
- Compact, ergonomic design
- Built-in handle

Item No.	Description	Pkg.
HUD8884719011	Incentive Spirometer, 2,500 mL	12/cs



HUD8884719011

Hudson RCI Voldyne® Pediatric Spirometer

- Volume measurement up to 2,500 mL accommodates pediatric inspired volumes
- Colorful, child-oriented labeling to facilitate patient evaluation and use
- Pediatric-sized mouthpiece
- Compact, ergonomic design
- Built-in handle; includes fun puzzle/coloring book

Item No.	Description	Pkg.
HUD8884719012	Pediatric Incentive Spirometer, 2,500 mL	12/cs



HUD8884719012

Hudson RCI TriFlo II™ Inspiratory Exerciser

- Wide flow rate range from 600 to 1200 cc/sec
- 3 chambers, 3 color-coded balls, one in each chamber
- Minimum flow imprinted on each chamber

Item No.	Description	Pkg.
HUD717301	Inspiratory Exercisor	12/cs



HUD717301

Peak Flow Meters

Personal Best® Peak Flow Meters

- Lightweight meter helps patients monitor their asthma
- Handheld device measures Peak Expiratory Flow Rate so patients can initiate preemptive therapy in order to help prevent serious asthma attacks
- Universal range allows the Personal Best peak flow meter to be used for all patients, from pediatric to adult
- Integrated case for storage
- Mouthpieces available separately

Item No.	Description	Pkg.
HS755ML	Personal Best Flow Meter, Full-Range, 60–810 LPM	12/cs
HS756ML	Personal Best Flow Meter, Low-Range, 50–390 LPM	12/cs
HS714200ML	Mouthpiece for Personal Best Peak Flow Meter, One-Way Valve, Adult	200/cs
HS715100ML	Mouthpiece for Personal Best Peak Flow Meter, Pediatric	100/cs
HS757100ML	Mouthpiece for Personal Best Peak Flow Meter, Oval, Adult	100/cs



HS755ML

AsthmaCheck Peak Flow Meter

AsthmaCheck offers the reliability of a highly accurate and portable peak flow meter due to its small, compact design.

- Easy-grip handle directs hands away from the airflow path and ensures accurate measurements every time.
- Oval mouthpiece for patient comfort

Item No.	Description	Pkg.
HS740ML	AsthmaCheck, Peak Flow Meter, 60–810 LPM	12/cs



HS740ML

AsthmaMD® Peak Flow Meter

- Portable, pocket-sized asthma management tool helps patients monitor lung performance and anticipate attacks
- Colorful integrated markers help patients easily monitor changes in air flow with large, visible scale

Item No.	Description	Pkg.
HUD1802	AsthmaMD Peak Flow Meter, 60–810 LPM	10/cs



HUD1802

Medline Peak Flow Meters

- Individually calibrated to ensure accuracy
- Durable for years of reliable performance
- The pocket model offers superb convenience and portability

Item No.	Description	Pkg.
HCSPF100	Portable Peak Flow Meter, Full Range, 60–720 LPM	12/bx
HCSPF220	Disposable Pediatric Mouthpiece, Adapter	20/bx
HCSPF262	Disposable One-way Valve Mouthpiece, Adult	100/cs



HCSPF220



HCSPF262



HCSPF100

Threshold PEP and IMT Threshold Systems

- Used for airway clearance, bronchial hygiene, or as an alternative to chest physical therapy
- Incorporates a flow-independent, one-way valve to provide consistent resistance and features an adjustable specific pressure setting (in cm H₂O) to be set by a healthcare professional
- When patients inhale through Threshold IMT, a spring-loaded valve provides resistance that exercises respiratory muscles through conditioning
- When patients exhale through Threshold PEP, the resistive load creates positive pressure that helps open the airways and allows mucus to be expelled during “huff” coughing (forced expiratory technique)

Item No.	Description	Pkg.
HS730ML	Threshold IMT	10/cs
HS735ML	Threshold PEP	10/cs



HS730ML



HS735ML



Adapters and Connectors

General	88–89
Elbows	90
Wyes	91
MDI Adapters	91–92
NeoNatal Adapters	92
ConchaTherm® Adapters	93
Drainage Systems	93
One-Way Valves	94
Temperature Monitoring	94
Water Traps	95
Tees	95–96
Disposable Nose Clips	96
Disposable Mouthpiece	96
Miscellaneous	97

Adapters and Connectors

Medline offers an extensive line of Hudson RCI® Adapters and Connectors, including elbows, wyes, tees, and much more, in a wide range of sizes and designs, to address your respiratory needs.



HUD1640 Tubing Adapter

- Standard 22 mm OD x 22 mm ID
- Packaged 50/cs



HUD1421 Universal Cuff Adapter

- Standard 22 mm ID connections
- Packaged 50/cs



HUD1769 Step-down Adapter

- 15 mm OD x 22 mm ID
- Packaged 50/cs



HUD1678 22-mm ISO Cuff Adapter

- 22 mm ID x 22 mm OD connections
- Packaged 50/cs



HUD1424 Swivel Oxygen Tubing Connectors

- Provide exceptionally smooth, 360° rotation
- Available with optional clip to secure tubing to clothing or bedding
- Male/male connection
- Individually packaged 20/cs



HUD1425 Swivel Oxygen Tubing Connectors

- Provide exceptionally smooth, 360° rotation
- Available with optional clip to secure tubing to clothing or bedding
- Male/female connection
- Individually packaged 20/cs



HUD1768 Step-down Adapter

- 15 mm OD x 22 mm OD
- Packaged 25/cs



HUD1746 Nebulizer Adapter

- 22 mm OD x 18 mm OD
- Packaged 50/cs



**HUD1663
Soft Cuff Adapter**

- Approximately 21 mm ID x 19 mm ID
- Made of flexible material
- Packaged 50/cs



**HUD1423
Tubing Adapter**

- 22 mm ID x 5 to 7 mm tapered nipple
- Can be used with trach mask (HUD1075)
- Packaged 50/cs



**HUD2555
Nut and Nipple Adapter**

- Hose adapter for standard DISS oxygen outlets
- Packaged 50/cs



**HUD2556
Nut and Nipple Adapter**

- Yellow hose adapter for standard DISS outlets
- Packaged 50/cs



**HUD1422
Multi-Adapter**

- Standard 15 mm ID x 22 mm OD
- Packaged 50/cs



**HUD1409
Large Tubing Connector**

- Standard 22 mm OD x 22 mm O.D. connections
- Packaged 50/cs



**HUD1642
Pressure Line Adapter**

- 22 mm OD x 22 mm ID connections
- Elbow swivels 360°
- Barbed elbow accepts 5 to 7 mm ID tubing
- Packaged 50/cs



**HCS676
Humidifier Connectors**

- Adapter for connecting an oxygen concentrator to a humidifier
- 15" length tube
- Packaged 50/cs



**HUD2557
Nut and Nipple Adapter**

- Clear hose adapter for standard DISS outlets
- Packaged 50/cs



**HUD1420
Tubing Connector**

- Couples oxygen supply tubing using 5 to 7 mm end connectors
- Packaged 50/cs



**HCS1790
Humidifier Connectors**

- Adapter for connecting an oxygen concentrator to a humidifier
- 21" length tube
- Packaged 50/cs

Elbows



HUD1623 Patient Elbow Adapter

- 15 mm OD x 15 mm ID/22 mm OD
- Packaged 50/cs



HUD1641 Elbow Adapter

- Standard 22 mm ID x 22 mm OD
- Packaged 50/cs



HUD1673 Ported Elbow Adapter

- Standard 22 mm ID x 22 mm OD connections
- Accepts most temperature probes
- Packaged 50/cs



HUD1624 Gas Sampling Adapter

- 15 mm OD x 15 mm ID/22 mm OD
- Female luer-lock connection swivels 360°
- Packaged 50/cs



HUD1632 Norman Mask Elbow

- 15 mm OD x 15 mm ID/22 mm OD
- Elbow accepts 5- to 7-mm ID tubing, swivels 360°
- Packaged 50/cs



HUD1682 CPAP Elbow

- For use on infant nasal CPAP sets and systems
- Packaged 50/cs



RSH395282 Barbed Elbow Connector

- 22 mm OD connection
- Packaged 50/cs

Wyes



HUD1660 Universal Wye Connector

- 22 mm OD tubing connections
- 15 mm ID/22 mm OD patient connection
- Airway pressure and temperature monitoring ports
- Packaged 50/cs



HUD1661 Wye Connector

- Same as HUD1660, but without pressure monitoring port
- Packaged 50/cs



HUD1662 Wye Connector

- Same as HUD1660, but without temperature and pressure monitoring ports
- Packaged 50/cs



HUD1674 Adult Near Parallel Wye

- 22 mm OD tubing connections
- 15 mm ID/22 mm OD patient connection
- No ports
- Packaged 50/cs

MDI Adapters



HUD1659 MDI Adapter

- 22 mm OD x 22 mm ID connections for in-line placement in 22 mm tubing
- Packaged 50/cs



HUD1751 MDI Adapter

- 15 mm OD x 15 mm ID connections for in-line placement in 15 mm tubing
- Packaged 50/cs



HUD1759 Dual Spray MDI Adapter

- Bidirectional nozzle creates a chamber effect and directs plumes in both directions simultaneously
- 22 mm OD x 22 mm ID connections
- Extended universal port
- Tethered cap
- Packaged 50/cs



SMDIA1000ML MDI Inline Adapter

- 22 mm ID → 22 mm OD/15 mm ID
- Can be used with or without dose counter
- Includes built-in tether strap
- Packaged 25/cs



SMDIA1100ML
MDI Inline Adapter

- 22 mm OD → 15 mm ID/ 22 mm ID
- Can be used with or without dose counter
- Includes built-in tether strap
- Packaged 25/cs



SMDIA1200ML
MDI Inline Adapter

- 22 mm OD → 22 mm OD/ 15 mm ID
- Can be used with or without dose counter
- Includes built-in tether strap
- Packaged 25/cs



SMDIA1300ML
MDI Inline Adapter

- 15 mm OD → 15 mm ID
- Can be used with or without dose counter
- Includes built-in tether strap
- Packaged 25/cs

Neonatal Adapters



HUD1763
Neonatal Probe Adapter

- Standard 10 mm OD connections
- Accepts most probes and thermometers
- Packaged 50/cs



HUD1402
Drive Line Adapter

- 0.152 in OD connection
- 0.092 in ID connection
- Packaged 100/cs



HUD1792
Neonatal Adapter Kit

- Includes 18 mm nebulizer adapter, neonatal tee, and 12 in., 10 mm corrugated tubing
- Packaged 50/cs



HUD1764
Neonatal Adapter

- 10 mm OD x 7.5 mm OD
- For connecting 10 mm ID and 7.5 mm ID orifices
- Packaged 100/cs



HUD1765
Neonatal Pressure Line Tee

- For connecting 1/8 in. ID tubing to 1/8 in. pressure line
- Packaged 100/cs



HUD1766
Neonatal Pressure Line Tee

- For connecting 1/8 in. ID tubing to 3/16 in. pressure line
- Packaged 100/cs

ConchaTherm® Adapters



HUD38400 ConchaTherm® Tubing Adapter

- Tapered fitting up to 7 mm
- 15.5 mm ID/19 mm OD
- Designed to be used with Concha® heated humidifier
- Packaged 30/cs



HUD38426 ConchaTherm® Oxygen Diluter

- Designed to be used with Concha® heated humidifier
- Allows delivery of heated humidity at specific oxygen concentrations
- Packaged 30/cs



HUD1645 Thermometer Adapter

- Standard 22 mm OD x 22 mm ID connections
- Port accepts Concha® temperature probes
- Packaged 50/cs

Drainage Systems



HUD1740 Aerosol Drainage System

- Disposable drainage bag with tee adapter accepts standard 22 mm ID corrugated tubing
- 750-cc bag capacity
- Self-locking plastic chain secures bag to a stable object
- Convenient drainage port
- Individually packaged, 50/cs



HUD1742 Aerosol Drainage System

- Disposable drainage bag with wye adapter accepts standard 22 mm ID corrugated tubing
- 750-cc bag capacity
- Self-locking plastic chain secures bag to a stable object
- Convenient drainage port
- Individually packaged, 50/cs

One-Way Valves



HUD1644 One-Way Valve

- Standard 22 mm OD x 22 mm ID connections
- Valve offers low flow resistance
- Capped monitoring port
- Packaged 50/cs



HUD1666 Unidirectional Valved Tee

- Assembled tee and two one-way valves allow gas to flow in only one direction
- Packaged 50/cs

Temperature Monitoring



HUD1635 Thermometer

- Clear housing with alcohol-based red liquid column
- Fahrenheit (50°-120°) and Celsius (10°-50°) scales
- Packaged 50/cs



HUD1637 Thermometer with Tee Adapter

- Combines Item No. HUD1635 and Item No. HUD1648
- Packaged 50/cs



HUD1647 Thermometer with Thermometer Adapter

- Combines HUD1635 and HUD1645
- Packaged 50/cs

Water Traps



HUD1650 Water Trap

- 70-cc capacity
- 22 mm OD x 22 mm OD tubing connectors
- Self-sealing lid to ensure circuit patency
- Packaged 50/cs



HUD1658 Neonatal Water Trap

- 20–30 cc capacity
- Self-sealing lid to ensure circuit patency
- Packaged 50/cs



HUD1679 O₂ Tubing Water Trap

- Collects moisture that can accumulate in O₂ tubing
- Approximately 8" long with a 1" diameter
- Packaged 250/cs

Nebulizer Tees



HUD1638 Nebulizer Tee Connector

- 22 mm OD x 18 mm ID / 22 mm OD x 22 mm OD
- Port accepts most small volume nebulizers
- Packaged 50/cs



HUD1639 Nebulizer Tee Connector

- 22 mm OD x 18 mm ID / 22 mm ID
- Port accepts most small volume nebulizers
- Packaged 50/cs

Aerosol Tees



HUD1648 Aerosol Tee Connector

- 22 mm OD x 15 mm ID / 22 mm OD x 22 mm OD
- Packaged 50/cs



HUD1077 Aerosol Tee Connector

- Anti-spill design minimizes aspiration of liquids
- 22 mm OD x 15 mm ID / 22 mm OD x 22 mm OD
- Packaged 50/cs



HCSA1008 Aerosol Tee

- Valved, 22 mm OD/OD
- Can be used to introduce aerosol into the circuit
- Spring-loaded valve opens and closes as nebulizer is inserted and removed
- Packaged 30/cs



HCSA1009 Aerosol Tee

- Valved, 22 mm OD/ID
- Can be used to introduce aerosol into the circuit
- Spring-loaded valve opens and closes as nebulizer is inserted and removed
- Packaged 30/cs

Valved Tees



HUD1743 Valved NebTee

- 22 mm OD x 18 mm ID x 22 mm ID connections
- Tee may be left in line when nebulizer is removed
- Port accepts most small-volume nebulizers
- Packaged 30/cs



HUD1744 Valved NebTee


- 22 mm OD x 18 mm ID x 22 mm OD connections
- Tee may be left in line when nebulizer is removed
- Port accepts most small-volume nebulizers
- Packaged 30/cs

Disposable Nose Clips



DYNJ750 Standard Nose Clip

- Packaged 100/bx



DYNJ751 Nose Clip with Marshmallow Pads

- Packaged 100/bx



DYNJ752 Large Nose Clip

- Packaged 50/bx



DYNJ755 Nose Clip and Mouthpiece

- Packaged 100/bx

Disposable Mouthpiece



DYNJ754 White Thermoplastic Mouthpiece

- Packaged 50/bx

Miscellaneous



HUD1668 Trach Tee Oxygenator

- Allows administration of oxygen to endotracheal or tracheostomy tubes
- Includes oxygen nipple adapter, aerosol tee, and 6-in. Corr-A-Flex® tube
- Packaged 50/cs



HUD1565 Mouthpiece

- Standard mouthpiece with 15 mm/22 mm OD connection
- Packaged 50/cs



HCS3222 Flex Connector

- Flexible adult connector, 15 mm x 22mm
- Used on a ventilator and connects to tubing
- Packaged 50/cs



HCS3215 Flex Connector

- Flexible pediatric connector, 15 mm x 15 mm
- Used on a ventilator and connects to tubing
- Packaged 50/cs



HUD1643 Exhalation Valve

- Includes drive line and hanger
- Packaged 50/cs



HUD2417 Comfort Flo® 15 mm Oxygen Tubing Adapter

- Connects Comfort Flo® High Flow Nasal Cannula to standard oxygen tubing during step-down therapy adapter
- Packaged 50/cs



HUD1669 Monitoring Adapter

- 15 mm OD x 15 mm ID/ 22 mm OD connectors
- Standard luer fitting
- Can be used with resuscitation bags
- Packaged 50/cs



HUD1667 Ventilator Monitoring Adapter Circuit

- Allows monitoring of ventilator weaning parameters in intubated patients
- Oxygen nipple adapter, aerosol tee, 6-in. tube, 15 mm x 22 mm adapter and two one-way valves
- Packaged 50/cs



Minor Capital

Fintertip Pulse Oximeters	99
Suction Regulators	100-101
Air & Oxygen Flow Meters	102-103
Oxygen Regulators	104
Oxygen Analyzers	105
Suction Aspirators	106

Fingertip Pulse Oximeters

- Get quick and precise SpO₂ and pulse readings in seconds
- Small and lightweight clip-on sensors automatically take readings when a finger is inserted
- Soft finger pad to enhance patient comfort
- Fit a wide range of finger sizes
- All models include AAA batteries, lanyard and warranty



Bluetooth Finger Pulse Oximeter

- Connects to Medline Health app via Bluetooth for quick readings of SPO₂ and pulse rate
- Includes carrying case
- OTC for everyday use
- One-year warranty

HCSM70T, 1 ea



OTC Fingertip Pulse Oximeter

- Small, lightweight finger clip with LCD display designed to quickly read SPO₂ and patient heart rate
- Complete with lanyard, instructions and 2 AAA batteries with 24-hour continuous use
- No prescription required to purchase

HCSM70R, 1 ea



Soft-Touch Finger Pulse Oximeter

- RX only for healthcare use
- One-year warranty

HCSM70S, 1 ea



Pediatric Finger Pulse Oximeter

- OTC for everyday use
- One-year warranty

HCSM70P, 1 ea

Note: For more Pulse Oximeter offerings please reach out to your Medline Representative.

Suction Regulators

The choice is clear

With a unique transparent and disposable back, Medline Regulators enhance efficiency, reduce costs and improve patient safety.

Identify contamination instantly

Up to 37% of suction regulators in hospitals are contaminated due to fluid breaches, and pathogens can spread back through the suction line within 48 hours.¹ The Medline Suction Regulator has a clear back plate offering a view of the internal components of the regulator, so it's easy to see when it's contaminated.

1. Suction Regulators: A Potential Vector for Hospital Acquired Pathogens by Keith S. Kaye, MD, MPH; Dror Marchaim, MD; Chester Smialowicz, MD; Lauren Bentley, MSBME



There's no time for downtime

The clear back plate provides an internal view of the regulator, making it possible to see fluid breaches. Plus, cleaning is simple: the unit is easily taken apart, allowing you to discard the contaminated back plate and replace with a new component in minutes.



Step 1:
Remove contaminated Regulator



Step 2:
Dispose



Step 3:
Reassemble with only 1 NEW component in just minutes



Clear back plate

Infection Control Measure

Remove from the outlet and identify contamination immediately

For specialty oxygen regulators (EMS, all-brass, etc.) please contact your **Medline Representative** or email ProductSupportRespiratory@Medline.com with your request.

Enhance efficiency, reduce costs and improve patient safety.

Continuous Models

- Continuous 3-mode (0-300mmHg)
- Pediatric 2-mode (0-160mmHg)
- Continuous high/surgical (0-760mmHg)

Continuous/Intermittent Models

- Adult continuous/intermittent (0-300mmHg)
- Pediatric continuous/intermittent (0-160mmHg)
- Timing cycle: 16 seconds on/8 seconds off
- Can be adjusted on the back of the unit



Analog




Digital



Suction regulator ordering information

RTXXXXXX



Vacuum Traps

Item No.	Description
RTTRAPMT	1/8" Male NPT Fitting
RTTRAPDH	Female DISS Handtight
RTTRAPDN	Female DISS Nut/Nipple



Vacuum/Safety Traps

Patient Fittings



DISS Male



Hose Barb

Examples of Wall Adapters



Chemetron



DISS Handtight



Ohmeda



Puritan-Bennett

For MRI compatible please add "-M" to the end of the part number (available for adult/analog units only).
For color coding, please specify in sales notes. Default color is white.

Air and Oxygen Flowmeters

5 YEAR warranty

- Complete line of flowmeters from neonatal to adult
- Available in standard Thorpe style and Dial style
- MRI options are available
- One of the smallest and lightest brass bodied flowmeters available
- Color coded inner tube and labeling



Dial/Click Style Flowmeter



Dual/Y-block Flowmeter Assembly



Thorpe Style



Swivel

Accessories



Swivels/Christmas Trees
HCS4700G (Green)
HCS4700Y (Yellow)



Power Takeoff

Examples of wall adapters




Chemetron



DISS Handtight




Ohmeda



Puritan-Bennett

Flowmeter ordering information

RTFMXXXXX(YY)*



For dual Y-block Flowmeter configurations: Replace RTFM prefix with RFMD

Medical Gas Hose Assemblies

Medline offers custom-configured medical gas hose assemblies in four easy steps:

- ¼" ID Non-Conductive
1. Select hose type
 2. Select the inlet fitting (gas service and type)
 3. Select the outlet fitting (gas service and type)
 4. Specify the hose length

RT = Hose Assembly

1. Hose Type

OXY = Oxygen	N2 = Nitrogen
VAC = Vacuum	CO ₂ = Carbon Dioxide
N2O = Nitrous Oxide	WAG = Waste Anesthetic Gas Disposal
AIR = Medical Air	

2. Inlet Fitting (Connects to the Source)

OM: OHIO Male	OF: OHIO Female Coupler
DC: DISS Male with Check	DN: DISS Nut/Nipple Female
CM: Chemetron Male	DH: DISS Handtight Female
PB: Puritan-Bennett Male Coupler	CF: Chemetron Female Coupler
	PF: Puritan-Bennett Female

3. Outlet Fitting (Connects to the Medical Device)


OM: OHIO Male	OF: OHIO Female Coupler
DC: DISS Male with Check	DN: DISS Nut/Nipple Female
CM: Chemetron Male	DH: DISS Handtight Female
PB: Puritan-Bennett Male Coupler	CF: Chemetron Female Coupler
	PF: Puritan-Bennett Female

4. Hose Length, in feet (e.g. 2 = 2-ft., 4 = 4-ft., 7 = 7-ft., etc.)

Example:

RTOXYDHDNC5 ¼" ID non-conductive oxygen hose assembly with a DISS handtight Female inlet and DISS Male outlet, 5' in length.

For custom configurations, contact your Medline Representative.



Oxygen Regulators

- For all medical and emergency applications
- Fitted with brass sleeves around key parts
- Laser-etched knob and body

Economy

Item No.	Description	Flow Range	Flow Increments	Pkg.
HCS8708M	870 CGA Connection	0–8 LPM	0.12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 8	1 ea
HCS8715M	870 CGA Connection	0–15 LPM	0.12, 0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15	1 ea
HCS8725M	870 CGA Connection	0–25 LPM	0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25	1 ea
HCS5408M	540 CGA Connection	0–8 LPM	0.12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 8	1 ea
HCS5415M	540 CGA Connection	0–15 LPM	0.12, 0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15	1 ea



HCS8715M

HCS5415M

Super Mini

Item No.	Description	Flow Range	Flow Increments	Pkg.
HCS8715MA	870 Connection	0–15 LPM	0.12, 0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15	1 ea



HCS8715MA

Accessories

Item No.	Description	Pkg.
HCS870W	Oxygen Cylinder Wrench, Small (CGA-870)	10/bg
HCS870WC	Oxygen Cylinder Wrench with Security Chain, Small (CGA-870)	10/bg
HCS870WLG	Oxygen Cylinder Wrench, Large (CGA-870)	25/bg
HCS166840	Brass Yoke Washer Seal for CGA-870 Oxygen Regulator	50/bg

EMS Regulators

- Easily transport ventilator patients
- Maximum inlet pressure: 3,000 PSI
- Chrome plated brass body
- Solid brass core

Item No.	50 PSI Power Take-off	Hose Barb	Liters Per Minute	Cylinder Style	Pin Index	Pkg.
CPRB34501C	1			CGA346	Air	1 ea
CPRB87251B	0	●	0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25	CGA870	Oxygen	1 ea
CPRB87251B2C	2	●	0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25	CGA870	Oxygen	1 ea
CPRB54251B2C	2	●	0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25	CGA540	Oxygen	1 ea
CPRB54402C	2			CGA540	Oxygen	1 ea
CPRB54501C	1			CGA540	Oxygen	1 ea



CPRB87251B

For specialty gas delivery products including blenders, analyzers, and more, please contact your **Medline Representative** or email ProductSupportRespiratory@Medline.com with your request.

Oxygen Analyzers and Sensors

Maxtec MaxO₂ ME Smart Alarms

Visual and audible smart alarms that can be independently adjusted to alert caregivers of a possible malfunction in the O₂ delivery system, which can provide additional time to quickly troubleshoot to avoid interruption in care.



CERR230P01

Oxygen Analyzers

- Easy, one-step calibration
- 2-year warranty

Item No.	Description	Pkg.
CERR230P01	MaxO ₂ ME Oxygen Monitor	1 ea
HCSR217P72	MaxO ₂ +AE Oxygen Analyzer	1 ea
HCSR217P62	MaxO ₂ Plus A Oxygen Analyzer	1 ea
CERR218P12	Handi Oxygen Analyzer	1 ea
CERR221P11	UltraMaxO ₂ Oxygen Analyzer	1 ea



HCSR217P62



CERR221P11

Sensors

Item No.	Description	Pkg.
CERR125P03002	Max-250E Sensor	1 ea
CERR140P02	Max-550E Sensor	1 ea
CERR113P10	Max-11 Sensor, FIO ₂	1 ea
CERR115P10	Max-13 Oxygen Sensor	1 ea
CERR109P45	Max-12 Oxygen Sensor	1 ea
HCSR125P02011	Oxygen Analyzer Sensor for HCSR217P6	1 ea
HCSR125P07	Max-13-250 Oxygen Sensor	1 ea
CERR107P20	Max-2 FIO ₂ Sensor	1 ea
CERR113P90	Max-1 FIO ₂ Sensor	1 ea
CERR109P48	Max-12A Oxygen Sensor	1 ea
CERR109P53	Max-12C Oxygen Sensor	1 ea
HCSR114P70	Max-16 Oxygen Sensor	1 ea



CERR125P03002



CERR140P02



HCSR114P70

Suction Aspirator and Disposable Canister Kits

- Easy-to-clean, simple-to-use portable suction aspirators are designed to withstand rigors of home care
- Devices feature convenient electrical cord storage, handles that help prevent breakage, and ultra-quiet motors
- Disposable 850 cc or 1,200 cc suction canister sets with float lids and 850 cc suction canister sets with filter lids
- HCS7121, HCS7851, and HCS7852 kits include 6" and 18" long suction tubing



HCS7001



HCS7000

Portable Suction Aspirators

Item No.	Description	Pkg.
HCS7000	Vac-assist Suction Aspirator with 800 cc Canister	1 ea
HCS7001	Vacu-Line Suction Aspirator with 850 cc Canister	1 ea

Canisters with Float Lids

Item No.	Description	Pkg.
HCS7850	Canister with Float Lid, 850 cc	12/cs
HCS7851	Canister with Float Lid and Tubing, 850 cc	12/cs
HCS7120	Canister with Float Lid, 1,200 cc	12/cs
HCS7121	Canister with Float Lid and Tubing, 1,200 cc	12/cs



Float Lid:
HCS7851



HCS7120

Canisters with Filter Lids

Item No.	Description	Pkg.
HCS7852	Canister with Filter Lid and Tubing, 850 cc	50/cs
HCS7853	Canister with Filter Lid, 850 cc	100/cs



Filter Lid:
HCS7852



HCS7850

Accessories

Item No.	Description	Pkg.
HCS7001SET	Replacement Suction Filter and Tubing Kit	1 ea
HCS4053	Hydrophobic Filter Set, Barb Both Ends	3/pk
HCS70131	Replacement Gauge Assembly, Includes Pressure Gauge, Pressure Adjustment Knob and Barb Tubing Adapter	1/ea



Tubing included with
HCS7851, HCS7852, and HCS7121

Index

A

Adapters:	
Elbow	47
Adapters and Connectors	88
Addipak Unit	
Dose Solutions	49
Aerosol:	
Masks	60, 78
AG Cuffill	44
Air Flowmeters	102
Aspirator, Suction	106
AsthmaCheck	
Peak Flow Meter	85
AsthmaMD	
Peak Flow Meter	85

B

Babi.Plus	
Bubble Pap Valve	27
Bacterial Filters	11
Bags,	
Resuscitation	37-39
Blender Buddy	
for MicroMax	32
Breathing Circuits	15
Bronchial Hygiene	83
Bubble Humidifiers	70
Bypass Heat and	
Moisture Exchanger	12

C

Canister Kits,	
Disposable	106
Cannulas:	
CO2	67-68
Catheters:	
Closed Suction	45
Open Suction	48
Sleeved Suction	48
Suction	51-52
CO2:	
Cannulas	67-68
Detectors, Colorimetric	40
Masks	69
Comfort Flo	
Humidification Systems	30
Comfort Flo Plus	29
Compressors	
Nebulizer	81-82
Concha Reservoirs	6
ConchaSmart and	
Concha Columns	6

CPAP:

Mask, Nasal	23
Masks	22-24
Full-Face	23
Systems	25-26

D

Detectors,	
Colorimetric CO2	40
Disposable Canister Kits	106

E

Elbow Adapter	47
Endotracheal	45, 47
Tubes	54
Exchangers:	
Heat and Moisture	13
Bypass	12

F

Filtered Heat and Moisture	
Exchangers	13
Filters:	
Bacterial/Viral	11
Finger Pulse Oximeter	99
Flowmeters:	
Air & Oxygen	102

G

Gas:	
Hose Assemblies	103

H

Heat and Moisture:	
Exchangers	13
Filtered	13
Heated High Flow	
Accessories	34
Hose Assemblies, Gas	103
Humidifiers, Bubble	70
Hyperinflation Systems	42

I

Infant Nasal Prong	
CPAP Cannulas	27
Interfaces, Non-Invasive	19
ISO-Gard Circuits	9
IV Pole Basket	34

M

Manual:	
Resuscitators	36-39
Masks:	
Aerosol	60, 78
CO2	69
CPAP	22-24
Full-Face	23
Nasal Gel	23
Full-Face	22
Oxygen	57-58
Pediatric Gel	23
Resuscitation	41
Maxtec MaxO2 ME	105
MaxVenturi	31
Medline Peak Flow Meters	86
MicroMax High Flow	32
MicroMax Low Flow	33
MRI-Conditional Transport	
Ventilator Circuit	16
Mucus Traps	55

N

Nebulizers:	
Compressors	81-82
Neptune Heated	
Humidifier	5
Non-Invasive	
Interfaces	19
Non-Invasive Ventilation (NIV):	
Adult Full Face Masks	20
Circuits	21

O

Open Suction Catheters	48
OptiChamber Diamond	79
Oximeters:	
Finger Pulse	99
Oxygen:	
Flowmeters	102
Masks	57-58
Oxygen Supply Tubing	66
OxyMask Aerosol	58
OxyMask ETCO2	69

P

Personal Best® Portable	
Peak Flow Meter	85
Pulse Oximeters	
Finger	99

R

Regulators	104
Resuscitation:	35
Accessories	41
Bags	37-39
Masks	41
Resuscitators:	
Manual	36-39

S

Sleeved Suction Catheters	48
Smart Stack Pole Stand	34
Solutions, Sterile	52
Specimen Traps	55
Suction:	
Aspirator	106
Catheters	45, 48, 51, 52
Suction Regulators	100
SuperSoft High	
Flow Oxygen Cannulas	65
SuperSoft Oxygen	
Cannulas	62

T

Tracheostomy	45, 47
Care	50
Trach-Flex Plus	29
Traps	
Mucus	55
Specimen	55
Trays	
Tracheostomy Care	50
Tubing	
Endotracheal	54

V

Ventilator:	
Circuit	16

Y

Y-Connector	47
-------------	----



**We make
healthcare
run better™**

Medline Industries, LP
 Three Lakes Drive, Northfield, IL 60093
 Medline United States | 1-800-MEDLINE (633-5463)
 medline.com | info@medline.com

Medline Canada
 1-800-268-2848 | medline.ca | canada@medline.com
Medline México
 01-800-831-0898 | medlinemexico.com | mexico@medline.com

Follow us    

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure. © 2023 Medline Industries, LP. IntelliCuff is the trademark of Hamilton Medical. Respironics and Sidestream are trademarks of Koninklijke Philips N.V. and its affiliates. All rights reserved. LiteTouch Respironics and LiteTouch are trademarks of Koninklijke Philips N.V. and its affiliates. All rights reserved. OptiChamber Respironics and OptiChamber are trademarks of Koninklijke Philips N.V. and its affiliates. All rights reserved. Respironics and Personal Best are trademarks of Koninklijke Philips N.V. and its affiliates, all rights reserved. Bleep, Dreamports, DreamWay, and DreamPreps are trademarks of Bleep LLC. AIRgel, Ascend, VeraSeal 2, Mojo, MiniMe, MiniMe 2, and Innova are trademarks of Sleepnet. AIRgel, Ascend, VeraSeal 2, Mojo, MiniMe, MiniMe 2, and Innova are trademarks of Sleepnet. Citrus II CPAP Mask Wipes is a trademark of Beaumont Products Inc. Go-PAP and O₂-Max are trademarks of Pulmodyne. FloCap is a trademark of Sunmed Medical. The Oxymask trademark is the property of Southmedic Inc., and is used under license by Medline Industries, LP. MKT19W4807401 / LIT085 / 55