

# SOL-CARE™ Safety Blood Collection Needle

with or without pre-attached holder



Discover our range of safety 'butterfly' needles for blood collection. Each device features a simple safety mechanism that should be activated immediately after blood withdrawal is complete, thereby reducing the risk of needlestick injury.

- Thin walled needle minimises patient discomfort
- Flexible wings comfortably secure device to patient's skin
- DEHP-free tubing is soft, flexible and kink-resistant
- Transparent design allows for clear visualisation of blood flow

**Simple, intuitive safety system  
minimises risk of needlestick injury:**

- Three simple ways to activate safety mechanism – using either one or two handed technique
- Audible click tells you when the needle is safely locked within the safety shield

With pre-attached holder



### SOL-CARE™ Safety Blood Collection Needle

Product configuration	Reference number current/ new*	Wing color	Needle outer diameter		Needle length		Tubing length	
Luer Adapter	110201030035	Orange	25G	0.5 mm	3/4"	19 mm	12"	305 mm
	110201030033	Blue	23G	0.6 mm	3/4"	19 mm	12"	305 mm
	110201030031	Green	21G	0.8 mm	3/4"	19 mm	12"	305 mm

Sales Unit: 50 pcs/Box; 200 pcs/Case

### SOL-CARE™ Safety Blood Collection Needle with Pre-attached Holder

Product configuration	Reference number	Wing color	Needle outer diameter		Needle length		Tubing length	
With pre-attached holder	110201030010	Orange	25G	0.5 mm	3/4"	19 mm	12"	305 mm
	110201030006	Blue	23G	0.6 mm	3/4"	19 mm	12"	305 mm
	110201030002	Green	21G	0.8 mm	3/4"	19 mm	12"	305 mm
	110201030012	Orange	25G	0.5 mm	3/4"	19 mm	7"	178 mm
	110201030008	Blue	23G	0.6 mm	3/4"	19 mm	7"	178 mm
	110201030004	Green	21G	0.8 mm	3/4"	19 mm	7"	178 mm

Sales Unit: 25 pcs/Box; 200 pcs/Case



CORP-BC-115-0813

#### For more information, please contact:

Sol-Millennium Medical Inc.  
 315 Shawnee North Drive, Suite 100,  
 Suwanee, GA 30024  
 T: 404-973-2200  
 e-mail: info@sol-m.com

[sol-m.com](http://sol-m.com)