


SOL-CARE™ Safety Winged Infusion Set





Discover our range of safety winged sets for short-term intravascular administration of medication.

Each device features a simple safety mechanism that should be activated immediately after infusion is completed, thereby reducing the risk of needlestick injury.

- Thin walled needle minimizes patient discomfort
- Flexible wings comfortably secure device to patient's skin
- DEHP-free tubing is soft, flexible and kink-resistant



SOL-CARE™ Safety Winged Infusion Set



Simple, intuitive safety system minimises risk of needlestick injury:

- Three simple ways to activate safety mechanism – using either one or two handed technique
- Audible click tells you when the needle is safely locked within the safety shield

CE 0123

FDA 501k Cleared

SOL-CARE™ Safety Winged Infusion Set

Product code	Color coding	Needle outer diameter		Needle length		Tubing length	
		Gauge	mm	Inches	mm	Inches	mm
193412SI	Cream	19G	1.08	3/4	19	12	305
213412SI	Deep Green	21G	0.8	3/4"	19	12	305
223412SI	Black	22G	0.7	3/4"	19	12	305
233412SI	Light Blue	23G	0.6	3/4"	19	12	305
253412SI	Orange	25G	0.5	3/4"	19	12	305

CORP-SIN-130-0920

For more information about SOL-CARE™ Safety Winged Infusion Set, please contact your local SOL-MILLENNIUM® representative.

Sol-Millennium Medical Inc.
 315 Shawnee North Drive, Suite 100, 30024 Suwanee, GA
 Tel: 404-973-2200
 info@sol-m.com
 sol-m.com

