



SOL-MILLENNIUM®  
Building a Healthier Tomorrow



Product Catalog

**Covid-19**

North America

## A FEW WORDS ABOUT US

SOL-MILLENNIUM Medical Group is a relatively new enterprise at just 10 years old, but is quickly building a reputation as a trusted, reliable and customer-focused partner to the healthcare industry.

Whilst we may be young, our reach is wide with presence across the globe. You will find us in Europe, North America, South America and SE Asia.



Our approach to product development is simple. We are particularly careful to ensure that the design of our products are intuitive and easy to use, thereby minimising the impact on application technique.

Our team is experienced and knowledgeable and we are here to help you operate efficiently, effectively and economically. Our new R&D facility in Switzerland has been increasing and enhancing our current comprehensive range of products.

At SOL-MILLENNIUM we know what we do best – we make good quality, affordable medical devices that are familiar, safe and easy to use.

# Table of contents

Medical Masks	06
Surgical Face Mask	07
Medical Surgical Mask, TieBack	08
Particulate Respirator KN95	09
Face Shield	10
Isolation Gowns	11
Wipes	12
Instant Hand Sanitizer	13
Alcohol Swabs	14
Infrared Thermometers	15
Fingertip Pulse Oximeter	16
Blood Collection products	17
Vaccination products	22



# BUILDING A HEALTHIER TOMORROW

It is a promise we made to our customers and society

SOL-MILLENNIUM® focuses on developing high quality and affordable safety medical devices that support best clinical practice – whilst also protecting healthcare workers and their patients from the risks associated to pathogens exposure.

Fully-automated production lines and technology lead to significant efficiency and cost reduction



## RESEARCH & DEVELOPMENT

Innovation runs through our veins – it has been our inspiration since the very start. We are passionate about the details and applying precision engineering techniques to produce ever better designs (we've patented over 100 so far).

Our R&D Centres are equipped with the latest design and laboratory technology, where our experienced team have the freedom to explore, imagine and create.


Close collaboration with frontline healthcare professionals and customers is critical to the way we work – by sharing knowledge and ideas together we ensure that our innovation meets everyday clinical needs.

## QUALITY CONTROL

Each and every product we manufacture meets internationally recognized and required standards of Quality Management

Passion, expertise and personal responsibility lie at the heart of our Quality Policy

A culture of continuous improvement means we never stop learning and advancing our know-how



We sell our products and services to hospitals, veterinary groups, medical supply companies, government health departments and scientific research facilities around the world.

## GLOBAL TERRITORIES: 500+ CUSTOMERS TODAY IN 40 COUNTRIES

### NORTH AMERICA

Canada  
USA  
Puerto Rico

### CENTRAL

Jamaica **AMERICA**  
Barbados  
Costa Rica  
Panama  
Trinidad  
Dominican Republic

### SOUTH AMERICA


Brazil

### ASIA-PACIFIC

Australia  
Taiwan  
Hong Kong  
Singapore  
Malaysia  
China (Continental)

### EUROPE, AFRICA AND MIDDLE EAST

UK	Austria	Romania
Ireland	Switzerland	Morocco
Finland	Poland	Saudi
Norway	Latvia	UAE
Sweden	Spain	Oman
Netherlands	Portugal	Kuwait
Belgium	Italy	Qatar



Headquarter: Chicago



Eight global sites



Employees based in 12 Countries




## OUR GLOBAL SUPPLY CHAIN



# Disposable Medical Procedure Face Mask (non-sterile)



- **Consists of a mask body, ear loops and a nose clip**  
The product is in compliance with EN 14683:2019 + AC Type I
- **Three layers:** outer layers are made of non-woven fabric and the medium layer is made of melt-blown fabric
- Bacterial filtration efficiency  $\geq 95\%$
- Intended to be worn by patients and other persons to reduce the risk of spread of infections
- Does not contain components made with natural rubber latex




Reference number	Case dimension	Gross weight
MSK003KU	530 x 380 x 400 mm	9.5 kg
MSK003X	520 x 380 x 300 mm	7.6 kg

**Sales Unit:** 50 pcs/Box; 2000 pcs/Case

\* Subject to availability in your country

# Disposable Surgical Face Mask (non-sterile)



Medical  
Supplies



- **Consists of a mask body, ear loops and a nose clip**

Complies with EN 14683:2019 Type IIR

- **Three layers:** outer layers are made of non-woven fabric and the medium layer is made of melt-blown fabric

- Bacterial filtration efficiency  $\geq 98\%$

- Mask is intended to be used by operation room personnel and other general healthcare workers to protect both patients and healthcare workers against transfer of microorganisms, blood and body fluids and airborne particulates

- Does not contain components made with natural rubber latex




Reference number	Case dimension	Gross weight
SMSK001	540 x 400 x 320 mm	8 kg
SMSK001KU	530 x 380 x 400 mm	9.5 kg

**Sales Unit:** 50 pcs/Box; 2000 pcs/Case

\* Subject to availability in your country

# Medical Surgical Mask, TieBack (non sterile)



Medical  
Supplies



- **Consists of a mask body, tieback and a nose clip**

Complies with EN 14683:2019 Type IIR

- **Three layers:** outer layers are made of non-woven fabric and the medium layer is made of melt-blown fabric


- Bacterial filtration efficiency  $\geq 98\%$

- Microbial cleanliness  $\leq 30$  CFU/g

- Splash resistance  $\geq 16$  kPa

- Mask is intended to be used by operation room personnel and other general healthcare workers to protect both patients and healthcare workers against transfer of microorganisms, blood and body fluids and airborne particulates

- Does not contain components made with natural rubber latex



**Ref. number:** SMSKT001

**Sales Unit:** 50 pcs/Box; 2000 pcs/Case

**Case dimension:** 520 x 380 x 440 mm

**Gross weight:** 7.6 kg

\* Subject to availability in your country



# Particulate Respirator KN95



## Consists of straps, Filter, Nose clip and face seal

- In compliance with norm NIOSH – 42 CFR 84
- Polypropylene non-woven fabric
- Particulate Filtration Efficiency  $\geq 95\%$
- Designed to help provide reliable respiratory protection for the wearer with 95 percent filtration efficiency against particulate aerosols free of oils.
- Does not contain components made with natural rubber latex



**Ref. number:** 1200H


**Sales Unit:** 20 pcs/Box; 480 pcs/Case

**Case dimension:** 525\*415\*465mm

**Gross weight:** 8 kg

\* Subject to availability in your country

# Face Shield (Safety Face Mask)



Medical  
Supplies



- **Consists of a protective cover made by polymer materials, elastic band and foam**

Complies with EN 166  
Disposable device

- **Face shield in PET 0.2-0.3 mm thick:**

Anti-spitting splash  
Soft and fitting  
High transparency

- Refractive powers EN 166: Optical Class 1  
Transmittance EN 167 > 74.4%  
Luminous Transmission ANSI/ISEA Z87.1-2015: > 85%

- Protect eyes and block body fluids, blood sprays or other splashes

- Does not contain components made with natural rubber latex



**Ref. number:** FS9132

**Sales Unit:** 50 pcs/Box; 2000 pcs/Case

**Case dimension:** 750 x 330 x 510 mm

**Gross weight:** 13.5 kg

\* Subject to availability in your country

# Isolation Gown (non-sterile)



Medical  
Supplies



## Basic coverage/fluid repellent

In compliance with ANSI/AAMI PB70 Level 2

Made of comfortable, strong and durable material, with overhead, waist belt, elastic cuff and full back

## Tri-laminate Spunbond - Meltblown - Spunbond, polypropylene fabric (SMS) 23 gsm

To protect against a risk of splashes from potentially infectious liquid and solid material, alcohol, blood, water, body fluids and dust

Does not contain components made with natural rubber latex



Reference number	Size	Case dimension	Gross weight
IG004	Large	485 x 270 x 320 mm	6.7 kg
IG004XL	Extra Large	485 x 270 x 320 mm	6.7 kg

**Sales Unit:** 10 units/ Box; 100 pcs/Case

\* Subject to availability in your country

# Alcohol Wipes (Tube)



Medical  
Supplies



- Saturated with 70% isopropyl alcohol (CAS 67-63-0)
- **It is a convenient multipurpose cleaner and general disinfectant** - ideal for most medical surfaces, devices, the offices and around the home
- Product can be used within 2 months after opening
- Does not contain components made with natural rubber latex



## Ref. numbers:

WTIA0001  
Practical dispenser tube (Ø85x150mm) contains 110 pcs of nonwoven wipes (200x130mm)

WTIA0003  
Practical dispenser tube (Ø100 x 165mm) contains 75 pcs of nonwoven wipes (420x145mm)

\* Subject to availability in your country

# SOL-M™ Instant Hand Sanitizer



Medical  
Supplies



- Ethyl Alcohol 75% Gel with moisturizers and vitamin E for antiseptic purpose (CAS: 64-17-5)
- **It is a convenient antiseptic and general disinfectant**  
The 236 and 500 ml volumes have pump delivery.
- Intended as antiseptic, for hand washing to decrease bacteria on the skin.  
Recommended for repeated use. Not recommended for infants
- Does not contain components made with natural rubber latex



Reference number	NDC Number (bottle)	Description	Case dimension	Gross weight	Case units
HSG59HU	73720-0059-1	Gel, 59 ml	250 x 195 x 220 mm	5.03 kg	72
HSG236HU	73720-0236-1	Gel, 236 ml	310 x 270 x 175 mm	6.14 kg	24
HSG500HU	73720-500-1	Gel, 500 ml	295 x 220 x 200 mm	6.29 kg	12

\* Subject to availability in your country

# Alcohol Swabs (sterile and non-sterile)




- Complies with USP standard
- Products are made of non-woven fabrics, saturated with 70% IPA. Individually foil wrapped for convenience
- **Four-layer wrapper** provides air-tight seal which prevents leakage and drying out. Softer and Thicker
- Does not contain components made with natural rubber latex



Reference number	Description	Sterility	Box/Case units
SIA1001	Alcohol Swabs Sterile Medium	sterile	100 / 10 000
SIA1002	Alcohol Swabs Sterile Large	sterile	100 / 1 000
SIA2001	Alcohol Swabs Non-Sterile Medium	non-sterile	100 / 10 000
SIA2002	Alcohol Swabs Non-Sterile Large	non-sterile	100 / 1 000

\* Subject to availability in your country



- CE marked and China SFDA Class 2 registered medical device
- Intended to measure the forehead temperature of human body
- Measurement range: 32.0°C~43.0°C (89.6°F~109.4°F)  
Accuracy:  $\pm 0.2^{\circ}\text{C}$   
Precision readable: 0.1°C
- **Response time: 1 second. Measurement distance: 0-5cm**



## SFDA Class 2 Medical Device

**Ref. number:** YT-1

**Sales Unit:** 40 pcs/Case

**Case dimension:** 575x 351 x 174 mm

**Gross weight:** 6.2 kg

Does not contain components made with natural rubber latex

\* Subject to availability in your country

# Fingertip Pulse Oximeter



Medical  
Supplies



- Intended for measuring the Pulse Oxygen Saturation, Pulse Rate, and Perfusion Index through patient's finger. It's applicable for spot-checking SpO2 and pulse rate of adult and pediatric at home and medical clinics
- SpO2 measuring range: 35%- 100% (accuracy: 2% from 70% to 100%)  
Pulse rate range 35bpm - 240 bpm (accuracy 2bpm / 2% whichever is greater)  
Perfusion index range 0%- 20%

## Features:

- Dual wavelength LED sensor transducer
- OLED displays SpO2, PR, PI, Pulse bar and waveform
- Special design with splash-proof and drop-resistant
- 4 direction display with big numbers
- Audible and visual alarm
- Low power consumption with 2 AAA alkaline batteries



FDA cleared



**Ref. number:** PC60B1

Does not contain components made with natural rubber latex

\* Subject to availability in your country



# SOL-CARE™ Safety Blood Collection Tube Holder Set



Blood Collection  
Devices



**SOL-CARE™ Safety Blood Collection Tube Holder Set** is a high-quality and high-tech blood collection device, designed to improve the performance of healthcare workers, while also offering maximum safety.



### Holder

Compatible with:

- standard and flashback multi-sample needles.
- standard blood collection tubes.



### Vacuum Activation Tube

Retraction technology to capture the needle after blood has been collected.

## FDA cleared

● **Minimum training requirement** - "Using the yellow activation tube is like adding a blood collection tube".

● **In-vein retraction** - effectively reduces the risk of needlestick injuries and blood exposure.

● **Speed of needle retraction** - lower risk of vein damage.

● Does not contain components made with natural rubber latex

**Ref. number:** 7001BCD

**Sales Unit:** 50 pcs/Box; 400 pcs/Case

\* Subject to availability in your country

# SOL-M™ Multi-Sample Needle



Blood Collection  
Devices



**SOL-M™ Multi-Sample Needle** makes it possible to collect several samples with a single vein puncture.

- **Low-angle bevel design allows easy penetration**

- **Color-coded caps and hubs** allow easy identification of needle gauge size

- Compatible with all standard holders

- Does not contain components made with natural rubber latex



FDA cleared

## SOL-M™ Multi-Sample Needle

Reference number	Needle outer diameter		Needle length	
	Gauge	mm	Inches	mm
110201020002	22G	0.7 mm	1"	25 mm
110201020003	22G	0.7 mm	1 1/2"	38 mm
110201020004	21G	0.8 mm	1"	25 mm
110201020005	21G	0.8 mm	1 1/2"	38 mm
110201020006	20G	0.9 mm	1"	25 mm
110201020007	20G	0.9 mm	1 1/2"	38 mm

**Sales Unit:** 100 pcs/Box; 1000 pcs/Case

\* Subject to availability in your country



## SOL-CARE™ Safety Multi-Sample Blood Collection Needle\*

is a safety-engineered tool for multi-sample venous blood collection. The device is designed with a safety mechanism covering the needle after use. In the activated position, the needle cover guards against accidental needlestick during normal handling and disposal.



FDA cleared

### SOL-CARE™ Safety Multi-Sample Needle

Reference number	Needle outer diameter		Needle length	
SMSN21125	22G	0.7 mm	1 1/4"	32 mm
SMSN22125	21G	0.8 mm	1 1/4"	32 mm

**Sales Unit:** 48 pcs/Box; 480 pcs/Case

FDA cleared

### SOL-CARE™ Safety Multi-Sample Needle with Pre-attached Holder

Reference number	Needle outer diameter		Needle length	
SMSNH22125	22G	0.7 mm	1 1/4"	32 mm
SMSNH21125	21G	0.8 mm	1 1/4"	32 mm

**Sales Unit:** 50 pcs/Box; 400 pcs/Case

Does not contain components made with natural rubber latex

\* Subject to availability in your country



## SOL-CARE™ Safety Blood Collection Needle

is a safety-engineered device intended for venous blood collection. The safety mechanism can be activated immediately after blood draw, helping to protect you against needlestick injury. It is also offered with a pre-attached holder.



FDA cleared

### SOL-CARE™ Safety Blood Collection Needle

Product configuration	Reference number current/ new*	Wing color	Needle outer diameter		Needle length		Tubing length	
Luer Adapter	110201030035	Orange	25G	0.5 mm	3/4"	19 mm	12"	305 mm
	110201030033	Blue	23G	0.6 mm	3/4"	19 mm	12"	305 mm
	110201030031	Green	21G	0.8 mm	3/4"	19 mm	12"	305 mm

**Sales Unit:** 50 pcs/Box; 200 pcs/Case

\* New reference number will be launched during 2020

### SOL-CARE™ Safety Blood Collection Needle with Pre-attached Holder

Product configuration	Reference number	Wing color	Needle outer diameter		Needle length		Tubing length	
With pre-attached holder	110201030010	Orange	25G	0.5 mm	3/4"	19 mm	12"	305 mm
	110201030006	Blue	23G	0.6 mm	3/4"	19 mm	12"	305 mm
	110201030002	Green	21G	0.8 mm	3/4"	19 mm	12"	305 mm
	110201030012	Orange	25G	0.5 mm	3/4"	19 mm	7"	178 mm
	110201030008	Blue	23G	0.6 mm	3/4"	19 mm	7"	178 mm
	110201030004	Green	21G	0.8 mm	3/4"	19 mm	7"	178 mm

**Sales Unit:** 25 pcs/Box; 200 pcs/Case

Does not contain components made with natural rubber latex

\* Subject to availability in your country



FDA cleared

**SOL-M™ Luer Adapter**

Reference number	Description	Box/Case (pcs)
LA21920	SOL-M™ Luer Adapter 20G	100/1000

FDA cleared



**SOL-M™ Blood Collection  
Tube Holder**



**SOL-M™ Tube Holder  
with Luer Lock Adapter\***



**SOL-M™ Blood Transfer Holder\***

Reference number	Description	Inner Bag/Case (pcs)
110201050001	SOL-M™ Blood Collection Tube Holder	250/1000
110201050003	SOL-M™ Tube Holder with Luer Lock Adapter*	200/800
BCTH002	SOL-M™ Blood Transfer Holder*	200/800 (Box/Case)

Does not contain components made with natural rubber latex

\* Subject to availability in your country

# SOL-CARE™ Safety Needle



Vaccination



**Designed for skin injection the SOL-CARE™ Safety Needle helps protecting from needlestick injuries after use.**

The SOL-CARE™ Syringe with Safety Needle combines SOL-M™ Syringe with SOL-CARE™ Safety Needle in one convenient solution.

Easy **one-hand technique** of activating protective arm.

**Audible click** indicates safety mechanism activation.

**Patented locking mechanism** ensures needle is secure and completely locked in protective arm.

**Two options** for protective arm activation: hard surface activation and thumb activation



Reference number	Product description	Gauge size	Outer diameter	Length		Shelf Box/ Shipping case (pcs)
SN2712	SOL-CARE™ Safety Needle 27G*1/2"	27G	0.4 mm	12.5 mm	1/2"	100/800
SN2612	SOL-CARE™ Safety Needle 26G* 1/2"	26G	0.45 mm	12.5 mm	1/2"	100/800
SN2558	SOL-CARE™ Safety Needle 25G*5/8"	25G	0.5 mm	16 mm	5/8"	100/800
SN2510	SOL-CARE™ Safety Needle 25G*1"	25G	0.5 mm	25 mm	1"	100/800
SN2515	SOL-CARE™ Safety Needle 25G*1 1/2"	25G	0.5 mm	40 mm	1 1/2"	100/800
SN2310	SOL-CARE™ Safety Needle 23G*1"	23G	0.6 mm	25 mm	1"	100/800
SN23125	SOL-CARE™ Safety Needle 23G*1 1/4"	23G	0.6 mm	30 mm	1 1/4"	100/800

Reference number	Product description	Gauge size	Outer diameter	Length		Shelf Box/ Shipping case (pcs)
12558SN	SOL-CARE 1ml Luer Lock Syringe	25G	0.5 mm	16 mm	5/8"	100/1200
32558SN	SOL-CARE 3ml Luer Lock Syringe	25G	0.5 mm	16 mm	5/8"	100/1200
32510SN	SOL-CARE 3ml Luer Lock Syringe	25G	0.5 mm	25 mm	1"	100/1200
32310SN	SOL-CARE 3ml Luer Lock Syringe	23G	0.6 mm	25 mm	1"	100/1200
32558SS	SOL-CARE 3ml Slip Tip Safety Syringe	25G	0.5 mm	16 mm	5/8"	100/1200
32510SS	SOL-CARE 3ml Slip Tip Safety Syringe	25G	0.5 mm	25 mm	1"	100/1200
32310SS	SOL-CARE 3ml Slip Tip Safety Syringe	23G	0.6 mm	25 mm	1"	100/1200

## FDA Cleared

Notes: Needle mounted on ref: 12558SN, 32310SN; needle tube length tolerances based on ISO 7864.

Does not contain components made with natural rubber latex

\* More sizes available on demand

# SOL-M™ Blunt Fill Needle

# SOL-M™ Hypodermic Needle



Vaccination



## SOL-M™ Blunt Fill Needle

designed to help reduce risk of needlestick injuries to healthcare workers when drawing up and preparing medication.

- Prevention from sharp injuries in hospitals and other health facilities.
- Eliminating unnecessary use of sharps.

**Ref. number:** 110022

SOL-M™ Blunt Fill Needle 18G\*1 1/2"

**Sales unit:** 100 pcs/Box; 1000 pcs/ Case



## SOL-M™ Hypodermic Needle

for injection and medication drawing up and preparation.

- Triple bevel cut and needle lubrication – minimizes patient discomfort.
- Color-coded hubs – easy gauge identification.
- Translucent polypropylene hubs – easy flashback confirmation.
- Luer Lock- and Slip Tip - compatible.

**FDA Cleared**



Reference number	Product description	Gauge size	Outer diameter	Length		Shelf Box/ Shipping case (pcs)
112705	SOL-M Hypodermic Needle 27G*1/2"	27G	0.4 mm	12.5 mm	1/2"	100/1000
112506	SOL-M Hypodermic Needle 25G*5/8"	25G	0.5 mm	16 mm	5/8"	100/1000
112510	SOL-M Hypodermic Needle 25G*1"	25G	0.5 mm	25 mm	1"	100/1000
112515	SOL-M Hypodermic Needle 25G*1 1/2"	25G	0.5 mm	40 mm	1 1/2"	100/1000
112310	SOL-M Hypodermic Needle 23G*1"	23G	0.6 mm	25 mm	1"	100/1000
111915	SOL-M Hypodermic Needle 19G*1 1/2"	19G	1.1 mm	40 mm	1 1/2"	100/1000
111910	SOL-M Hypodermic Needle 19G*1"	19G	1.1 mm	25 mm	1"	100/1000
111815	SOL-M Hypodermic Needle 18G*1 1/2"	18G	1.2 mm	40 mm	1 1/2"	100/1000
111810	SOL-M Hypodermic Needle 18G*1"	18G	1.2 mm	25 mm	1"	100/1000

Does not contain components made with natural rubber latex

\* More sizes available on demand

# SOL-CARE™ Safety Syringe

## SOL-CARE™ TB Safety Syringe




Vaccination



### SOL-CARE™ Safety Syringe

is a retracting needle manual safety syringe. It is designed to aid the prevention of needle stick injuries by manually retracting the needle into the barrel. A patented locking mechanism hold the plunger in position once fully retracted.

SOL-CARE™ Safety Syringe is available with exchangeable needle and without needle. It is compatible with standard Luer needles.



Reference number	Product description	Shelf Box/ Shipping case (pcs)
120006IM	SOL-CARE 3ml Luer Lock Safety Syringe w/o Needle	100/800
120007IM	SOL-CARE 5ml Luer Lock Safety Syringe w/o Needle	100/800

FDA Cleared

### SOL-CARE™ TB Safety Syringe with Fixed Needle

are designed to have a low dead space, to reduce medication waste.



Reference number	Product description	Gauge size	Outer diameter	Length		Shelf Box/ Shipping case (pcs)
100019IM	SOL-CARE 1ml TB Safety Syringe w/Fixed Needle	27G	0.4 mm	12.5 mm	1/2"	100/1000
100018IM	SOL-CARE 1ml TB Safety Syringe w/Fixed Needle	25G	0.5 mm	16 mm	5/8"	100/1000
100034IM	SOL-CARE 1ml TB Safety Syringe w/Fixed Needle	25G	0.5 mm	25 mm	1"	100/1000
100078IM	SOL-CARE 3ml Luer Lock Safety Syringe w/Exch Needle	25G	0.5 mm	16 mm	5/8"	100/800
100079IM	SOL-CARE 3ml Luer Lock Safety Syringe w/Exch Needle	25G	0.5 mm	25 mm	1"	100/800
100077IM	SOL-CARE 3ml Luer Lock Safety Syringe w/Exch Needle	23G	0.6 mm	25 mm	1"	100/800

FDA Cleared

100/800

Does not contain components made with natural rubber latex

\* More sizes available on demand



# SOL-GUARD™ TB Safety Syringe with Fixed Needle



**SOL-GUARD™** is a safety-engineered Syringes with Fixed Needle for injections. It is designed to reduce the risk of needlestick injuries to patients and healthcare workers before and during the injection process.

**Transparent Locking Safety Shield with Easy Grip Ridges** - protects during transportation

**Transport Position**

**Low Dead Space** - reduces medication waste



**FDA Cleared**

Reference number	Product description	Gauge size	Outer diameter	Length		Shelf Box/ Shipping case (pcs)
200018SG	SOL-GUARD 1ml TB Safety Syringe w/Fixed Needle 25G*5/8	25G	0.5 mm	16 mm	5/8"	100/800
200088SG	SOL-GUARD 1ml Safety Syringe Flu Tray 25G*1	25G	0.5 mm	25 mm	1"	25/1000

Does not contain components made with natural rubber latex

\* More sizes available on demand



## SOL-CARE™ and SOL-M™ Syringe Trays

offer syringes with needle conveniently packaged in bulk sterile fashion, in a rigid tray. This solution is ideal when multiple injection are planned in a series, like for a vaccination campaign.



**FDA Cleared**

SOL-CARE™ Safety Syringe Flu Tray						
Reference number	Capacity	Gauge size	Outer diameter (mm)	Length		Tray / Shipping Case (pcs)
100088IM	1 ml	25G	0.5 mm	25 mm	1"	25/1000

SOL-M™ Allergy and Mixing Syringe Trays						
Reference number	Capacity	Gauge size	Outer diameter (mm)	Length		Tray / Shipping Case (pcs)
181027T	1 ml	27G	0.4 mm	12.5 mm	1/2"	25/1000
181025T	1 ml	25G	0.5 mm	12.5 mm	1/2"	25/1000
181023T	1 ml	23G	0.6 mm	12.5 mm	1/2"	25/1000

SOL-M™ Luer Lock Syringe with Exchangeable Needle - Tray						
Reference number	Capacity	Gauge size	Outer diameter (mm)	Length		Shelf Box / Shipping Case (pcs)
1832510T	3 ml	25G	0.5 mm	25 mm	1"	25/1000
1832310T	3 ml	23G	0.6 mm	25 mm	1"	25/1000

Does not contain components made with natural rubber latex

\* More sizes available on demand

# SOL-M™ Syringe with needle



**A selection of 3 piece syringes** conveniently packaged with pre-assembled hypodermic needle for vaccination injection

**High-transparency barrel** - makes it easy to see syringe content

**Silicone lubrication** - ensures smooth plunger movement

**Safe plunger backstop** - prevents accidental withdrawal of plunger

**Large and stable plunger grip** - ensures stability and comfort during injection

**Triple bevel cut and needle lubrication** – minimizes patient discomfort

**Translucent polypropylene hubs** – easy flashback confirmation

**Luer Slip, Luer Lock options**



**FDA Cleared**

SOL-M™ Luer Lock Syringe with Exchangeable Needle

Reference number	Capacity	Gauge size	Outer diameter (mm)	Length		Shelf Box / Shipping Case (pcs)
1832558	3 ml	25G	0.5 mm	16 mm	5/8"	100/800
1832510	3 ml	25G	0.5 mm	25 mm	1"	100/800
1832310	3 ml	23G	0.6 mm	25 mm	1"	100/800

Does not contain components made with natural rubber latex

\* More sizes available on demand



**A selection of 3-piece syringes** to be used for vaccination injection together with SOL-CARE™ Safety Needle and SOL-M™ Hypodermic Needle

**High-transparency barrel** - makes it easy to see syringe content

**Silicone lubrication** - ensures smooth plunger movement

**Safe plunger backstop** - prevents accidental withdrawal of plunger

**Large and stable plunger grip** - ensures stability and comfort during injection

**Luer Slip, Luer Lock options**



**FDA Cleared**

SOL-M™ Luer Lock Syringe without Needle		
Reference number	Capacity	Shelf Box / Shipping Case (pcs)
P180001PP	1 ml	100/800
P180003	3 ml	100/800
P180005	5 ml	100/800

SOL-M™ Slip Tip Syringe without Needle		
Reference number	Capacity	Shelf Box / Shipping Case (pcs)
1800011	1 ml	100/800

Does not contain components made with natural rubber latex

\* More sizes available on demand

# SOL-M™ 1 ml Syringe with Dose Sparing feature



## Specially designed 1 ml syringe for accurate dose administration of medication

The „spine” molded on the gasket fits the Luer connection space, reducing the dead-volume and therefore the waste space.

The syringe can be used to inject vaccine with with SOL-CARE™ Safety Needle and SOL-M™ Hypodermic Needle.

**High-transparency barrel** - makes it easy to see syringe content

**Silicone lubrication** - ensures smooth plunger movement

**Large and stable plunger grip** - ensures stability and comfort during injection



**FDA Cleared**

SOL-M™ Slip Tip Syringe without Needle LDS			
Reference number	Tip	Capacity	Shelf Box / Shipping Case (pcs)
180011LDS*	concentric	1 ml	100/800

\* LDS: with dose sparing feature

Does not contain components made with natural rubber latex

\* More sizes available on demand



**SOL-MILLENNIUM**<sup>®</sup>  
Building a Healthier Tomorrow

**Sol-Millennium Europe Sp. z o.o.**

Twarda 18; PL-00105 Warsaw

T: +48 22 202 67 23; F: +48 22 202 67 19

[www.sol-m.com](http://www.sol-m.com)