

## 1, 2 & 3-Gallon TakeAway Recovery Systems

1-Gallon #81000, 2-Gallon #82000, 3-Gallon #83000

**Small quantity generators in need of a cost-effective, DOT-Approved, and regulatory-compliant disposal solution for sharps will find the TakeAway Recovery System the perfect fit for their needs.**

The TakeAway Recovery System is an innovative system designed for simplified and environmentally friendly treatment of used needles, syringes, and other used healthcare products.

### SYSTEM INCLUDES

- Sharps collection container
- Prepaid UPS return shipping box
- UPS tracking document
- Proper treatment and online documentation for proof of recovery

### ADVANTAGES

- Cost-effective alternative to pickup service, reducing costs by as much as 50% or more
- No contracts or hidden fees - order systems as needed
- Approved for sharps and other used healthcare products, keeping you in compliance
- Easy to package for return, saving you time for other important job responsibilities
- Online tracking and proof of recovery via SharpsTracer, available 24/7, giving you complete cradle-to-grave protection



1-Gallon #81000

3-Gallon #83000

## 1-Gallon TakeAway Recovery System

### Product Specifications and Shipping Data

#### #81000 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	8.75" L	5.5" W	7" H
Single Mailback Box:	9.75" L	6.75" W	7.75" H
Shipping Box:	10.75" L	7" W	8" H
UPC Code:	634188810009		
Max Return Weight:	10 lbs		
Units / Pallet:	132		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	3168 cases		

## 2-Gallon TakeAway Recovery System

### Product Specifications and Shipping Data

#### #82000 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	8.75" L	5.5" W	10.75" H
Single Mailback Box:	9.75" L	6.75" W	11.75" H
Shipping Box:	10.75" L	7" W	12" H
UPC Code:	634188820008		
Max Return Weight:	17 lbs		
Units / Pallet:	144		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	3456 cases		

## 3-Gallon TakeAway Recovery System

### Product Specifications and Shipping Data

#### #83000 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	8.75" L	5.5" W	16.75" H
Single Mailback Box:	9.75" L	6.75" W	17.75" H
Shipping Box:	10.75" L	7" W	18" H
UPC Code:	634188830007		
Max Return Weight:	25 lbs		
Units / Pallet:	80		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	1920 cases		

## Wall Mount Brackets for 1-, 2-, and 3-Gallon Sharps Collection Containers

#50004, #50005, #50014

Conveniently mount 1-, 2-, and 3-gallon sharps containers with three bracket styles: basket, front, and side.



For more information, please call **800.772.5657** or visit us online at **www.sharpsinc.com**.