

20-Gallon Sharps Recovery System for Trace Chemo

20-Gallon #40020

The 20-Gallon Sharps Recovery System for Trace Chemo/ Cytotoxic Drugs is the ideal solution for disposing of trace chemo waste.

The Sharps Recovery System is designed for regulatory-compliant containment, transport and treatment of needles, empty vials and syringes and other items contaminated with trace amounts of chemotherapy or other cytotoxic drugs or chemicals. The leak-proof, puncture-resistant container accommodates your smaller existing sharps containers (no loose sharps allowed).

SYSTEM INCLUDES

- Collection container
- Prepaid USPS return shipping box
- Tracking form
- Proper treatment and online documentation of treated waste

ADVANTAGES

- Cost-effective alternative to pickup service, reducing costs by as much as 50% or more
- No contracts or hidden fees — order systems as needed
- Approved for sharps, other medical waste and empty vials/ syringes of trace amounts of chemotherapy or cytotoxic drug or chemical waste*
- Easy to package for return, saving you time for other important job responsibilities
- Option to schedule free USPS priority pickup online at your convenience
- Online tracking and proof of treatment via SharpsTracer available 24/7, giving you complete cradle-to-grave protection

*Contact Sharps Compliance at 800.772.5657 for proper disposal of sharps and other items contaminated with acutely hazardous (p-listed) or bulk hazardous waste.



20-Gallon #40020

20-GALLON SYSTEM ACCOMMODATES THE FOLLOWING SHARPS CONTAINERS

1-Quart:	1.5-Quart:	5-Quart:	1-Gallon:	2-Gallon:	3-Gallon:
up to 28	up to 18	up to 6	up to 8	up to 5	up to 3

20-Gallon Sharps Recovery System

Product Specifications and Shipping Data

#40020 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	21.5" L	21.25" W	18.5" H
Single Mailback Box:	22" L	22" W	20" H
Shipping Box:	22" L	22" W	20" H
UPC Code:	634188000660		
Max Return Weight:	35 lbs		
Units / Pallet:	12 units		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	288 cases		

Rolling Base

#80005

Easily maneuver larger systems throughout your facility with the steel Rolling Base, perfect for the 20-Gallon System. The rolling base features four heavy-duty casters and locking wheels.



For more information, please call **800.772.5657** or visit us online at www.sharpsinc.com.