


 | *Grip-Lok*[®]

Safe. Simple.
Secure.

DESIGNED TO LOCK YOUR LINE.

Grip-Lok[®] securement solutions is a portfolio of versatile adhesive-based devices that promotes patient comfort and secures a variety of tubes, lines, and catheters.





Grip-Lok's "Line Lock design" allows clinicians to open and close the device to reposition the line without removing the securement device from the skin.

+ WRAP-AROUND ADHESIVE DESIGN

+ PHYSICAL BARRIER BETWEEN TUBE AND SKIN

DESIGNED TO LOCK YOUR LINE

Unlike tape, which only adheres to the top of the tube, Grip-Lok's wrap-around design provides security around the entire tube, while also providing a physical barrier between the tube and the skin.



TAPE

LINE SECUREMENT

Using Grip-Lok to secure tubes and lines allows clinicians to standardize line securement for consistency across all units.



TAPE



GRIP-LOK®

SKU #	DESCRIPTION	QTY	SIZE RANGE	CATHETER / DEVICE MATERIALS
3200S	SMALL SECUREMENT DEVICE	100 / BX	4.5 - 13.5 FR	PVC, POLYETHYLENE, POLYURETHANE
3300M	MEDIUM SECUREMENT DEVICE	100 / BX	9 - 24 FR	PVC, POLYETHYLENE, POLYURETHANE
3400L	LARGE SECUREMENT DEVICE	100 / BX	16 - 40 FR	PVC, POLYETHYLENE, POLYURETHANE
3300MWA	MEDIUM WIDE SECUREMENT DEVICE	100 / BX	16 - 40 FR	SILICONE, PVC, POLYETHYLENE, POLYURETHANE

