



IV Clear[®]

The Essentials
of Infection
Prevention,
All in One.

IV Clear[®] is the only antimicrobial vascular access dressing that offers:

- ✓ Complete surface, dual antimicrobial action and protection
- ✓ Full site visibility
- ✓ Silicone adhesive to preserve skin

covalon
TECHNOLOGIES LTD.

Broad Spectrum Kill

Antimicrobials are essential in combatting infection in vascular access, but not all antimicrobials are considered equal.

- ✓ IV Clear® incorporates two antimicrobials across the entire surface of the dressing, offering **full surface coverage** to protect against bacterial migration by preventing regrowth for up to 7 days
- ✓ IV Clear® uses chlorhexidine and silver to provide broad spectrum protection to kill on average **99.99% of microbes** associated with CRBSIs
- ✓ IV Clear® uses **safe and effective** concentrations to protect against irritant contact dermatitis

Total Visibility

In infection prevention, early detection is crucial for ensuring positive patient outcomes.



- ✓ IV Clear® offers complete transparency at and around insertion site for easy daily assessment

Minimizing Skin Injuries

Damaging the skin not only causes pain, but it can delay healing and increase the risk of infection.

- ✓ IV Clear® uses a silicone adhesive that preserves the skin barrier function and minimizes skin injuries


Availability

Material Number/ Part Number	Size (cm)
 TWBD1020 TWBD1014	6cm x 7cm 10cm x 12cm
 TWBD1016 TWBD1012	6cm x 7cm 10cm x 12cm




References

- Maene, B. Hidden Costs of medical tape-induced skin injuries. (2013). Wounds UK, Vol 9(1): 46-49.
- 3M. (2016). Tegaderm CHG Dressing Data Sheet. Retrieved from <https://multimedia.3m.com/mws/media/4883820/tegaderm-chg-dressing-datasheet-1657-1658-1659.pdf>
- Johnson & Johnson. (2018). BIOPATCH Disk Core Brochure. Retrieved from https://www.jnjmedicaldevices.com/sites/default/files/user_uploaded_assets/pdf_assets/2018-10/032669-

Area of Antimicrobial Coverage



IV Clear® provides the greatest area of antimicrobial coverage in order to minimize bacterial regrowth and migration.*

Brands	Area of Antimicrobial Coverage ^{1,2} (most comparable size)
 IV Clear®	120cm ²
 Competitor 1 CHG gel pad	12cm ²
 Competitor 2 CHG disk	5cm ²

www.covalon.com
customerservice@covalon.com
1.888.236.4944

covalon
TECHNOLOGIES LTD.