




Trifurcated Pharmacy Transfer Tubing with Three Large Bore Spike Leads

Our custom Trifurcated Lead Pharmacy Fluid Transfer Set features three large volume intake tubing with a universal spike, silicone tubing insert, and a male luer lock at the distal end.

- For use with Baxa Repeater[®] Pump
- Three large volume intake tubing
- Universal vented spike
- Silicone tubing insert
- DEHP-Free
- Latex-Free
- Pyrogen-Free
- Male luer lock at distal end
- Sterile



Product # TCBPTT063

Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.