



Sterile Isopropyl Alcohol Wipers

Our sterile isopropyl alcohol wipers are 70% saturated and are low linting. Composed of a 55/45 cellulose blend, each 9" x 9" wiper is 65gsm, with 30 wipers per zip pack, and 24 packs per case.

- Gamma irradiated to a 10⁻⁶ sterility assurance level
- Saturated with isopropyl alcohol, 70% USP
- Alcohol is filtered twice through 0.2µm filters and then terminally sterilized
- 9" x 9" size
- Extremely low linting, 65gsm cellulose blend
- Double bagged for added cleanroom protection



Product # TCBWALC30

Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.



Sterile Cleanroom Wipers

Our sterile cleanroom wipers are sterile gamma irradiated to a 10⁻⁶ sterility assurance level. They are constructed of an extremely low linting 68gsm 55/45 cellulose blend. They are appropriate for critical areas up to and including ISO 5 cleanroom environments. The Cleanroom Wipers feature maximum absorbencies of 20-25 mL, and are available in 9" x 9" and 12" x 12" sizes, and available in different quantities.

Product # TCBWIP09

300 wipers per pack, 12 packs per case, a total of 3600 wipers per case

Product # TCBWIP09SP-20

20 wipers per bag, 15 bags per bundle, 12 bundles per case, a total of 3600 wipers per case

Product # TCBWIP12

150 wipers per pack, 8 packs per case, a total of 1200 wipers per case

Product # TCBWIP12SP-20

20 wipers per bag, 8 bags per bundle, 8 bundles per case, a total of 1280 wipers per case



Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBWIP09
Product # TCBWIP09SP-20
Product # TCBWIP12
Product # TCBWIP12SP-20



Sterile Lab Gowns for Compounding

Our sterile SMS 40gsm weight white anti-static treated cleanroom gowns are designed for non-hazardous drug compounding. The gowns are low linting and latex-free. They feature a stand-up collar, elastic cuffs and thumb loops. They come individually packaged 1 pc/sterile bag with 50 gowns to a case.

- Low linting
- Thumb loops
- Stand up collar
- ISO-certified breathable fabric
- Latex-Free
- Available in sizes Small, Regular, XL and XXXL



Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBA40ST-S
Product # TCBA40ST-R
Product # TCBA40ST-XL
Product # TCBA40ST-XXXL



Sterile Cleanroom Coverall

Our Sterile Cleanroom Coverall, composed of 60gsm weight anti-static material, is designed to meet USP 800 requirements. The Coverall is breathable with PE/PP and offers Level 3 protection against chemo drugs. They come individually packaged, 1pc/sterile bag, 20 bags/case.



- Polypropylene spunbond fabric with a polyethylene exterior coating
- ASTM tested
- ISO-certified breathable fabric
- Level 3 protection
- Low linting
- Thumb loops
- Seam seals
- Zippered with adhesive front cover
- Elastic cuffs
- Latex-Free
- Available in sizes Small, Regular, XL and XXXL

Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBACV54ST-S
Product # TCBACV54ST-R
Product # TCBACV54ST-XL
Product # TCBACV54ST-XXXL



Sterile Lab Gowns for Hazardous Drug Compounding

Our Sterile Cleanroom Gown is disposable and composed of 52gsm weight anti-static material designed to meet USP 800 requirements. The Gown is breathable with PE/PP and offers Level 3 protection against chemo drugs. They come individually packaged, 1pc/sterile bag, 50 bags/case.

- Polypropylene spunbond fabric with a polyethylene exterior coating
- ASTM tested
- ISO-certified breathable fabric
- Level 3 protection
- Low linting
- Thumb loops
- Seam seals
- Stand up collar
- Elastic cuffs
- Latex-Free
- Available in sizes Small, Regular, XL and XXXL



Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBA54ST-S
Product # TCBA54ST-R
Product # TCBA54ST-XL
Product # TCBA54ST-XXXL



Micro Bore IV Extension Set

Our Micro Bore IV Extension Sets are equipped with male and female luer locks and a slide clamp. The sets are DEHP-Free and Latex-Free, and are packaged in quantities of 50 per box.

Product # TCBEXT001

Male and female luer locks, slide clamp, 7" length, DEHP-Free, Latex-Free, 0.21mL priming volume

Product # TCBEXT003

Male and female luer locks, slide clamp, 15" length, DEHP-Free, Latex-Free, 0.35mL priming volume

Product # TCBEXT004

Male and female luer locks, slide clamp, 30" length, DEHP-Free, Latex-Free, 0.6mL priming volume



Product # TCBEXT001
Product # TCBEXT003
Product # TCBEXT004

Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.



Ambulatory Pump Pouch

Our ambulatory pump pouches offer extreme durability and quality at low costs, featuring high quality, 100% polyester (600 X 3000) fabric. The pouch interiors feature a polyurethane backing and reinforced seams for added strength and resilience. Our pouches offer easy-to-use Velcro™ style strap systems to secure pump and cassette/bag. The unique pump pocket design is universal in nature, enabling the pouch/backpack to accommodate all brands of ambulatory pumps.

Product # TCBP100

Universal Pump Pouch, shoulder swivel, convertible to fanny waist belt, loop to secure on belt, window, pump pocket and inner wrap accommodates 50-100mL cassette or medication/solution bag

Product # TCBP510

Universal Pump Pouch, Elite, two adjustable straps included (fanny waist and shoulder swivel), window/privacy flap, ancillary zip, pump pocket and inner wrap accommodates 250-500mL medication/solution bag

Product # TCBP1100

Universal Ambulatory Pump Pouch, Elite, adjustable shoulder swivel strap (fanny waist strap additional), window/privacy flap, ancillary zip, pump pocket and inner wrap accommodates 750-1250mL medication/solution bag

Product # TCBP3200

Universal Ambulatory Pump Backpack, Elite, padded adjustable straps, enhanced solution fastening system, window/privacy flap, ancillary zip, pump pocket and inner wrap accommodates 1500-3500mL medication/solution bag

Product # STRAP-530

Additional fanny pack belt



Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBP100
Product # TCBP510
Product # TCBP1100
Product # TCBP3200
Product # STRAP-530



Pharmacy Transfer Tubing

Our custom Single Lead Pharmacy Fluid Transfer Set features a universal spike, silicone tubing insert, and a male luer lock at the distal end.

- Compatible with Baxa Repeater® pump
- Single lead
- Universal vented spike
- Silicone tubing insert
- DEHP-Free
- Latex-Free
- Pyrogen-Free
- Male swivel luer lock at distal end



Product # TCBPTT001

Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.



PharmAssist[®] Compatible Solution Transfer Tubing

Our custom PharmAssist[®] compatible transfer tubing pump sets are designed for use with the Excelsior PharmAssist[®] pump and are constructed for high volume use.

- Universal vented spike
- Silicone tubing
- DEHP-Free
- Latex-Free
- One pinch clamp
- Male luer lock at distal end
- High volume capacity
- Pyrogen-Free



Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBPTPS1



Pharmacy Transfer Tubing (Weighted)

Our custom Single Lead Pharmacy Fluid Transfer Set features large volume intake tubing, a weighted male luer lock intake connector, and silicone tubing insert.

- For use with Baxa Repeater[®] Pump
- Large volume intake tubing
- Weighted male luer lock intake connector
- Silicone tubing Insert
- DEHP-Free
- Latex-Free
- Pyrogen-Free
- Sterile



Product # TCBPTT013

Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.



Pharmacy Transfer Tubing with Large Bore MLL

Our custom Single Lead Pharmacy Fluid Transfer Set features a universal male luer lock and silicone tubing insert.

- For use with Baxa Repeater[®] Pump
- Large bore intake tubing
- Universal male luer lock
- Silicone tubing insert
- DEHP-Free
- Latex-Free
- Pyrogen-Free
- Sterile



Product # TCBPTT021


Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.



Bifurcated Pharmacy Transfer Tubing with Two Large Bore Spike Leads

Our custom Bifurcated Lead Pharmacy Fluid Transfer Set features two large volume intake tubing with a universal spike, silicone tubing insert, and a male luer lock at the distal end.

- For use with Baxa Repeater[®] Pump
- Two large volume intake tubing
- Universal vented spike
- Silicone tubing insert
- DEHP-Free
- Latex-Free
- Pyrogen-Free
- Male luer lock at distal end
- Sterile



Product # TCBPTT062


Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.



Trifurcated Pharmacy Transfer Tubing with Three Large Bore Spike Leads

Our custom Trifurcated Lead Pharmacy Fluid Transfer Set features three large volume intake tubing with a universal spike, silicone tubing insert, and a male luer lock at the distal end.

- For use with Baxa Repeater[®] Pump
- Three large volume intake tubing
- Universal vented spike
- Silicone tubing insert
- DEHP-Free
- Latex-Free
- Pyrogen-Free
- Male luer lock at distal end
- Sterile



Product # TCBPTT063

Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.



Baxter FloGard® Compatible IV Administration Sets

Our Baxter FloGard® Compatible IV Administration sets are individually blister-packaged in quantities of 50 per box.



TCBINF621: Universal vented/non-vented 10 drop/mL spike, DEHP-Free, PVC fluid path, 100" length, one distal Y-Site, roller clamp, slide clamp, Pyrogen-Free, Latex-Free



TCBINF6419: Universal vented/non-vented 10 drop/mL spike, DEHP-Free, PVC fluid path, 102" length, one Clave® Y-Site, roller clamp, slide clamp, Pyrogen-Free, Latex-Free



TCBINF6480: Universal vented/non-vented 10 drop/mL spike, DEHP-Free, PVC fluid path, 102" length, one Clave® Y-Site, GVS Speedflow® 0.2µm Filter, roller clamp, slide clamp, Pyrogen-Free, Latex-Free



TCBINF6481: Universal vented/non-vented 10 drop/mL spike, DEHP-Free, PVC fluid path, 102" length, one Clave® Y-Site, GVS Speedflow® 1.2µm Filter, roller clamp, slide clamp, Pyrogen-Free, Latex-Free



TCBINF6519: Universal vented/non-vented 10 drop/mL spike, DEHP-Free, PVC fluid path, 97" length, two Clave® Y-Sites, backcheck valve, roller clamp, slide clamp, Pyrogen-Free, Latex-Free



TCBINF6537: Universal vented/non-vented 10 drop/mL spike, DEHP-Free, PVC fluid path, 110" length, three Clave® Y-Sites, backcheck valve, roller clamp, slide clamp, Pyrogen-Free, Latex-Free



TCBINF6571: Universal vented/non-vented 10 drop/mL spike, DEHP-Free, PVC fluid path, 105" length, two Clave® Y-Sites, backcheck valve, GVS Speedflow® 0.2µm Filter, roller clamp, slide clamp, Pyrogen-Free, Latex-Free

Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBINF621
Product # TCBINF6419
Product # TCBINF6480
Product # TCBINF6481
Product # TCBINF6519
Product # TCBINF6537
Product # TCBINF6571