



TCBWALC6004

Sterile Non-Woven Polyester Presaturated Wipes

- + Our wipes ensure consistent, repeatable application for enhanced contamination control.
- + Convenience: pre-saturated with isopropyl alcohol 70% USP and 30% deionized water.
- + 100% non-woven polyester: 60 g/m²
- + Gamma irradiated to a 10⁻⁶ sterility assurance level.
- + Secure packaging: Peel-and-reseal design prevents evaporation when stored correctly.
- + Double bagged for added cleanroom protection.



Product #	Size	Quantity
TCBWALC6004	4" x 4"	60 wipes/pouch; 4 pouches/inner bag; 2 inner bags/outer bag; 3 outer bags/case; 24 pouches/case, 1440 wipes/case

