



2022 Enteral Feeding Devices Product Catalog



Vesco

(*vēscō*)

Latin to English Translation
to feed, to eat



The Droplet

Our “Nike Swoosh” symbolic of a drop of nutrition artistically flowing into carefully coiled tubing



Vision

Worldwide leader delivering the broadest portfolio of the highest value enteral feeding devices



Feeding Innovation

Our promise to deliver breakthrough enteral products to better serve our customers





VESCO

MEDICAL™



High-Quality Devices

- ISO 13485 Certified Quality Management System.
- CE Marked, FDA Registered & 510(k) Premarket Clearance.
- Made without natural rubber latex, DEHP, BPA or other phthalates plasticizers.



Unbeatable Service

- Friendly, skilled and experienced customer care specialists.
- Full portfolio of ENFit products direct or through distributors.
- Responsive when and where you need us.



Competitive Pricing

- Nimble, agile and driven to keep costs down and pass on savings.
- High-quality products without the premium prices.

ENFit Tip Syringes

ENFit Tip Syringes


- Single Use syringes compliant with ISO 80369-3 Standard
- For Feeding, Flushing & Liquid Medication Administration
- Sterile and Clean/Non-Sterile

ENFit Standard Tip Syringes (5, 10, 20, 35 & 60mL)


- Single use enteral syringes for use in initial tube placement, flushing/hydration, bolus feeding administering medication and irrigation.
- Printed graduated length markings in mL.
- One cap included in each blister pack.
- Made without natural rubber latex, DEHP, Bisphenol A (BPA) or other phthalates plasticizers.
- Syringes are sold in clean/non-sterile and sterile configurations. Sterile devices are sterilized via (EO for ethyl-o-xylene). Do not re-sterilize.

ENFit Low Dose Tip (LDT) Syringes (0.5, 1 & 3mL)

- Designed to minimize process changes in pharmacy and bedside administration.
- Demonstrated dose accuracy administering smaller volumes compliant with GEDSA*.
- Internal lumen consistent with traditional male syringe design reduces dead space and priming volume.
- Made without natural rubber latex, DEHP, Bisphenol A (BPA) or other phthalates plasticizers.
- Device is sterilized (EO for ethyl-o-xylene). Do not re-sterilize.




ENFit
Standard Tip



ENFit
Low Dose Tip (LDT)

ENFit Syringe Item Numbers



ENFit Tip Syringes

Size (mL)	Single Use - Sterile		Clean/Non-Sterile	
	Item #	Units/Case	Item #	Units/Case
1 (LDT)	VED-601EO	800	N/A	N/A
3 (LDT)	VED-603EO	800	N/A	N/A
5	VED-605EO	800	VED-605	500
10	VED-610EO	800	VED-610	500
20	VED-620EO	600	VED-620	500
35	VED-635EO	400	VED-635	500
60	VED-660EO	400	VED-660	180
60 Thumb Ring	N/A	N/A	VED-661	180
60 Thumb Ring + 500mL Container	N/A	N/A	VED-665	180

ENFit Syringe Pharmacy Bulk Pack

Size (mL)	Item #	Units/ Case	Amber Item #	Units/ Case
0.5mL (LDT)	VED-600B	4000	VED-600AB	4000
1 (LDT)	VED-601B	4000	VED-601AB	4000
3 (LDT)	VED-603B	4000	VED-603AB	4000
5	VED-605B	2000	VED-605AB	1000
10	VED-610B	1000	VED-610AB	1000
20	VED-620B	1000	VED-620AB	1000
35	VED-635B	800	VED-635AB	800
60	VED-660B	400	VED-660AB	400

Nasogastric Feeding Tubes

Nasogastric Polyurethane (PU) Feeding Tubes

- Single or Dual ENFit Ports
- French Sizes: 5, 6, 8, 10, 12, 14 & 16
- Optional Stylet and/or Weighted Tip
- Lengths: 18, 36, 45 & 55 Inches

Features

- Single or Dual ENFit Port for feeding and administering medication.
- Premium medical-grade kink resistant polyurethane tubing.
- 4 side ports and distal end port improves flow and reduces likelihood of tube clogging.
- Water-activated lubricant aids in tube placement and allows for aseptic lubrication.
- Radiopaque for x-ray confirmation of tube placement.
- Centimeter graduated length markings.
- Available with and without flow through stylet for ease of placement and confirmation.
- Available non-weighted or non-magnetic tungsten weighted (closed tip).
- Not made with natural rubber latex, DEHP, Bisphenol A (BPA) or other phthalate plasticizers.



Single & Dual ENFit Port NG Tubes



Important Safety Information

- Use only administration/feeding sets with ENFit connectors. Feeding sets with different connectors are not compatible with Vesco feeding tubes.
- Device is non sterile.
- Tube has the potential to misconnect with small bore connectors of other healthcare applications.
- Single patient use and Rx only.

NG Tube Item Numbers

Single & Dual Port Polyurethane (PU) NG Tubes					
French Size	18" (45cm) Item #	36" (91cm) Item #		45" (114cm) Item #	Units/ Case
	Single Port	Single Port	Dual Port	Dual Port	
5	VED-14505	VED-19105	N/A	N/A	10
6	VED-14506	VED-19106	VED-29106	VED-21406	10
8	VED-14508	VED-19108	VED-29108	VED-21408	10
10	N/A	VED-19110	VED-29110	VED-21410	10
12	N/A	VED-19112	VED-29112	VED-21412	10
14	N/A	VED-19114	VED-29114	VED-21414	10
16	N/A	N/A	VED-29116	N/A	10

Single & Dual Port Polyurethane (PU) NG Tubes with STYLET				
French Size	36" (91cm) Item #		45" (114cm) Item #	Units/ Case
	Single Port	Dual Port	Dual Port	
6	N/A	VED-29106S	VED-21406S	10
8	VED-19108S	VED-29108S	VED-21408S	10
10	VED-19110S	VED-29110S	VED-21410S	10
12	VED-19112S	VED-29112S	VED-21412S	10
14	VED-19114S	VED-29114S	VED-21414S	10
16	N/A	VED-29116S	N/A	10

Dual Port Polyurethane (PU) NG Tubes, WEIGHTED with STYLET					
French Size	Weight (Grams)	36" (91cm) Item #	45" (114cm) Item #	55" (140cm) Item #	Units/ Case
6	0.55	VED-29106SW	VED-21406SW	N/A	10
8	0.75	VED-29108SW	VED-21408SW	VED-24008SW	10
10	1.33	VED-29110SW	VED-21410SW	VED-24010SW	10
12	2.07	VED-29112SW	VED-21412SW	VED-24012SW	10

Pediatric Nasogastric Tubes

Pediatric Polyurethane (PU) NG Tubes with ENFit

- For use up to 30 days
- French Sizes: 4, 5, 6, 8, 10, 12
- Lengths: 16, 24 and 35 Inches (40, 60 & 90cm)

Features

- Clear and transparent tube for monitoring fluid content
- Single ENFit port for feeding, flushing and administering medication.
- Indicated for long term use up to 30 days.
- Premium medical-grade kink resistant polyurethane tubing.
- 2 side ports improves flow and reduces likelihood of tube clogging.
- Radiopaque for x-ray confirmation of tube placement.
- Device is sterilized (EO for ethyl-o-xylene). Do not re-sterilize.
- Not made with natural rubber latex, DEHP, Bisphenol A (BPA) or other phthalate plasticizers.
- Centimeter graduated length markings.



Pediatric
NG Tubes

Important Safety Information

- Use only administration/feeding sets with ENFit connectors. Feeding sets with different connectors are not compatible with Vesco feeding tubes.
- Devices are sterilized (EO for ethyl-o-xylene). Do not re-sterilize.
- Tube has the potential to misconnect with small bore connectors of other healthcare applications.
- Single patient use.

Pediatric NG Tube Item Numbers



New

Pediatric Polyurethane NG Tubes - Sterile

French Size	16" (40cm) Item #	24" (60cm) Item #	35" (90cm) Item #	Units/Case
4	VED-84004EO	VED-86004EO	N/A	40
5	VED-84005EO	VED-86005EO	N/A	40
6	N/A	VED-86006EO	VED-89006EO	40
8	N/A	VED-86008EO	VED-89008EO	40
10	N/A	VED-86010EO	VED-89010EO	40
12	N/A	N/A	VED-89012EO	40



ENFit Pharmacy Items

ENFit Medicine Bottle Adapters (Multiple Sizes)

- Variety of sizes intended to meet most common medication bottles.
- Fits inside bottle opening to allow continued use of child proof caps.
- Ideal for aseptic technique.
- Made without natural rubber latex, DEHP, Bisphenol A (BPA) or other phthalates plasticizers.

ENFit Transition Connectors

- Stepped Transition Connector for backwards compatibility to legacy feeding tubes.
- Available as a stand-alone item by the case.

ENFit Tamper Evident Syringe Caps (Fits All Sizes)


- Provides security to syringe contents during transportation, prior to administration and enhances medication safety protocols.
- Made without natural rubber latex, DEHP, Bisphenol A (BPA) or other phthalates plasticizers.

ENFit Syringe Tip Cap (Fits All Sizes)


- Pharmacy friendly purple cap designed to aid in single handed placement.
- Ideal for aseptic technique and secure for transport.
- Made without natural rubber latex, DEHP, Bisphenol A (BPA) or other phthalates plasticizers.

ENFit Medicine Straw (50, 100 & 150mm)


- Connects to an ENFit Tip Syringe to easily facilitate the withdrawal of medication or feed.
- Ideal for aseptic technique.
- Device is sterilized (EO for ethyl-o-xylene). Do not re-sterilize.



Transition Connectors



ENFit Medicine
Bottle Adapters



ENFit Medicine Straw

Important Safety Information

- Use ENFit tip syringes only with ENFit feeding tube ports. Feeding tubes with different open female ports are not compatible with Vesco ENFit tip syringes.
- Single-patient use.

Pharmacy Item Numbers



ENFit Medicine Bottle Adapters			
Size #	Inner Diameter (mm)	Item #	Units/Case
1	12 - 15	VED-690	100
2	14.5 - 16.5	VED-691	100
3	17 - 20.5	VED-692	100
4	19 - 23	VED-693	100
5	24 - 26	VED-694	100
6	27-30	VED-695	100
7	31-34	VED-696	100

New

Ancillary Items	Item #	Units/Case
Transition Connectors	VED-001	300
2 Inch (5cm) ENFit Medicine Straw	VED-670	100
4 Inch (10cm) ENFit Medicine Straw - Sterile	VED-671EO	150
6 Inch (15cm) ENFit Medicine Straw - Sterile	VED-672EO	150
ENFit Syringe Screw-on-Cap	VED-697	4000
ENFit Syringe Pop-on-Cap (Dark)	VED-699P	4000
ENFit Syringe Pop-on-Cap (Lite)	VED-699PF	4000
ENFit Tamper Evident Cap	57-220	1000

New




Extension Feeding Sets

Universal Low Profile Feeding Tube Extension Sets

- Bolus Feed, Single ENFit Port with Straight or Right Angle available in 12" (30cm) and 24" (60cm)
- Continuous Feed, Dual ENFit Port available in 12" (30cm) and 24" (60cm)

Features

- Compatible with MIC Key® and MINI One® and all other low profile gastrostomy feeding tubes.
- Dual ENFit Ports with 90 degree connector allows for continuous feeding, easy flushing and administration of liquid medication without disconnecting set.
- Single ENFit port & straight connector for bolus feeding, flushing and medication delivery.
- Not made with natural rubber latex, Bisphenol A (BPA) or other phthalate plasticizers.



Low Profile
Extension Sets

Important Safety Information

- Use only administration/feeding sets with ENFit connectors. Feeding sets with different connectors are not compatible with Vesco feeding tubes.
- See individual product labeling for directions and important safety information.
- Single patient use.

Extension Set Item Numbers



Universal Low Profile Extension Sets			
Set Description	12" (30cm) Item #	24" (60cm) Item #	Units/Case
Continuous - Dual Port, Right Angle	VED-1000	VED-1001	10
Bolus - Single Port, Straight	VED-1002	VED-1003	10
Bolus - Single Port, Right Angle	VED-1004	VED-1005	10




Gravity Feeding Sets

Gravity Feeding Sets with ENFit

- 500 & 1000mL Gravity Bag Sets
- ENPlus Spike and 40mm Screw Cap Closed System Gravity Sets
- Large Bore 1000mL Gravity Bag Set

Features

- Gravity Bag Sets
 - Snap cap with rigid neck for easy filling and closure.
 - Easy-to-read graduations and easy view translucent material.
- Closed System Sets
 - Screw cap allows for flexibility in feeding containers with a 40mm opening.
 - Convenient suspension bag allows for use with smaller containers (baby bottles) for easy hanging.
- Roller clamp and 14 drip chamber for easy priming and consistent flow.
- Large Bore Gravity Sets designed to feed thicker formulas faster with 30° tapered bag for quicker, more complete feed.
- Not made with natural rubber latex, DEHP, Bisphenol A (BPA) or other phthalate plasticizers.
- Device intended to deliver liquid nutrition formulas or water to a patient's feeding tube.



1000mL
Gravity Feed Set

Important Safety Information

- To avoid risk of a tubing misconnection, ONLY connect the distal connector to an ENTERAL feeding tube. This device is compatible with a flexible female connector on the patient's enteral feeding tube.
- Not for IV use.
- This device should be replaced at least every 24 hours.
- Do not store for prolonged periods at extreme temperatures.
- See individual product labeling for directions and important safety information.
- Single patient use.

Gravity Set Item Numbers



Gravity Feed Sets	Transition Set (w/TC) Item #	NEW Item# (No TC)*	Units/Case
40mm Screw Cap Set with ENFit	VED-020	VED-022	30
ENPlus Spike with ENFit	VED-026	VED-028	30
500mL Bag Set (LL) with ENFit	VED-038	VED-040	30
1000mL Bag Set (LL) with ENFit	VED-032	VED-034	30
1000mL TOP FILL Bag Set (PVC) with ENFit	VED-045	VED-046	30
1000mL LARGE BORE Bag Set with ENFit New	VED-048	VED-049	30

* New Item # refers to the same gravity set as previous without the Transition Connector.



Gastrostomy Feeding Tubes

Dual ENFit Port Replacement Gastrostomy Feeding Tubes

- French Sizes: 12, 14, 16, 18, 20, 22 and 24

Features

- Dual ENFit Port allows for easy flushing and administration of liquid medication without disconnecting set.
- Medical-grade silicone tubing with clear shaft and graduation marks in 1cm increments.
- Rounded skin disk for retention.
- Gastric balloon for internal retention with rounded tip, open distal end and radiopaque tip marker.
- Device is sterilized (EO for ethyl-o-xylene). Do not re-sterilize.
- Not made with natural rubber latex, Bisphenol A (BPA) or other phthalate plasticizers.



Single ENFit Port Replacement Gastrostomy Feeding Tubes

- French Sizes: 16, 18, 20 and 22

Features

- Biocompatible-grade silicone tubing with clear shaft and graduation marks in 2cm increments.
- Rounded skin disc with bumpers for retention.
- Gastric balloon for internal retention with rounded flush tip and open distal end.
- Device is sterilized (EO for ethyl-o-xylene). Do not re-sterilize.
- Not made with natural rubber latex, Bisphenol A or other phthalate plasticizers.



Important Safety Information

- Use only administration/feeding sets with ENFit connectors.
- Single patient use.

G-Tube Item Numbers



Dual ENFit Port Gastrostomy Tubes - Sterile			
French Size	Balloon Size (mL)	Item #	Units/Case
14	5	VED-214	5
16	5	VED-216	5
18	20	VED-218	5
20	20	VED-220	5
22	20	VED-222	5
24	20	VED-224	5

Single ENFit Port Gastrostomy Tubes - Sterile			
French Size	Balloon Size (mL)	Item #	Units/Case
16	10	VED-116	2
18	15	VED-118	2
20	15	VED-120	2
22	15	VED-122	2




Drainage Bags

Enteral Drainage Bags with ENFit

- 250, 500, 1000 and 1500mL Drainage Bag Sets
- 1000mL and 1500mL Drainage Bag with 75cm proximal tubing

Features

- Designed specifically for drainage in enteral feeding systems with male and female ENFit connectors compliant with ISO 80369-3 Small Bore Connector Standards.
- Large, clear tubing allows for visualization of gastric contents and facilitates gastric decompression or drainage.
- Eliminates the need for non-enteral drainage products such as urinary bags.
- Tethered cap & thumb clamp.
- Larger sizes (1000mL & 1500mL) equipped with 30 inches of extension tubing.
- Smaller sizes (250mL & 500mL) have 2 inches of proximal tubing.
- Not made with natural rubber latex, DEHP or other phthalate plasticizers.



Important Safety Information

- Use only with enteral feeding tubes with ENFit connectors.
- Device is sterilized (EO for Ethyl-o-xylene). Do not re-sterilize.
- Single patient use.

Drainage Bag Item Numbers



New

Drainage Bag Sets - Sterile	Item #	Units/Case
250mL PVC ENFit Drainage Bag with 2" (5cm) tubing	VED-009EO	30
500mL PVC ENFit Drainage Bag with 2" (5cm) tubing	VED-010EO	30
1000mL PVC ENFit Drainage Bag with 30" (75cm) tubing	VED-012EO	30
1500mL PVC ENFit Drainage Bag with 30" (75cm) tubing	VED-015EO	30



Vesco Medical

Tube feeding experts delivering high quality products with unbeatable service at a competitive price. Our portfolio is All ENFit All The Time!

Group Purchasing Organization Agreements:

HealthTrust (HT31810) MHA (MS60238N)

Premier (PP-DI-1468) Vizient (MS6216)

www.VescoMedical.com

info@vescomedical.com

(855) 239-0131

GEDSA

ENFit is a federally registered trademark of GEDSA used with permission.
©2022 Vesco Medical, LLC MKT 1128-8-02



For more information regarding new ISO Standard ENFit connectors, visit www.stayconnected.org



1039 Kingsmill Parkway
Columbus, OH 43229 USA