

AirLife[®] heated circuit and chamber for adults

Fact sheet and performance specifications



AirLife[®]



AirLife high-performance, dual-limb heated circuit for adults (Product codes: AH102 and AH180)

Duration of use:	Up to 30 days
Circuit length:	5'
Heater platform compatibility:	Compatible with Fisher & Paykel MR850 heater base
Weight:	AH102: 141 g AH180: 205 g
DEHP/BPA/latex information:	AH102/AH180: Manufactured without DEHP, BPA or latex
Resistance to flow @ 30 LPM:	AH102: Inspiratory limb: 0.02 hPa/L/min (<i>cm H₂O/L/min</i>) AH180: Inspiratory limb: 0.02 hPa/L/min (<i>cm H₂O/L/min</i>) Expiratory limb: 0.02 hPa/L/min (<i>cm H₂O/L/min</i>)
Compressible volume:	AH102: 629 mL AH180: 1,474 mL
Compliance @ 60 hPa:	AH102: 1.8 mL/hPa (<i>mL/cm H₂O</i>) AH180: 2.2 mL/hPa (<i>mL/cm H₂O</i>)
Flow rate:	AH102/AH180: > 3 L/min
Rated flow:	AH102: 60 L/min AH180: 30 L/min

AirLife humidification chamber, auto-feed (Product code: AH290)

Duration of use:	Up to 14 days
Heater platform compatibility:	Compatible with Fisher & Paykel MR850 heater base
Compatibility with HFOV:	Compatible with HFOV
Interface connections:	ISO 5356-1 (22 mm male)
Compressible volume:	245 mL
Resistance to flow @ 60 LPM:	0.40 cm H ₂ O
Compliance:	0.6 mL/cm H ₂ O
Maximum operating pressure:	13.2 kPa
Maximum peak flow:	200 L/min for 3 seconds

For more information, please call Customer Service at **800.323.9088** or visit carefusion.com/AirLife

CareFusion
Vernon Hills, IL

carefusion.com

© 2016 CareFusion Corporation or one of its subsidiaries. All rights reserved. AirLife, CareFusion and the CareFusion logo are trademarks or registered trademarks of CareFusion Corporation or one of its subsidiaries. AL8006 (1216/3000)

