

Choose Vyaire™ Medical solutions for your skilled nursing facility needs.

Vyaire Medical is breathing new life into respiratory care.

While our name may be new to you, our products have been helping you save lives and provide care for over 65 years. Today, we manufacture and market more than 27,000 unique products for the diagnosis, treatment and monitoring of respiratory conditions across the continuum of care.

Vyaire was born of innovative brands that include AirLife®, Vital Signs®, Bear Medical Systems, Pulmonetic Systems, Sensormedics™, Viasys™ and many others; and we bring decades of industry leadership, clinical experience and unrivaled expertise.

Now, as a single global “breathing company,” we are focused on improving patient outcomes and increasing value for our customers.

LTV® ventilator series




The LTV® ventilator series, with greater than 85,000 units sold globally, has been providing ventilator support for pediatric to adults for since 1999. The LTV® ventilators provide portable, advanced ventilation for adult and pediatric patients at home or a post-acute care facility. At only 14 lbs, the ventilator provides a wide range of ventilation therapies to meet demanding patient needs, including volume control, pressure control, pressure support and spontaneous breath types.

Noninvasive ventilation (NIV)



Noninvasive ventilation (NIV) full face and nasal masks offer clinicians a reliable option, while providing patients comfort and security. The adjustable four-point headgear and flexible swivel connector enable ease in fit and adjustment. Masks feature a lightweight silicone design that contours to the patient's face for a comfortable fit, while the quick release clip eliminates the need to refit after removal. The odorless, easy-to-clean non-vented masks are available in three sizes to fit a wide range of patients.

Oxygen delivery



Our oxygen therapy disposables help meet nearly every oxygen therapy need across the care continuum. We offer an extensive range of AirLife® nasal cannulas, oxygen masks and accessories .

Suction kits



Optimize suctioning while reducing trauma. Our Tri-Flo™ catheters and catheter kits offer high suctioning efficiency in many catheter sizes and configurations to fulfill both clinician and patient needs.

Trach care kits



Sterile, single-use tracheostomy care and cleaning kits include complete trach cleaning trays and trach care kits. The trays include three rigid basins for cleaning, soaking and rinsing; pre-sewn trach dressing; and various other supplies for optimizing trach care.

HMEs



Lightweight heat and moisture exchangers (HMEs) and heat moisture exchanger with filter (HMEF) provide low resistance to flow and are treated to provide bacteriostatic and hygroscopic (HCH) benefits. Our portfolio also includes bypass heat and moisture exchangers as well as the Edith HME, which meets requirements for adult and pediatric patients in various critical care environments.

Modudose™ unit dose solutions



These easy-to-use unit dose solutions contain sterile water, normal saline or half-normal saline solution and are available in different sizes to accommodate your inhalation and irrigation needs.

Closed suction



Our closed suction catheters and Verso® adapters facilitate access to patient airways while helping to maintain continuous ventilation.

PureSom headgear and chinstrap



Keep a cool head with these CPAP accessories. PureSom headgear and chinstraps are made with Breathe-O-Prene®, a breathable material that allows air to flow, keeping your patients cool and comfortable. Our universal headgear can be used with our CPAP series and other vendor masks.

For more information, please contact us at **+1.833.327.3284**.

VYAIRE MEDICAL

26125 N. Riverwoods Blvd.
Mettawa, IL

vyaire.com

For U.S. distribution only.

Trademarks are the property of their respective owners.

© 2018 Vyair. Vyair, the Vyair Logo and all other trademarks are property of Vyair Medical, Inc. | CF/805/17/0138 (0118)

vyaire[™]
M E D I C A L