



WalkMed

TECHNOLOGIES



350VL Ambulatory Infusion Pump

For over 30 years, WalkMed has been a trusted provider of reliable and cost-effective ambulatory infusion equipment. We offer patient-focused solutions, including the WalkMed 350VL infusion pump, reservoir bags, administration sets, and carrying pouches, all designed for ease of use.

Continuous Infusion Therapy You Can Depend On



Flow Rate:

0.10 – 19.99 mL/hr
20.0 – 30.0 mL/hr

Total Volume Limit:

1 – 1,999 mL

Power:

9 volts Battery (Standard
or Rechargeable)

Typical Battery Capacity:

1 mL/hr – 18.7 Days
10 mL/hr – 2.7 Days
30 mL/hr – 0.7 Days

Dimensions:

4.4 x 4.0 x 1.5 in
(including front cover)

Weight: 12 oz

(including battery and
front cover)

Delivery Accuracy:

+/- 5%

Infusion Mode:

Continuous

Indications for Use:

The WalkMed 350VL pump is indicated for intravenous, subcutaneous, arterial, enteral, and epidural infusion of:
antibiotics, analgesics, chemotherapeutic agents, and other medications or fluids requiring precisely controlled infusion rates.



24 Hour Patient Helpline
(800) 578 - 0555