

By your side  
and beyond



New from Cardinal Health

Cardinal Health™  
Kangaroo OMNI™ ENTelliSet™  
Thick Formula Feeding Set

Discover the future of enteral feeding

The first and only feeding set of its kind, **Kangaroo OMNI™ ENTelliSet™ Thick Formula Feeding Set** delivers higher viscosity formulas. Thick formulas are defined as enteral fluids of smooth consistency that would be categorized as level 2, 3 or 4 drinks within the International Dysphagia Diet Standardisation Initiative (IDDSI) framework. These versatile sets are DEHP-, latex- and BPA-free.



E5FDR

Explore the advanced capabilities of Kangaroo OMNI™  
ENTelliSet™ Thick Formula Feeding Sets



**Hydration**  
Maintain tube patency and hydration when used with the Kangaroo OMNI™ ENTelliSet™ Thick Formula Feeding Set with Flush Bag



**Innovative**  
Delivers more viscous formulas



**Easy to use**  
Designed for use with ENFit™ connections and color-coded, patented cassettes that load with one hand



Ordering information

Catalog No.	Description	Volume	Ship case
E5FDR	Kangaroo OMNI™ ENTelliSet™ Thick Formula Feeding Set	500 mL	30
E5FFR	Kangaroo OMNI™ ENTelliSet™ Thick Formula Feeding Set with Flush Bag	500 mL	30
E10FDR	Kangaroo OMNI™ ENTelliSet™ Thick Formula Feeding Set	1,000 mL	30
E10FFR	Kangaroo OMNI™ ENTelliSet™ Thick Formula Feeding Set with Flush Bag	1,000 mL	30

For more information or to order, contact your sales representative or visit [cardinalhealth.com/kangaroo](http://cardinalhealth.com/kangaroo)