

By your side  
and beyond



New from Cardinal Health

# Cardinal Health™ Kangaroo OMNI™ Feeding Sets

Count on us for advanced nutritional delivery solutions

Our DEHP-, latex- and BPA-free **Kangaroo OMNI™ Feeding Sets** are designed to deliver commercially available formula to meet nutritional needs.



B5FD

Discover the value of Kangaroo OMNI™ Feeding Sets



**Accurate**

Pair with the Kangaroo OMNI™ Enteral Feeding Pump for +/- 5% pump accuracy



**Versatile**

Deliver a variety of enteral fluids, including commercially available formulas



**Easy to use**

Designed for use with ENFit™ connections and color-coded, patented cassettes that load with one hand



Ordering information

Catalog No.	Description	Volume	Ship case
<b>B5FD</b>	Kangaroo OMNI™ Feeding Set	500 mL	30
<b>B10FD</b>	Kangaroo OMNI™ Feeding Set	1,000 mL	30
<b>BSPFD</b>	Kangaroo OMNI™ Feeding Set with ENPlus Spike	-	30

For more information or to order, contact your sales representative or visit [cardinalhealth.com/kangaroo](http://cardinalhealth.com/kangaroo)