

250 mL

List No. 7951-12

RX ONLY

0

0

50

50

100

100

150

150

200

200

EMPTY CONTAINER

*SINGLE-USE
I.V. CONTAINER
FOR ADMIXING
COMPATIBLE MEDICATIONS.*

FLUID PATH AND AREAS BENEATH
UNDISTURBED PROTECTIVE COVERS
ARE STERILE AND NONPYROGENIC IN INTACT UNIT
PACKAGE. USE ASEPTIC TECHNIQUE. MIX CONTAINER
CONTENTS THOROUGHLY. ALL GRADUATIONS ARE
APPROXIMATE. LABEL APPROPRIATELY. BEFORE
ADMINISTRATION USE, CHECK FOR MINUTE LEAKS
BY SQUEEZING CONTAINER FIRMLY. DO NOT USE
IF CONTAINER IS DAMAGED. ADHERE TO STORAGE
REQUIREMENTS OF ADDED MEDICATIONS. MUST
NOT BE USED IN SERIES CONNECTIONS.

MADE IN COSTA RICA.

PRODUCT INQUIRIES SHOULD BE DIRECTED TO
ICU MEDICAL, INC., LAKE FOREST, ILLINOIS, 60045, USA



Not made with
natural rubber latex



CONTAINS DEHP
Recyclable as Type-3 Plastic

IM-4363 (06/2018)

icumedical